



---

# **Solute Attributes and Molecular Interactions Contributing to “U-Shape” Retention on Fluorinated HPLC Stationary Phases**

**David S. Bell and A. Daniel Jones  
The Pennsylvania State University  
Department of Chemistry**




# Introduction

---

- **LC/MS rapidly becoming the cornerstone of analytical chemistry**
  - Pharmaceutical
  - Environmental
  - Agricultural.....
- **Many efforts aimed at improving the technique**
  - Improve sensitivity
  - Greater speed
  - More universal
- **Much work centered on interfaces and instrumentation**
- **Less attention paid to LC**
- **Lack of retention/separation leads to**
  - Poor quantification
  - Problematic qualitative analysis


# Introduction

“U-Shape Retention”



# Introduction

- **Retention at high organic modifier percentages may:**
  - Increase sensitivity in LC/MS experiments through facilitated desolvation.





# Introduction

---

- **Retention at high organic modifier percentages may also:**
  - Speed up analyses.
  - Induce retention otherwise not obtained in RP.
  - Provide alternative mechanisms of interaction.
- **Needham (*J. Chromatogr. A*, 869 (2000) 159)**
  - Examined several stationary phase chemistries
  - Determined that pentafluorophenylpropyl (PFPP)-bonded silica provided greatest retention with good peak shape and reproducibility
- **Certain analytes under certain conditions using certain stationary phases.....**



# Introduction


---

- **Paramount to understand**
  - Solute attributes
  - Molecular interactions
  
- **In this research we sought to better understand “U-Shape” retention by:**
  - Investigating the relationship of retention and percent organic for several pharmaceutical acids, bases and neutrals on PFPP compared to C18.
  - Performing several studies aimed at elucidating the dominant molecular interactions responsible for the observed behavior.

# Introduction

---

- **PFPP offers:**
  - Dipole-dipole interactions
  - Pi-pi interactions
  - Charge-transfer interactions
  - Others?



# Acidic, Basic and Neutral Retention Profiles


---

- **Experimental Design:**

- 6 analytes representing pharmaceutical acids, bases and neutrals
- Retention monitored from 40 to 90% acetonitrile at two pH levels (4 and 6.7) on Discovery<sup>®</sup> HS F5 (PFPP) and Discovery HS C18
- Aqueous component-10 mM ammonium acetate

# Acidic Analyte Profiles on PFPP


Retention Profiles ( $k'$ ) of Acidic Probes on PFPP at pH 6.7



Retention ( $k'$ ) of acidic probes x ibuprofen, ☎ aspirin, + naproxen, ✉ ketoprofen, 📠 piroxicam and 🍷 diclofenac using PFPP from 40% to 90% acetonitrile under pH 6.7 conditions.

# Basic Analyte Profiles on PFPP


Retention Profiles ( $k'$ ) of Basic Probes on PFPP at pH 6.7








Retention ( $k'$ ) of basic probes ☐ amitriptyline, 📞 nortriptyline, 📺 diphenhydramine, x verapamil, ✉ alprenolol and ☎ lidocaine using PFPP from 40% to 90% acetonitrile under pH 6.7 conditions

# Neutral Analyte Profiles on PFPP


Retention Profiles ( $k'$ ) Neutral Probes on PFPP at pH 6.7



Retention ( $k'$ ) of neutral probes + hydrocortisone,  hydrocortisone acetate,  progesterone,  corticosterone,  cortisone acetate and  prednisone using PFPP from 40% to 90% acetonitrile under pH 6.7 conditions

# Basic Analyte Profiles on C18

Retention Profiles ( $k'$ ) of Basic Probes on C18 at pH 6.7




Retention ( $k'$ ) of basic probes amitriptyline, nortriptyline, diphenhydramine, verapamil, alprenolol and lidocaine using C18 from 40% to 90% acetonitrile under pH 6.7 conditions


# Selectivity Analysis

- Secondary amines show relatively greater increase in the normal-phase region
- Except for verapamil, amines are located on structural “arms”
- Basic moiety in verapamil hindered
- Lidocaine shows evidence of a pKa dependence

Amitriptyline



Verapamil





# Retention Profiling Conclusions

---

- **Only basic analytes exhibit appreciable “U-Shape” retention**
- **Not all bases**
- **Not all to the same extent - selectivity**
- **Also observed to a small degree on C18**
  - Commonality is the silica surface
  - Known to interact with bases via ion-exchange
- **Pointed to ionic interactions with surface silanols**


# **“U-Shape” Retention Profile Dependence on MP Ionic Strength**

---

- **As % organic increased, ionic strength decreased**
- **What if the ionic strength was held constant?**
- **Experimental:**
  - Basic analytes nortriptyline and amitriptyline run from 50 to 90% acetonitrile on PFPP keeping the buffer concentration at 10 mM throughout

# “U-Shape” Retention Profile Dependence on MP Ionic Strength

Comparison of Amitriptyline and Nortriptyline Retention on PFPP Phase: Decreasing vs. Constant Buffer Concentration



- Retention at all % acetonitrile levels is attenuated
- Little “U-Shape” character observed at constant ionic strength

# Contribution of Bonded Phase to Retention

---


- **Compared basic analyte retention on bare silica to PFPP**
  - 2 mM ammonium acetate (pH 6.7) at 90% acetonitrile

Compound	k' Silica	k' PFPP
amitriptyline	4.63	12.71
nortriptyline	6.92	18.98
diphenhydramine	4.69	10.98
verapamil	1.76	5.40
alprenolol	6.09	12.59
lidocaine	0.09	0.46

# Contribution of Bonded Phase to Retention

kappa-kappa plot: Retention of Basic Analytes on PFPP Versus Bare Silica

- good correlation indicative of similar dominant interactions



# Dependence of Retention on Ionic Strength at Constant % Acetonitrile

---

- **Two-site model of retention** (*Yang, et. al., J Chromatogr. A 996 (2003) 13*)

$$\log k' = \log(k'_{RP} + B_{IEX}/[C^+]_m)$$

Reversed-Phase Contribution  
(independent of salt concentration)

Ion-Exchange Component  
(inversely proportional to salt concentration)

Plot of  $\log k'$  vs  $\log[C^+]_m$

slope = -1, if exclusively ion-exchange


slope =  $\sim 0$ , if RP dominates

# Dependence of Retention on Ionic Strength at Constant % Acetonitrile

Amitriptyline Retention ( $\log k'$ ) vs. Ammonium Acetate Concentration ( $\log[\text{mM}]$ ) on PFPP at 85% Acetonitrile

- Slope of  $-0.6659$  indicates more than just ion-exchange


- Note linear dependence of retention on buffer conc.



# Dependence of Silanol pKa on Bonding Chemistry


---

- For IEX silanols must be ionized
- Neue (*J. Chromatogr. A* 925 (2001) 49) used quaternary ammonium ion (bretylum) retention as a function of pH to estimate silanol pKa values.
- Bare silica, PFPP and C18 were subjected to similar experiment
  - Ammonium ion concentration held constant at 25 mM
  - pH varied from 2 to 8
  - Bretylum ion retention monitored



# Dependence of Silanol pKa on Bonding Chemistry

Retention ( $k'$ ) of Bretylium Ion as a Function of pH



pH 2 to 8 on bare silica, PFPP and C18

- Bonded-phase alters silanol acidity
- Bonded phase inhibits interaction
- Explains why PFPP and not C18



# Summary

- **Only basic analytes exhibit appreciable “U-Shape” retention on PFPP**
  - Also observed on C18 to a lesser extent
  - Selectivity and good peak shape are obtained
  - Selectivity appears to be a function of pKa and silanol/base accessibility
- **The “normal-phase” region of the profile was shown to be a function of MP ionic strength**
- **Based on the two-site model of retention, both RP and IEX contribute to retention (**hydrophobically-assisted ion-exchange**)**
- **Further supported by the preferential retention of bases on PFPP compared to silica**
- **Estimated pKa values for C18 and PFPP also support the existence of IEX mechanisms**





# Conclusions

---

- **Ability to retain basic analytes at high organic modifier concentrations offers potential to increase LC/MS sensitivity**
- **Alternative mechanisms may provide unique options in method development**
- **Bonded-phase chemistry plays an important role in both RP and IEX**



# Conclusions

---

- **PFPP presents new opportunities to manipulate retention and selectivity**
  - Buffer concentration
  - pH (pKa of analytes and surface silanol groups)
- **Knowledge should lead to:**
  - greater utilization
  - more robust methods
  - further advances in stationary phase chemistries
  - expansion to other forms of separation

# Acknowledgements

---

- **A. Daniel Jones (Penn State)**
- **Keith Duff (Supelco) and Shane Needham (Alturas Analytics)**
- **Supelco**
- **Exygen Research**
- **Organizers**

