



Mobius[®] Power MIX 2000 Single-Use Mixer

High performance mixing technology and ergonomic design features.

The **Mobius® Power MIX 2000** combines high performance mixing technology with design features that make it easy to use. The impeller design and motor are based on our magnetically coupled NovAseptic® mixing technology, a proven mixing technology in stainless steel tanks. The system efficiently mixes the most challenging buffers, media and biopharmaceutical ingredients.

Features and benefits:

- Powerful magnetically coupled impeller designed to efficiently mix the most challenging buffer, media and biopharmaceutical ingredients
- Broad working volume, from 300 L to 2200 L, for processing flexibility
- Double doors that open to one half the mixer's diameter and easily accessible off-set drain for easy bag loading
- Optional integrated hoist for powder delivery
- Integrated load cells and temperature sensor
- Two in-process pH measurement options: Probe Port for easy insertion of a re-usable standard PG 13.5 threaded probe or installed, pre-calibrated and pre-sterilized Hamilton OneFerm[®] Single-Use pH Sensor
- Closed, sterile, and zero deadleg sampling directly from the mixing bag
- Backed by industry-leading technical support and services





Specifications

| Specifications | | | |
|---|--|--|--|
| Materials of Construction | | | |
| Carrier | | | |
| Tank | 304L stainless steel | | |
| Stand | 304L stainless steel | | |
| Latch | 18-8 stainless steel | | |
| Casters | 304L stainless steel | | |
| Dimensions LXW | 63.75 x 63.75 inches (1619,3 x 1619,3 mm) | | |
| Height without hoist | 86.7 in (2210mm) | | |
| Height with hoist | 133.2 in (3384mm) | | |
| Weight | lb (kg) | | |
| With hoist and lid | 2100 lb (953kg) | | |
| Fluid volume in jacket (approx.) | 8 gal (30L) | | |
| Heat transfer area in jacket (approx.) | 57.7 (5.4) sq. ft. (5.36 sq. m.) | | |
| Maximum working pressure | 150 psi | | |
| Stainless Steel Jacketed Mixers ASME U-1 code stamped | Optional | | |
| Integrated Load Cells | Mettler-Toledo 0970 RingMount Load Cells | | |
| Ingress protection rating | IP68 | | |
| RTD Temperature Sensor OMEGA RTD 809 | Standard | | |
| Storage | 0 - 40 °C non-condensing environment | | |
| Operating temperature | 4 - 60 °C non-condensing environment | | |
| Motor | | | |
| Weight lb (kg) | 30.4 lb (13.8 kg) | | |
| Electrical (max) | 0.55 kW | | |
| Design speed (max) | 380 rpm | | |
| Maximum operating speed has to be validated for e | each process and is dependent on volume, viscosity, etc. | | |
| Cable length | 12 ft (3.65m) | | |
| Ingress protection rating | IP55 | | |
| Regulatory information | CSA, CE, UL | | |
| Control Box | | | |
| Materials of construction | 304L stainless steel | | |
| Weight lb (kg) | 21 lbs (9.5kg) | | |
| Electrical | 120V, 60Hz,15A or 230V, 50Hz, 10A | | |
| Cable length in (cm) | 15 ft (4.57m) | | |
| Ingress protection rating | IP 65 | | |
| Regulatory information | CE, UL Compliant | | |
| Display Box (with temperature sensor and load | cell display) | | |
| Materials of construction | 304L stainless steel | | |
| Weight | Permanently mounted to carrier | | |
| Electrical | 120 / 230 V | | |
| Cable length in (cm) | 15 ft (4.57m) | | |
| Ingress protection rating | IP55 | | |
| Regulatory information | CE, UL compliant | | |
| | | | |

Specifications

| Materials of Construction | |
|---------------------------|--|
| Mixer Container Assembly | |
| Mix container | PureFlex™ Plus Film |
| Adapters | SterilEnz™ LDPE |
| Luer plugs | Kynar [®] Polyvinylidene Fluoride |
| Connectivity | CPC [®] Polysulfone |
| Tubing | Silicone |
| Pinch clamp | PVDF |
| Impeller | Neodymium magnets encased in polypropylene |
| Powder port | Polyethylene |
| TC gasket | Silicone |
| TC clamp | Nylon |
| TC end cap | Polypropylene |
| Cup | HDPE |
| Probe port | HDPE |
| O-ring | Silicone |
| Plug | Polypropylene |
| Needle-free sample port | HDPE |
| O-ring | Silicone |

Ordering Information

| | | Top Ports | | | |
|--|---------------|-----------------------|----------------------|--------|--------------|
| Description | Certification | 4 inch Powder port | ³ ⁄4 inch | 1 inch | Part Code |
| Container Assembly | | | | | |
| Mobius® Power MIX 2000 Assembly, Bronze, PureFlex™ Plus | Bronze | 1 | 1 | 2 | MXRA2000BPL1 |
| Mobius® Power MIX 2000 Liner, Bronze, PureFlex™ Plus | Bronze | n/a | n/a | n/a | MXRL2000BPL1 |
| Mobius [®] Power MIX 2000 Assembly, Gold, PureFlex [™] Plus | Gold | 1 | 1 | 2 | MXRA2000GPL1 |

Please note: All mixer assemblies are provided as one packs.

Ordering Information

| Description | Qty/Pk | Catalog No. |
|--|--------|--------------|
| Mobius [®] Power MIX 2000 Carriers | | |
| 2000 L jacketed carrier, hoist, temp sensor, load cells, ASME cert. | | MXRJ2000HTLA |
| 2000 L jacketed carrier, hoist, temp sensor, load cells | | MXRJ2000HTL |
| 2000 L jacketed carrier, temp sensor, load cells, ASME cert. | | MXRJ2000TLA |
| 2000 L jacketed carrier, temp sensor, load cells | | MXRJ2000TL |
| Mobius [®] Power MIX Electronics | | |
| Motor | | MXRMTR1K2K |
| Control box, 120V | | MXRCTRLB2K |
| Control box, 230V | | MXRCTRLB2K2 |
| Electrical power cords (country specific) | | |
| North America, 120 V (15A) | | MXRPCNA1K2K |
| North America, 120 V (20A) | | MXRPCNA125 |
| Europe, 230 V (10A) | | MXRPCEU |
| Argentina, Australia, Colombia, New Zealand, Paraguay, Uruguay | | MXRPCAACNPU |
| Denmark | | MXRPCDN |
| Israel | | MIXRPCIS |
| Chile, Ethiopia, Italy, 250V (10A) | | MXRPCCEI |
| Bangladesh, Burma, Pakistan, South Africa, Sri Lanka | | MXRPCBBPSS |
| Switzerland | | MXRPCSZ |
| Display have with tampageture senses and land call displays i | | |

Display box with temperature sensor and load cell displays is standard and provided with the carrier.

Please note: With the exception of NA, two identical power cords must be purchased - one power cord for the Control Box and one power cord for the Display Box.

In NA, the two power cords that need to be purchased are MXRPCNA1K2K and MXRPCNA125.

| Spare Parts | |
|--|---------------|
| Quick-connect coupling, male, 34 inch NPT | MIXSST0109442 |
| Quick-connect coupling, female, ³ / ₄ inch NPT | MIXSST0109441 |

To place an order or to receive technical assistance

In the U.S. and Canada, call toll-free **1-800-645-5476**

For other countries across Europe and the world, please visit: emdmillipore.com/offices

For Technical Service, please visit: emdmillipore.com/techservice



EMDMillipore.com

MilliporeSigma and the vibrant M are trademarks and Mobius is a registered trademark of Merck KGaA, Darmstadt Germany. All other trademarks are property of their respective owners. © 2016 EMD Millipore Corporation. All Rights Reserved.

Lit#SP1826EN00 Ver. 2.0 PS-16-12836 09/2016