

Identification of Glycyrrhizic Acid in Licorice According to USP and Quantification Using the TLC Explorer

Valanka D'Silva

R&D APAC Lab, Jigani, Bangalore, India

Abstract

A thin-layer chromatographic method for the identification of glycyrrhizic acid in licorice was established in accordance with the United States Pharmacopeia (USP) monograph. The TLC Explorer system was utilized for improved digital documentation, automated assessment, and quantitative analysis, thereby broadening the use of TLC for this application. Effective chromatographic separation and documentation were achieved, and glycyrrhizic acid was successfully detected in licorice samples. The results demonstrated that the TLC Explorer system provides a reliable platform for TLC-based analysis and documentation.

Introduction

Plants are known to produce various secondary metabolites that serve a wide variety of medicinal functional beneficial to humans. Licorice (*Glycyrrhiza glabra*, *G. uralensis*) has been extensively used in Eastern and Western medicine for the treatment of conditions, including the common cold, fever, liver ailments, gastric ulcers, asthma, bronchitis, Addison's disease, and rheumatoid arthritis. Additionally, its use as a laxative, antitussive, and expectorant has also been documented.¹

Glycyrrhizic acid (**Figure 1**), also known as glycyrrhizin, is a triterpenoid saponin found in the roots, rhizomes, and stolons of licorice plant. In the gastrointestinal tract, glycyrrhizin is metabolized to glycyrrhetic acid, which has been identified as the primary contributor to the pharmacological and biological activities associated with licorice.²




Figure 1. Chemical structure of glycyrrhizic acid (glycyrrhizin).

Thin-Layer Chromatography (TLC) is designated in the United States Pharmacopeia (USP) monograph for licorice as an identification procedure.³ TLC-based approaches are widely cited in pharmacopeial identity tests. High-Performance Thin-Layer Chromatography (HPTLC), an advanced form of TLC commonly paired with automation, has been established as a robust, reliable, rapid, and cost-effective technique for both qualitative and quantitative analysis of medicinal compounds. Chromatographic fingerprints generated by this technique can be visualized and archived as electronic images.^{4,5}

In the present study, the TLC identification test for glycyrrhizic acid prescribed in the USP monograph was executed using the TLC Explorer documentation system (**Figure 2**). The system enabled not only standardized digital documentation but also a quantitative assessment of glycyrrhizic acid content to evaluate compliance with the USP requirement for licorice (NLT 2.5% on the dried basis). Three commercially available licorice brands were examined.



Figure 2. TLC Explorer.

The TLC Explorer documentation system enables digital and automated evaluation of TLC plates, thereby enhancing the efficiency and accuracy of thin layer chromatography analyses. Three illumination modes are provided through LED light sources—white light (VIS), UV-A (366 nm), and UV-C (254 nm) — to facilitate the detection and fast measurement of the compounds of interest. The system software offers special features like automated track recognition, the ability to evaluate multiple plates simultaneously, and background signal correction. These capabilities allow accurate TLC imaging suitable for video-densitometric measurements, supporting quantitative analysis and consistent data interpretation. (Read more at [SigmaAldrich.com/tlc-explorer](https://www.sigmaaldrich.com/tlc-explorer)).

Experimental

Reagent Preparation

Mobile phase: Mix butyl alcohol, acetic acid glacial and water in a ratio of 7:1:2, v:v:v.

Diluent: Mix ethyl alcohol and water in a ratio of 7:3, v:v.

Standard Preparation

Standard solution 1 for identification (5 mg/mL): Weigh and dissolve 50 mg of glycyrrhizic acid RS in 10 mL of diluent.

Standard solution 2 for quantification (0.125 mg/mL): Weigh and dissolve 50 mg of glycyrrhizic acid RS in 100 mL of diluent. Dilute 2.5 mL of this solution to 10 mL with diluent.

Sample Preparation

Samples: Three different commercially available licorice brands.

Test solutions for identification I + II + III: Weigh and add 2 g of pulverized licorice to 10 mL of diluent. Heat the solution by shaking on a water bath for 5 minutes, cool and filter through a 0.45 µm PVDF syringe filter.

Test solutions for quantification IV + V +VI: Weigh and add 500 mg of pulverized licorice to 10 mL of diluent. Heat the solution by shaking on a water bath for 5 minutes, cool and filter through a 0.45 µm PVDF syringe filter. Dilute 1 mL of this filtrate to 10 mL with diluent.

Instrument Parameters

Table 1. TLC conditions

TLC parameters	
Plate:	TLC aluminum sheets, Silica gel 60 F ₂₅₄ 20 x 20 cm (1.05554)
Sample application:	2 µL of each solution
Mobile phase:	Butyl alcohol:acetic acid glacial:water (7:1:2 v:v:v)
Chamber conditions:	Twin trough unsaturated chamber
Migration distance:	10 cm
Drying:	Air-drying
Detection:	UV at 254 nm

Results

Identification

The identification of glycyrrhizic acid in licorice, performed according to the USP monograph and visualized at 254 nm using the TLC Explorer system, is shown in **Figure 3**. A dark purple zone corresponding to glycyrrhizic acid was observed for both the standard solution and the test solutions under UV 254 nm, in addition to other sample constituents. The retention factors (R_f) obtained from these analyses are summarized in **Table 2**.




Figure 3. TLC chromatogram showing the identification of glycyrrhizic acid under UV 254 nm using the TLC Explorer.

Table 2. Chromatographic data obtained for the identification of glycyrrhizic acid under UV 254 nm using the TLC Explorer

Tracks no.	Solution name	R_f of Glycyrrhizic acid
1 & 2	Standard solution 1	0.381
3 & 4	Test solution I	0.389
5 & 6	Test solution II	0.388
7 & 8	Test solution III	0.391

Quantification

The quantification of glycyrrhizic acid in three commercial licorice brands was carried out at 254 nm (Figure 4) using the TLC Explorer system with video densitometry. A one-point calibration was applied using the 0.125 mg/mL standard solution, which corresponds to a sample concentration equivalent to 2.5% under the conditions employed. Table 3 summarizes the obtained chromatographic results showing that all three licorice samples contained glycyrrhizic acid at levels exceeding the USP acceptance criterion (NLT 2.5%).



Figure 4. TLC chromatogram for the quantification of glycyrrhizic acid under UV 254nm as documented by the TLC Explorer system.

Table 3. Chromatographic data and quantitative results for glycyrrhizic acid in three licorice samples obtained by video densitometry using the TLC Explorer at 254nm

Tracks no.	Solution name	R_f Glycyrrhizic acid	Peak area	% Content
1 & 2	Standard solution 2	0.394	55.4	2.50
3 & 4	Test solution IV	0.401	153.8	6.94
5 & 6	Test solution V	0.406	95.8	4.32
7 & 8	Test solution VI	0.406	130.1	5.87

Conclusion

The USP monograph-specified identification test for glycyrrhizic acid in licorice was carried out using thin-layer chromatography, and the chromatographic assessment and documentation were performed with the TLC Explorer system. The principal spots observed in the chromatograms of the test solutions exhibited R_f values consistent with the principal spot in the standard chromatogram, fulfilling the requirements of the monograph. In addition to the identity test, the method was extended to a quantitative evaluation of glycyrrhizic acid using the TLC Explorer's video-densitometric capabilities. The glycyrrhizic acid content of the examined licorice samples was determined to range from 4.32% to 6.94%, exceeding the USP acceptance threshold of not less than 2.5%.

Overall, the study demonstrated that the TLC Explorer system provided reliable chromatographic documentation, track identification, R_f calculation, and video-densitometric quantification for the analysis of glycyrrhizic acid in licorice.

References

1. Isbruckner RA and Burdock GA. Risk and safety assessment on the consumption of Licorice root (*Glycyrrhiza* sp.), its extract and powder as a food ingredient, with emphasis on the pharmacology and toxicology of glycyrrhizin. *Regulatory Toxicology and Pharmacology*. 2006;46(3):167-192. <https://doi.org/10.1016/j.yrtph.2006.06.002>
2. Chauhan S, Gulati N, Nagaich U. Glycyrrhizic acid: extraction, screening and evaluation of anti-inflammatory property. *Ars Pharm*. 2018;59(2):61-67. <http://dx.doi.org/10.30827/ars.v59i2.7513>
3. USP monograph: Licorice. https://doi.org/10.31003/USPNF_M45050_02_01
4. Attimarad M, Ahmed KK, Aldhubaib BE, Harsha S. High-performance thin layer chromatography: A powerful analytical technique in pharmaceutical drug discovery. *Pharm Methods*. 2011;2(2):71-75. <https://doi.org/10.4103/2229-4708.84436>
5. Srivastava MM. An overview of HPTLC: A modern analytical technique with excellent potential for automation, optimization, hyphenation, and multidimensional applications. In: *High-Performance Thin-Layer Chromatography (HPTLC)*. Springer Berlin Heidelberg; 2011. p 3-24. https://doi.org/10.1007/978-3-642-14025-9_1

Related Products

Description	Cat. No.
TLC Analysis & Documentation	
TLC Explorer, digital TLC analysis and documentation device	1.52610
TLC aluminum sheets, Silica gel 60 F ₂₅₄ 20 x 20 cm, Pk. 2	1.05554
Solvents, Reagents, Standards, and Consumables	
1-Butanol, suitable for HPLC, ≥99.7%	34867
Acetic acid glacial, for chromatography	6.18665
Water, suitable for HPLC	270733
Ethyl alcohol, pure, ≥99.5% (GC), EMPARTA® ACS reagent, for extraction and analytical purposes	1.07017
Millex™ HV Durapore (PVDF) 0.45µm	SLHV033N
Glycyrrhizic acid, Pharmaceutical Secondary Standard; Certified Reference Material	PHR1516

Find more information on our TLC portfolio at

[SigmaAldrich.com/tlc](https://www.sigmaaldrich.com/tlc)

More details on the TLC Explorer Documentation System are available at

[SigmaAldrich.com/tlc-explorer](https://www.sigmaaldrich.com/tlc-explorer)

MilliporeSigma
400 Summit Drive
Burlington, MA 01803

To place an order or receive technical assistance

Order/Customer Service: [SigmaAldrich.com/order](https://www.sigmaaldrich.com/order)

Technical Service: [SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice)

[SigmaAldrich.com](https://www.sigmaaldrich.com)

