


A Case Study: 3-Step Process for Efficient mAb Purification Using Different Commercially Available Chromatography Resins

Melissa Holstein*, Kristen Cotoni, Nanying Bian, EMD Millipore Corporation, 80 Ashby Road, Bedford, MA 01730

Summary

This study showcases a portfolio of commercially available biopharmaceutical chromatography resins designed for the efficient purification of monoclonal antibodies. A 3-step purification process has been implemented which showed effective removal of the main contaminants, low ligand leakage, and high yields over the entire process. Eshmuno® A affinity chromatography resin was evaluated as the first step in the process. The Protein A elution pool was further purified using cation exchange chromatography. Two cation exchange resins with different selectivities were compared. The final purification step consisted of anion exchange chromatography. The product yield, HCP removal, leached Protein A removal, and aggregate content were evaluated at each step in the process. Final process yields ranged from 74 to 81 percent. This work highlights the comparable purification capabilities of a panel of EMD Millipore resins designed for efficient mAb purification.



Results

mAb purity was evaluated based on Host Cell Proteins (HCP) and leached Protein A after each purification step. The cation exchange elution fractions were pooled in order to maximize purity while maintaining at least 80% step yield.

Comparison of Leached Protein A Levels

Leached Protein A: Repligen Protein A Kit 9000-1

Yield: Based on UV absorbance at 280 nm


3G Immunoenzymetric Assay


mAb Feed HCP: 281,190 ppm

HCP: Cygnus Technologies F550 CHO HCP ELISA Kit

this study. The clarified feed was sterilized using 0.22

µm Stericup® filters immediately prior to use.





Purity and Yield for the Two Process Trains

Neutralized the Protein A elution pool to pH 5 using 2 M Tris

• Elution: 0.1 M acetic acid, pH 3

1	HCP (ppm)	HCP log removal (cumulative)	Protein A (ppm)	Protein A log removal (cumulative)	Step Yield (%)
PrA - Eshmuno® A	212	3.1	35	N/A	97.8
CEX - Fractogel® SO ₃	81	4.1	7.8	0.65	86.1*
AEX - Eshmuno® Q	ND	> 5	8.8	0.60	87.3
Overall					73.5

Anion Exchange (Flow-through mode):

• Eshmuno® Q resin: 25 mM Tris, pH 7.3


2	HCP (ppm)	HCP log removal (cumulative)	Protein A (ppm)	Protein A log removal (cumulative)	Step Yield (%)
PrA - Eshmuno® A	283	3.0	35	N/A	97.8
CEX - Eshmuno® CPX	12	4.4	10	0.55	89.0*
AEX - Eshmuno® Q	2	> 5	12	0.62	93.0
Overall					81.0

ND = Not Detectable (<1 ppm); * After pooling fractions

Conclusion

A 3-step purification process has been implemented which showed effective removal of the main contaminants, low residual leached Protein A, and high yields over the entire process.

The two different combinations of commercially available biopharmaceutical resins from EMD Millipore used in this case study gave comparable yields and final mAb impurity profiles.



Corresponding Author: melissa.holstein@emdmillipore.com

EMD Millipore, ChromaSorb and the M mark are trademarks of Merck KGaA, Darmstadt, Germany. Eshmuno, Fractogel, Opticap, Millipore Express, Millistak+ and Stericup are registered trademarks of Merck KGaA, Darmstadt, Germany. 08/13 Lit. No. PS4218EN00 Rev.B © 2013 EMD Millipore Corporation, Billerica MA USA. All rights reserved.

EMD Millipore is a division of Merck KGaA, Darmstadt, Germany