




Analysis of Banned and Other Neonicotinoids from Plant Material

Back in 2018, the representatives of the Member States of the European Union decided to ban three products (clothianidin, imidacloprid, thiamethoxam) of the neonicotinoid family of insecticides deemed dangerous for bees on all outdoor crops permitting only greenhouse use.

On Monday, January 15, 2020, the European Commission announced that they would not renew the license for a fourth neonicotinoid, thiacloprid, based on a scientific opinion from the European food Safety Authority (EFSA). One of the major concerns was the impact on groundwater.

To help you ensure accurate testing of these pesticides, we offer a comprehensive portfolio of analytical standards and certified reference materials for the whole range of neonicotinoids, their main metabolites and isotope-labeled products.

The below application demonstrates the extraction and cleanup of these banned compounds and others of the neonicotinoid class from dandelion blossoms using the QuEChERS methodology.



For the QuEChERS approach, Supel[™] QuE tubes were used for extraction (Acetate) and cleanup (PSA/C18) prior to analysis of seven neonicotinoid pesticides. For the LC-MS analysis, an Ascentis[®] Express C18 Fused-Core[®] particle column was chosen to achieve high efficiency at relatively low backpressure compared to sub-2 µm UHPLC columns. This allowed the analysis to be performed on a standard pressure HPLC system with sufficient efficiency.



Conditions

| sample/ matrix | 3 g pulverized dandelion blossoms, homogenized in 10 mL water | |
|-----------------------|--|--|
| extraction process | add 25 mL of acetonitrile; add contents of Supel [™] QuE Acetate extraction tube; shake immediately for 1 minute; centrifuge at 3400 rpm for 5 min | |
| clean-up process | transfer 1 mL of the acetonitrile layer into a Supel [™] QuE PSA/C18 cleanup tube; shake for 1 minute; centrifuge at 3500 rpm for 3 minutes; (draw off 700 µL of supernatant, evaporate to dryness at 50 °C under nitrogen; reconstitute in 200 µL of 50:50 0.1% formic acid:0.1% formic acid in methanol) | |
| column | Ascentis® Express C18, 10 cm x 3.0 mm I.D., 2.7 μm particles | |
| column temp. | 25 °C | |
| mobile phase | [A] 0.1% formic acid in water; [B] 0.1% formic acid in methanol | |
| gradient | 30% B from 0 to 5 min; to 100% B in 0.2 min; held at 100% B for 5.3 min; to 30% B in 0.5 min; held at 30% B for 5 min | |
| flow rate | 500 µL/min | |
| pressure | 3800 psi (262 bar) | |
| injection | 3 μL | |
| detector | MS, ESI(+), MRM, m/z: 203.2/129.2, 271.2/225.0, 292.1/211.0, 256.0/175.2, 250.0/132.0, 223.2/126.0, 253.0/125.8 | |

Materials

| Product | Description |
|---------|---|
| 55234-U | Supel [™] QuE Acetate extraction tube, AOAC 2007.01, pk of 50 |
| 55288-U | Supel [™] QuE PSA/C18 cleanup tube, 2mL, AOAC 2007.01, pk of 100 |
| 52248-U | Supel [™] QuE Empty 50 mL Centrifuge Tube with Lid, pk of 50 |
| 53814-U | Ascentis [®] Express C18, 10 cm x 3.0 mm I.D., 2.7 µm particles |
| 1.00029 | Acetonitrile hypergrade for LC-MS LiChrosolv® |
| 1.06035 | Methanol hypergrade for LC-MS LiChrosolv [®] . |
| 1.15333 | Water for chromatography (LC-MS Grade) LiChrosolv $^{\mbox{\tiny 8}}$ |
| 5.33002 | Formic acid 98% - 100% for LC-MS LiChropur^ $\!\!^{\scriptscriptstyle (\!8\!)}$ |

For more information about neonicotinoids visit us at

SigmaAldrich.com/neonicotinoids

Standards

| Product | Description | Package Size | | |
|--------------|---|-----------------|--|--|
| Acetamiprid | | | | |
| 33674 | N-(6-Chloro-3-pyridylmethyl)-N-cyano-N- methylacetamidine Pestanal® analytical standard | 100 mg | | |
| 51625 | N-(6-Chloro-3-pyridylmethyl)-N-cyano-N- methylacetamidine TraceCERT [®] certified reference material | 50 mg | | |
| 39246 | d_3 -Acetamiprid, analytical standard | 50 mg | | |
| 32979 | Acetamiprid-N-desmethyl PESTANAL® analytical standard | 10 mg | | |
| Clothianidin | | | | |
| 33589 | Clothianidin PESTANAL® analytical standard | 100 mg | | |
| 68965 | Clothianidin TraceCERT® certified reference material | 50 mg | | |
| 56816 | $d_{\rm 3}\mbox{-}Clothianidin \mbox{PESTANAL}^{\rm @}$ analytical standard | 50 mg | | |
| Imidacloprid | | | | |
| 37894 | Imidacloprid PESTANAL® analytical standard | 100 mg | | |
| 68694 | TraceCERT [®] certified reference material | 50 mg | | |
| 46341 | Imdiacloprid solution, 100µg/mL, PESTANAL [®] , analytical standard | 2 mL | | |
| 34170 | d ₄ -Imidacloprid PESTANAL [®] analytical standard | 10 mg | | |
| 34534 | $Imidacloprid-olefin \ PESTANAL^{\circledast} \ analytical \ standard$ | 10 mg | | |
| 37052 | Desnitro-imidacloprid hydrochloride PESTANAL® analytical standard | 25 mg | | |
| Thiacloprid | | | | |
| 37905 | Thiacloprid PESTANAL® analytical standard | 100 mg | | |
| 14783 | Thiacloprid TraceCERT [®] , certified reference material | 50 mg | | |
| 33897 | Thiacloprid-amide PESTANAL® analytical standard | 100 mg | | |
| 30673 | Thiacloprid-(<i>thiazolidin ring-d</i> ₄) analytical standard | 10 mg | | |
| Thiamethoxam | | | | |
| 37924 | Thiamethoxam PESTANAL® analytical standard | 100 mg | | |
| 68901 | Thiamethoxam TraceCERT [®] , certified reference material | 50 mg | | |
| 38176 | $d_{3}\mbox{-}Thiamethoxam \mbox{PESTANAL}^{\mbox{\tiny (B)}}$ analytical standard | 25 mg | | |
| 73348 | N-Desmethylthiamethoxam Pestanal® analytical standard | 50 mg | | |
| Metabolites | | | | |
| 19386 | 6-Hydroxypyridine-3-carboxylic acid analytical standard | 100 mg | | |
| 68678 | 6-Chloropyridine-3-carboxylic acid analytical standard | 100 mg | | |
| 31534 | 2-Imidazolidone PESTANAL® analytical standard | 250 mg | | |
| Others | | | | |
| 46077 | Nitenpyram PESTANAL® analytical standard | 100 mg | | |
| 32499 | Dinotefuran PESTANAL® analytical standard | 50 mg | | |
| 69091 | Dinotefuran reference material | 50 mg | | |

MilliporeSigma

400 Summit Drive Burlington, MA 01803

To place an order or receive technical assistance

Order/Customer Service: SigmaAldrich.com/order Technical Service: SigmaAldrich.com/techservice Safety-related Information: SigmaAldrich.com/safetycenter

SigmaAldrich.com

© 2020 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. MilliporeSigma, the vibrant M, Ascentis, Supel, PESTANAL, TraceCERT, LIChrosolv and LiChropur are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

