


### **Trusted Results**

SupraSolv® High Purity Solvents for Gas Chromatography







## SupraSolv® and SupraSolv® hypergrade Solvents for gas chromatography

As the world's leading supplier of high-purity solvents, we offer a full range of products for every gas chromatography application in the laboratory - including highly sensitive pesticide and dioxin analyses. Our SupraSolv® and SupraSolv® hypergrade solvents are developed specifically for sensitive detection processes in residue and environmental analysis. They cover all areas of application, and provide the highest level of reliability for your analytical results.

To ensure purity and suitability, we employ only the latest manufacturing processes. SupraSolv® solvents are recovered during special distillation cuts, and suitability testing involves a variety of detectors and highly concentrated solvents.

As a result, these high-purity products support you in countless ways during your daily work - with individual specifications that are tailored to their specific area of application.

### SupraSolv® and SupraSolv® hypergrade solvents benefits

- The most comprehensive application range due to the largest retention time range
- Analytical reliability due to the highest possible purity and a minimal signal-to-noise ratio
- Time and cost savings due to the best possible batch consistency, thus avoiding analysis repetition



### Table of Contents

| gas chromatography  | 2 |
|---|---|
| SupraSolv® and SupraSolv® hypergrade solvents for all your gas chromatography needs | 4 |
| SupraSolv® solvents for gas chromatography  | 6 |
| SupraSolv $^{^{(\!g)}}$ headspace solvents for headspace gas chromatography         | 8 |
| SupraSolv® hypergrade solvents for organic trace analysis1                          | 0 |
| SupraSolv® hypergrade solvents for organic trace analysis1                          | 2 |
| Optimum packaging and withdrawal systems1   | 4 |
| Ordering information 1  | 6 |
| Column Selection by GC Technique & Application1                                     | 9 |
| Reference substances for gas chromatography2  | 1 |

### Quality for the widest range of applications

Our solvents for gas chromatography offer the highest quality for the widest range of applications. The retention time window during which specified contaminants are lowest is broader than that of any competitor product.

### Quality for the best batch consistency

Every batch that leaves our premises is tested to the same high standards of quality. This is your guarantee of consistently reliable analytical results.

### **Quality for reliable analyses**

Our solvents for gas chromatography provide reliable, consistent analytical results without the need for costly purification and repeat analysis. As a result, they make your work easier, more efficient and more economical.


# SupraSolv® and SupraSolv® hypergrade solvents for all your gas chromatography needs


### FID FLAME IONIZATION DETECTOR

- BTX (Benzene, Toluene, Xylene) highly volatile aromatic hydrocarbons in sewage, ground-water, juices, canned fish etc.
- Hydrocarbon-oil index in water
- Determination of emissions in car cockpit material

#### ECD ELECTRON CAPTURE DETECTOR

- Pesticide analysis in fruits and vegetables
- Acrylamide in e.g. potato chips, crisps and crisp bread
- · Polychlorinated biphenyls (PCB) in water and sludge
- DDT (preserver and insecticide) in milk, fish, meat, fruits etc.
- Highly volatile halogenated hydrocarbons in water
- Nitrate in lettuce, radish etc.






### SupraSolv® solvents for gas chromatography

With gas chromatography, only solvents with the highest levels of purity are suitable for sample preparation tasks such as the extraction and concentration of the extracts before injection. SupraSolv® solvents are developed specially for this highly sophisticated application area.

Our comprehensive portfolio of GC solvents offers the right product for your specific application and detection method. SupraSolv® ECD and FID is specially developed and tested for ECD (Electron Capture Detector) and FID (Flame Ionization Detector). Typical applications include the determination of polychlorinated biphenyls (PCB) in water and soil or pesticides in fruits and vegetables. SupraSolv® MS is dedicated for use in gas chromatography coupled with mass spectrometric detection. This method is of increasing importance and used e.g. for the analysis of dioxins and furans (PCDD/PCDF) in food and water samples or for the determination of PAH (polycyclic aromatic hydrocarbons) in food. Both SupraSolv® qualities are carefully tested for the specific detectors and offer a clear baseline and minimal signal-to-noise ratio within a specified retention time range. Therefore SupraSolv® solvents help you achieve consistently accurate, reliable and reproducible results.

### **EPA Method 508: Determination of** chlorinated pesticides in water, standard chromatogram

Classical pesticide analysis according EPA Method 508 is employed for the qualitative and quantitative determination of pesticides in food and environmental samples. This method uses gas chromatography coupled with ECD. For sample preparation, the solvents Dichloromethane Methanol, Methyl tert-butyl ether (MTBE) or Acetone are used. Due to their particular suitability for GC-ECD, as well as their high purity and minimal interference signals within the relevant retention time range, SupraSolv® solvents will help you to achieve consistently accurate, reliable and reproducible results in pesticide analysis. Furthermore, the specified ECD retention time range of SupraSolv® ECD and FID covers all analytes of interest for this application, resulting in best application security.



- 1. Etridiazole
- 2. Chlorneh
- Propachlor
- 4. Trifluralin 5. a-HCH
- 6. Hexachlorobenzene
- 7. B-HCH
- 8. Ω-HCH 9. γ-HCH
- 10. Chlorothalonil
- 11. Heptachlor
- 12. Aldrin
- 13. DCPA
- 14. Heptachlor epoxide
- 15. γ-Chlordane

- 16. Endosulfan I
- 17. a-Chlordane
- 18. Dieldrin
- 19. 4,4'-DDE
- 20. Endrin 21. Endosulfan II
- 22. Chlorobenzilate
- 23. 4,4'-DDD
- 24. Endrin aldehyde
- 25. Endosulfan sulfate 26. 4,4'-DDT
- 27. Methoxychlor
- 28. cis-Permethrin
- 29. trans-Permethrin


#### SupraSolv® - the reliable solution

SupraSolv® has minimal interference signals in the relevant retention time window (Fig. 2). This ensures reliable, reproducible and accurate analysis results.

Thanks to outstanding batch consistency, SupraSolv® also saves you time and money by making repeat analyses a thing of the past.


#### SupraSolv® ECD specification

Outstanding analytical capabilities form the basis for providing you with comprehensive quality information - the specifications document our quality level and give you the reliability you need for your day-to-day laboratory decisions.



### Use of a non-specific solvent and competitor comparison

Both competitor chromatogram and the chromatogram of a non-GC-specific solvent (n-Hexane reagent grade) exhibit highly unstable baselines and many unidentifiable contaminant peaks. The competitor chromatogram also shows very low batch consistency. The bottom line: No clear analytical results, a risk of misinterpretation - and expensive, time-consuming repeat analyses.



[Hz] 18,000 14,000 10,000 6,000 2,000 0 2.5 7.5 20[min] 10 12.5 15 17.5 Batch chromatogram - Reference chromatogram

Figure 3B GC-ECD, batch chromatogram, n-Hexane GC grade, competitor  ${\sf A}.$ 

## SupraSolv® headspace solvents for headspace gas chromatography

Headspace gas chromatography is a precise, well-accepted method for the analysis of residual solvents in drug substances and products. It is recommended as the preferred method of analysis for this application by the European Pharmacopoeia (Chapter 2.4.24) and the United States Pharmacopoeia (Chapter 467).

The ICH (International Conference on Harmonization of Technical Requirements for Registration of Pharmaceuticals for Human Use) Guideline Q3C »Impurities: Guideline for Residual Solvents« divides all residual solvents into three classes according to their harmfulness for human health, and defines permissible maximum concentrations in actives, excipients and drug products. Both the European


### **Extract of specification**

ICH = International Conference on Harmonisation of Technical Requirements for Registration of Pharmaceuticals for Human Use


and the United States Pharmacopoeia refer to this guideline. Accurate analysis with headspace gas chromatography demands the use of very pure solvents with extremely low concentrations of the defined residual solvents. By specifying for SupraSolv® headspace the concentrations of all residual solvents of the three defined classes in the ICH guideline, we offer a precise purity window for this application — for unique, application-orientated quality. Since we also perform a headspace application test on each batch, every delivery gives you the reliability, accuracy and analytical safety you need.

**SupraSolv**<sup>®</sup> **headspace** solvents are specially designed for the analysis of residual solvents

Every residual solvent of class 1 acc. ICH  $\leq$  1 µg/g Every residual solvent of class 2 acc. ICH  $\leq$  10 µg/g Every residual solvent of class 3 acc. ICH  $\leq$  50 µg/g



**Figure 4** Chromatogram of DMSO Headspace SupraSolv® 101900 without addition compared to a chromatogram of DMSO Headspace SupraSolv® 101900 with 0.8 ppm benzene.



**Figure 5** Chromatogram of DMF Headspace SupraSolv® 100202 without addition compared to a chromatogram of DMF Headspace SupraSolv® 100202 with 0.8 ppm benzene.

according to Ph Eur and USP. We have developed them in close cooperation with an experienced headspace laboratory, and manufacture them using

special production processes. As a result, these high purity products ensure reliable, accurate analytical results.

### **Application: Quantification of residual solvents in an API**




Figure 6 Quantification of residual solvents in an AFI using Dimethyl sulfoxide (DMSO) SupraSolv $^{\otimes}$  for headspace gas chromatography (101900).

#### **Chromatographic conditions**

| Column      | fused silica capillary column, DB 1, ID 0.32 mm, film 5 $\mu m$ | length 30 m,   |
|-------------|---|----------------|
| Pressure    | 0.6 bar / 8 psi (Helium)  |                |
| Injection   | splitless, 150 °C   |                |
| Headspace c | onditions   |                |
|             | thermostating temperature                                       | 80 °C          |
|             | transfer and needle temperature                                 | 130 °C         |
|             | thermostating time  | 30 min         |
|             | pressurisation  | 1.0 min        |
|             | injection time  | 0.04 min       |
|             | withdrawal time   | 0.2 min        |
|             | High pressure   | 2 bar / 28 psi |
| Detection   | FID, 250 °C   |                |
| Temperature | 50 °C for 5 min, with 8 °C/min up t<br>hold 240 °C for 5 min    | o 240 °C,      |
| Method      | Quantification of residual solvents in                          | n an API       |

#### **Chromatographic data**


| No. | Compound                | Time [min] | Area   |
|-----|-------------------------|------------|--------|
| 1   | Methanol                | 2.0        | 12361  |
| 2   | Ethanol                 | 3.8        | 399048 |
| 3   | 2-Propanol              | 5.4        | 2368   |
| 4   | tert-Butyl methyl ether | 9.0        | 34637  |
| 5   | Diisopropyl ether       | 10.5       | 43000  |
| 6   | Tetrahydrofuran         | 11.4       | 14083  |
| 7   | Toluene                 | 11.5       | 11502  |

## SupraSolv® hypergrade solvents for organic trace analysis

SupraSolv® hypergrade solvents offer a unique solution for all applications. The specification is even broader and higher than that of SupraSolv® hypergrade solvents: the specified retention time range is larger (so even low-boiling substances can be reliably detected), and the permissible concentration of interference signals within the retention time range is also lower. We recommend SupraSolv® solvents for all areas that demand the highest levels of reliability in analytical results — for example, environmental analyses. Intensive resear ch — combined with ongoing product development — not only ensures reliability in standard applications, but also permits easier, more precise analyses in new fields, such as determining the Hydrocarbon-oil index of water and soil samples.

### Pesticide residue analysis in apple juice with GC-MS and SupraSolv® hypergrade Dichloromethane

As the consumption of fruit-based soft drinks is significant in vulnerable groups, these products should be closely monitored regarding pesticides. Classical pesticide residue analysis is still performed with GC-ECD (see page 6) and SupraSolv® solvents n-Hexane, Ethyl acetate, Dichloromethane or Acetone. An alternative and faster method using QuEChERS sample preparation uses GC-MS instead. This method reduces manual effort, improves analytical safety, and extends the range of detectable pesticides. The extraction agent with the best dissolution properties for pesticides is Dichloromethane SupraSolv® (106454).




**Figure 7** Sample chromatogram (TIC), apple juice spiked. Sample preparation via liquid-liquid-extraction with EXtrelut® NT 20. Eluting solvent: Dichloromethane SupraSolv® hypergrade (1.06454.1000).

### Hydrocarbon-oil index analysis (C10-C40) in water in accordance with DIN EN ISO 9377-2 (H53, 2001) using SupraSolv® hypergrade

The specifications of SupraSolv® hypergrade solvents have been fine-tuned to permit even the determination of the hydrocarbon-oil index in water (DIN EN ISO 9377-2 [H53, 2001]). SupraSolv® petroleum benzine is used as the extraction agent for this application. Its minimal signal-to-noise ratio allows the employment

of even steeper temperature gradients. For you, this means that your analysis times are reduced to a minimum, while the quality of the results remains unchanged. This allows your sample throughput to be greatly increased: an advantage from which you will profit daily.








## SupraSolv® hypergrade solvents for organic trace analysis

No matter which gas chromatography method you use, and regardless of whether you are analyzing soil or water samples: With SupraSolv® hypergrade solvents, you only need to use a single quality. SupraSolv® hypergrade

solvents are specified for GC-ECD and GC-FID, and also for mass spectroscopy (MS), which is rapidly growing in importance for the structure determination and quantification of sample components (Fig. 9a - 11).







|   | GC-ECD  | GC-ECD pesticide analysis  | GC-FID  | GC-MS   |
|---|---|--|---|---|
| Specifications at a glance                                | Dichloromethane to 1,2,4-Trichlorobenzene (Tetrachloromethane standard) | 1,2,4-Trichlorobenzene<br>to Decachlorobiphenyle<br>(Lindane standard) | n-Undecane to<br>n-Tetracontane<br>(n-Tetradecane standard) | n-Undecane to<br>n-Tetracontane; scan<br>range 30 - 600 amu<br>(n-Tetradecane standard) |
| SupraSolv® solvents for gas chromatography ECD and FID    | -   | max. 3 pg/mL   | max. 3 ng/mL  | -   |
| SupraSolv® solvents for gas chromatography MS             | -   | -  | -   | max. 3 ng/mL  |
| SupraSolv® hypergrade solvents for organic trace analysis | max. 1 ng/mL  | max. 2 pg/mL   | max. 2 ng/mL  | max. 2 ng/mL  |

## Optimum packaging and withdrawal systems

### Quality for high-grade packaging

Our quality standards apply to each individual package of SupraSolv® solvents. We place great value on providing a large selection of application-orientated package sizes, from 1 liter glass bottles up to 10 and 30 liter returnable barrels made of stainless steel.

#### **Glass bottles**

Optimum characteristics for handling, storage and transport. Safe footprint, low center of gravity, optimum emptying. Safety screw cap S40.

#### Stainless steel barrels

Optimum material characteristics (avoidance of interactions between solvents and packaging material). These stackable, returnable barrels ensure optimum emptying, and can be combined with a variety of withdrawal systems.

#### 10 Liter



#### 30 Liter





### Sustainable environmental protection

Stainless steel barrels are unbreakable, and help to minimize packaging waste and environmental pollution risks. By using them, suppliers and users can proactively help to protect our environment. Our stainless steel barrels are returnable, and remain the property of Merck KGaA Darmstadt, Germany throughout their life cycle. After use, empty barrels must be returned. We will then ensure they are properly cleaned, checked and refilled.

### Important safety advice

Withdrawal of flammable liquids should only be made from vessels that have been properly earthed as well as the withdrawal system itself. This can be done e.g. using our antistatic device (Cat. No. 1.07070.0001).

Our withdrawal systems have been developed and optimized for the use with Supelco® containers and solvents. We therefore disclaims any warranty or liability for the operability of its withdrawal systems in connection with containers or solvents from other manufacturers. We reserve the right to refrain from the delivery of withdrawal systems if the respective order does not indicate that each withdrawal system will be used in combination with appropriate Supelco® solvents and containers.

We inform and advise our customers to the best of our knowledge and ability but without any engagement or liability on our part. Our customers must obey all existing laws and regulations. This also applies in respect of any protected rights of third parties. Our information and advice does not eliminate the need for our customers to check, on their own responsibility, that our products are suitable for the purpose envisaged.

### Quality for reliable and safe dispensing

Our specially designed withdrawal systems for stainless steel containers safeguard the high quality of our solvents during dispensing. To prevent any detectable contamination, solvents, packaging and dispensing systems are optimally matched. Our withdrawal systems also provide the highest levels of safety when used in daily work.

### Withdrawal system for manual pressure build-up in stainless steel barrels

System compounds: exchangeable dip tube for 10 and 30 liter stainless steel barrels, clamp for outlet tube, ball valve, pump ball, three-way stopcock.

Advantages: independent from gas supply, enables simple and safe filling of smaller bottles or containers, can be used in the laboratory, central storage possible.

### Withdrawal system for inert gas pressurizing in stainless steel barrels

System compounds: 2"-thread adapter with two rapid-action connections, spiral gas feed tube for pressurizing, stainless steel-coated PTFE tube with rapid-action connector and threaded connector, self-closing stainless steel nozzle with large handle.

Advantages: flexible tubing (gas: 180 cm usable length, product: 80 cm length), enables simple and safe filling of smaller bottles or containers, central storage and supply possible.





### Ordering information

### SupraSolv® solvents for gas chromatography ECD and FID

|   | Product                   | Purity (GC)<br>min. [%] | Evap. residue<br>max. [mg/L] | Water max.<br>[%] | Color max.<br>[Hazen] | Content /<br>Packaging | Cat. No.     |
|---|---------------------------|-------------------------|------------------------------|-------------------|-----------------------|------------------------|--------------|
| Α | Acetone                   | 99.8                    | 3.0                          | 0.05              | 10                    | 1 L GL                 | 1.00012.1000 |
|   |                           |                         |                              |                   |                       | 2.5 L GL               | 1.00012.2500 |
|   |                           |                         |                              |                   |                       | 4 L GL                 | 1.00012.4000 |
|   |                           |                         |                              |                   |                       | 30 L ST                | 1.00012.9030 |
|   | Acetonitrile              | 99.8                    | 3.0                          | 0.05              | 10                    | 1 L GL                 | 1.00017.1000 |
|   |                           |                         |                              |                   |                       | 2.5 L GL               | 1.00017.2500 |
|   |                           |                         |                              |                   |                       | 4 L GL                 | 1.00017.4000 |
| В | tert-Butyl methyl ether   | 99.8                    | 3.0                          | 0.02              | 10                    | 1 L GL                 | 1.01995.1000 |
|   |                           |                         |                              |                   |                       | 2.5 L GL               | 1.01995.2500 |
| С | Chloroform, stabilized    | 99.8                    | 5.0                          | 0.01              | 10                    | 1 L GL                 | 1.02432.1000 |
|   |                           |                         |                              |                   |                       | 2.5 L GL               | 1.02432.2500 |
|   | Cyclohexane               | 99.8                    | 3.0                          | 0.01              | 10                    | 1 L GL                 | 1.02817.1000 |
|   |                           |                         |                              |                   |                       | 2.5 L GL               | 1.02817.2500 |
|   |                           |                         |                              |                   |                       | 4 L GL                 | 1.02817.4000 |
|   |                           |                         |                              |                   |                       | 10 L ST                | 1.02817.9010 |
| D | Dichloromethane,          | 99.8                    | 5.0                          | 0.01              | 10                    | 1 L GL                 | 1.06054.1000 |
|   | stabilized                |                         |                              |                   |                       | 2.5 L GL               | 1.06054.2500 |
|   |                           |                         |                              |                   |                       | 4 L GL                 | 1.06054.4000 |
|   |                           |                         |                              |                   |                       | 10 L ST                | 1.06054.9010 |
|   | Diethyl ether, stabilized | 98.0                    | 3.0                          | 0.05              | 10                    | 1 L GL                 | 1.00931.1000 |
|   |                           |                         |                              |                   |                       | 2.5 L GL               | 1.00931.2500 |
|   |                           |                         |                              |                   |                       | 4 L GL                 | 1.00931.4000 |
|   | N,N-Dimethylformamide     | 99.8                    | 3.0                          | 0.05              | 10                    | 1 L GL                 | 1.10983.1000 |
|   |                           |                         |                              |                   |                       | 2.5 L GL               | 1.10983.2500 |
| E | Ethanol                   | 99.8                    | 3.0                          | 0.01              | 10                    | 1 L GL                 | 1.02371.1000 |
|   |                           |                         |                              |                   |                       | 2.5 L GL               | 1.02371.2500 |
|   | Ethyl acetate             | 99.8                    | 3.0                          | 0.02              | 10                    | 1 L GL                 | 1.10972.1000 |
|   |                           |                         |                              |                   |                       | 2.5 L GL               | 1.10972.2500 |
|   |                           |                         |                              |                   |                       | 4 L GL                 | 1.10972.4000 |
|   |                           |                         |                              |                   |                       | 10 L ST                | 1.10972.9010 |
|   |                           |                         |                              |                   |                       | 30 L ST                | 1.10972.9030 |
| Н | n-Hexane                  | 98.0 *                  | 3.0                          | 0.01              | 10                    | 1 L GL                 | 1.04371.1000 |
|   |                           |                         |                              |                   |                       | 2.5 L GL               | 1.04371.2500 |
|   |                           |                         |                              |                   |                       | 4 L GL                 | 1.04371.4000 |
|   |                           |                         |                              |                   |                       | 10 L ST                | 1.04371.9010 |
|   |                           |                         |                              |                   |                       | 30 L ST                | 1.04371.9030 |
| I | Isohexane                 | 99.8                    | 3.0                          | 0.01              | 10                    | 2.5 L GL               | 1.04340.2500 |
|   | Isooctane                 | 99.8                    | 3.0                          | 0.01              | 10                    | 1 L GL                 | 1.15440.1000 |
|   |                           |                         |                              |                   |                       | 2.5 L GL               | 1.15440.2500 |
| М | Methanol                  | 99.8                    | 3.0                          | 0.1               | 10                    | 1 L GL                 | 1.06011.1000 |
|   |                           |                         |                              |                   |                       | 2.5 L GL               | 1.06011.2500 |
|   |                           |                         |                              |                   |                       | 4 L GL                 | 1.06011.4000 |

|   | Product           | Purity (GC)<br>min. [%] | Evap. residue<br>max. [mg/L] | Water max.<br>[%] | Color max.<br>[Hazen] | Content /<br>Packaging | Cat. No.     |
|---|-------------------|-------------------------|------------------------------|-------------------|-----------------------|------------------------|--------------|
| Р | n-Pentane         | 99.8                    | 3.0                          | 0.02              | 10                    | 1 L GL                 | 1.00882.1000 |
|   |                   |                         |                              |                   |                       | 2.5 L GL               | 1.00882.2500 |
|   | Petroleum benzine | -                       | 3.0                          | 0.01              | 10                    | 1 L GL                 | 1.01772.1000 |
|   | (40 - 60 °C)      |                         |                              |                   |                       | 2.5 L GL               | 1.01772.2500 |
|   |                   |                         |                              |                   |                       | 4 L GL                 | 1.01772.4000 |
|   |                   |                         |                              |                   |                       | 10 L ST                | 1.01772.9010 |
|   |                   |                         |                              |                   |                       | 30 L ST                | 1.01772.9030 |
|   | 2-Propanol        | 99.8                    | 3.0                          | 0.1               | 10                    | 1 L GL                 | 1.00998.1000 |
|   |                   |                         |                              |                   |                       | 2.5 L GL               | 1.00998.2500 |
| Т | Toluene           | 99.8                    | 3.0                          | 0.03              | 10                    | 1 L GL                 | 1.08389.1000 |
|   |                   |                         |                              |                   |                       | 2.5 L GL               | 1.08389.2500 |
|   |                   |                         |                              |                   |                       | 4 L GL                 | 1.08389.4000 |
|   |                   |                         |                              |                   |                       | 10 L ST                | 1.08389.9010 |
| W | Water             | -                       | 5.0                          | -                 | 10                    | 1L GL                  | 1.02699.1000 |
|   |                   |                         |                              |                   |                       | 2.5 L GL               | 1.02699.2500 |

GL = glass bottle I ST = stainless steel barrel I \* = sum of hexane isomers + methyl cyclopentane (GC)  $\geq$  99.8 % GC-ECD (retention range 1,2,4-Trichlorobenzene to Decachlorobiphenyle individual signals (Lindane standard)) :, 3 pg/mL GC-FID (retention range n-Undecane to n-Tetracontane individual signals (n-Tetradecane standard)) :, 3 ng/mL

### SupraSolv® solvents for gas chromatography MS

|   | Product                     | Purity (GC)<br>min. [%] | Evap. residue<br>max. [mg/L] | Water<br>max. [0%] | Color<br>max. [Hazen] | Content /<br>Packaging | Ord. No.     |
|---|-----------------------------|-------------------------|------------------------------|--------------------|-----------------------|------------------------|--------------|
| Α | Acetone                     | 99.8                    | 3.0                          | 0.05               | 10                    | 1 L GL                 | 1.00658.1000 |
|   |                             |                         |                              |                    |                       | 2.5 L GL               | 1.00658.2500 |
|   | Acetonitrile                | 99.8                    | 3.0                          | 0.05               | 10                    | 1 L GL                 | 1.00665.1000 |
|   |                             |                         |                              |                    |                       | 2.5 L GL               | 1.00665.2500 |
| С | Cyclohexane                 | 99.8                    | 3.0                          | 0.01               | 10                    | 1 L GL                 | 1.00667.1000 |
|   |                             |                         |                              |                    |                       | 2.5 L GL               | 1.00667.2500 |
| D | Dichloromethane, stabilized | 99.8                    | 5.0                          | 0.01               | 10                    | 1 L GL                 | 1.00668.1000 |
|   |                             |                         |                              |                    |                       | 2.5 L GL               | 1.00668.2500 |
| Е | Ethyl acetate               | 99.8                    | 3.0                          | 0.02               | 10                    | 1 L GL                 | 1.00789.1000 |
|   |                             |                         |                              |                    |                       | 2.5 L GL               | 1.00789.2500 |
| Н | n-Hexane                    | 98.0 *                  | 3.0                          | 0.01               | 10                    | 1 L GL                 | 1.00795.1000 |
|   |                             |                         |                              |                    |                       | 2.5 L GL               | 1.00795.2500 |
| М | Methanol                    | 99.8                    | 3.0                          | 0.1                | 10                    | 1 L GL                 | 1.00837.1000 |
|   |                             |                         |                              |                    |                       | 2.5 L GL               | 1.00837.2500 |
| Т | Toluene                     | 99.8                    | 3.0                          | 0.03               | 10                    | 1 L GL                 | 1.00849.1000 |
|   |                             |                         |                              |                    |                       | 2.5 L GL               | 1.00849.2500 |
| W | Water                       | -                       | 5.0                          | -                  | 10                    | 1L GL                  | 1.03702.1000 |
|   |                             |                         |                              |                    |                       | 2.5 L GL               | 1.03702.2500 |

GL = glass bottle I \* = sum of hexane isomers + methyl cyclopentane (GC)  $\geq$  99.8 %

GC-MS (retention range n-Undecane to n-Tetracontane; scanning area 30 - 600 amu individual signals (n-Tetradecane standard))  $\leq$  3 ng/mL

### SupraSolv® headspace for the analysis of residual solvents according to ICH, Ph Eur and USP

|   | Product                | Purity (GC)<br>min. [%] | Evap. residue<br>max. [mg/L] | Water<br>max. [%] | Color<br>max. [Hazen] | Content /<br>Packaging | Cat. No.     |
|---|------------------------|-------------------------|------------------------------|-------------------|-----------------------|------------------------|--------------|
|   | 1-Methyl-2-pyrrolidone | 99.8                    | -                            | 0.05              | 10                    | 0.5 L GL               | 1.02497.0500 |
|   |                        |                         |                              |                   |                       | 1 L GL                 | 1.02497.1000 |
|   |                        |                         |                              |                   |                       | 2.5 L GL               | 1.02497.2500 |
| В | Benzylalcohol          | 99.5                    | -                            | 0.1               | 10                    | 1 L GL                 | 1.02695.1000 |
|   |                        |                         |                              |                   |                       | 2.5 L GL               | 1.02695.2500 |
| D | N,N-Dimethylacetamide  | 99.8                    | 3.0                          | 0.05              | 10                    | 0.5 L GL               | 1.00399.0500 |
|   | N,N-Dimethylacetamide  | 99.8                    | 3.0                          | 0.05              | 10                    | 1 L GL                 | 1.00399.1000 |
|   | N,N-Dimethylformamide  | 99.8                    | 3.0                          | 0.05              | 10                    | 1 L GL                 | 1.00202.1000 |
|   |                        |                         |                              |                   |                       | 2.5 L GL               | 1.00202.2500 |
|   | Dimethyl sulfoxide     | 99.8                    | 3.0                          | 0.05              | 10                    | 1 L GL                 | 1.01900.1000 |
|   |                        |                         |                              |                   |                       | 2.5 L GL               | 1.01900.2500 |
| W | Water                  | -                       | 5.0                          | -                 | -                     | 1 L GL                 | 1.00577.1000 |
|   |                        |                         |                              |                   |                       | 2.5 L GL               | 1.00577.2500 |

GL = glass bottle

Every residual solvent of class 1 acc. CH  $\leq$  1 µg/g Every residual solvent of class 2 acc. CH  $\leq$  10 µg/g Every residual solvent of class 3 acc. CH  $\leq$  50 µg/g

### SupraSolv® hypergrade solvents for organic trace analysis

|   | Product           | Purity (GC)<br>min. [%] | Evap. residue<br>max. [mg/l] | Water max.<br>[%] | Color<br>max. [Hazen] | Content /<br>Packaging | Cat. No.     |
|---|-------------------|-------------------------|------------------------------|-------------------|-----------------------|------------------------|--------------|
| D | Dichloromethane   | 99.9                    | 3.0                          | 0.005             | 10                    | 1 L GL                 | 1.06454.1000 |
| Н | n-Hexane          | 99.0*                   | 3.0                          | 0.005             | 10                    | 1 L GL                 | 1.04369.1000 |
|   |                   |                         |                              |                   |                       | 2.5 L GL               | 1.04369.2500 |
| Р | n-Pentane         | 99.9                    | 3.0                          | 0.01              | 10                    | 1 L GL                 | 1.07288.1000 |
|   |                   |                         |                              |                   |                       | 2.5 L GL               | 1.07288.2500 |
|   | Petroleum benzine | -                       | 3.0                          | 0.005             | 10                    | 1 L GL                 | 1.16740.1000 |
|   | (40 - 60°C)       |                         |                              |                   |                       | 2.5 L GL               | 1.16740.2500 |

GL = glass bottle

GC-ECD (retention range Dichloromethane to 1,2,4-Trichlorobenzene individual signals (Tetrachloromethane standard))  $\leq 1$  ng/mL

GC-ECD (retention range 1,2,4-Trichlorobenzene to Decachlorobiphenyle individual signals (Lindane standard))  $\leq 2 \text{ pg/mL}$ 

GC-FID (retention range n-Undecane to n-Tetracontane individual signals (n-Tetradecane standard))  $\leq 2 \text{ ng/mL}$ 

GC-MS (retention range n-Undecane to n-Tetracontane; scanning area 30 - 600 amu individual signals (n-Tetradecane standard) ≤ 2 ng/mL



<sup>\*</sup> Sum of hexane isomers + methylcyclopentane (GC) > 99.9 %

## Column Selection by GC Technique & Application

We have developed an extensive line of special purpose columns designed for industry specific applications. These columns are manufactured to deliver high resolution, great analyte response, low bleed, and long column life; allowing you to achieve the analytical performance you require. Look at our easy-to-read phase selection charts highlight choices for applications that are independent of any industry. Simply locate the application to identify a recommended column phase.

The stationary phase also dictates the minimum and maximum temperatures at which a column can be used. Therefore, it is critical to ensure the selected stationary phase can withstand the temperature requirements of the GC method.

### **Fast GC Applications**

|                                    | SLB®-1ms** | Equity®-1** | SLB®-5ms** | Equity®-5** | SPB®-624 | %*®TOOON | Equity®-1701** | Omegawax®** | SUPELCOWAX® 10** | SLB®-1L59 | SP®-2560 | SLB®-1L82 | SLB®-1L100 | SLB®-IL111i |
|------------------------------------|------------|-------------|------------|-------------|----------|----------|----------------|-------------|------------------|-----------|----------|-----------|------------|-------------|
| Volatiles                          |            |             |            |             | •        | •        |                |             |                  |           |          |           |            |             |
| Semivolatiles                      |            |             | •          |             |          |          |                |             |                  |           |          |           |            |             |
| Pesticides                         |            |             | •          | •           |          |          | •              |             |                  |           |          |           |            |             |
| PCBs                               |            |             | •          | •           |          |          | •              |             |                  | •         |          | •         |            | •           |
| Fuels by Pattern Recognition       | •          | •           | •          |             |          |          |                |             |                  |           |          |           |            |             |
| Aromatics                          |            |             |            |             |          |          |                |             | •                | •         |          |           | •          | •           |
| Sulfur Compounds                   |            |             |            |             |          |          |                |             |                  | •         |          |           |            |             |
| Biodiesel: FAME Profile            |            |             |            |             |          |          |                | •           |                  |           |          |           |            | •           |
| Solvents                           |            |             | •          |             |          |          |                |             | •                |           |          |           |            |             |
| Aromatic Amines (Anilines)         |            |             |            |             |          |          |                |             |                  | •         |          |           |            |             |
| Omega 3 and Omega 6 FAMEs          |            |             |            |             |          |          |                | •           |                  | •         |          |           |            |             |
| cis/trans FAME Isomers             |            |             |            |             |          |          |                |             |                  |           | •        |           |            | •           |
| Nitrosamines                       |            |             | •          |             |          |          |                |             |                  | •         |          |           |            |             |
| Essential Oils                     | •          | •           | •          |             |          |          |                |             | •                |           |          |           |            |             |
| Drugs of Abuse                     |            |             | •          |             |          |          |                |             |                  |           |          |           |            |             |
| General Purpose Non-Polar          | •          | •           |            | •           |          |          |                |             |                  |           |          |           |            |             |
| General Purpose Intermediate Polar |            |             |            |             |          |          | •              |             |                  |           |          |           |            |             |
| General Purpose Polar              |            |             |            |             |          |          |                |             | •                |           |          |           |            |             |

<sup>\*\*</sup> Indicated columns are available in Intuvo format

### **GCxGC Applications**

|                                 | SLB®-5ms** | Equity®-5** | SPB®-5** | PTA-5 | SAC-5 | SLB®-1ms** | SPB®-1** | SUPELCOWAX® 10** | SLB®-IL60i** | SLB®-1L59 | SLB®-1L61 | SLB®-1L111i | SP®-2380** | SLB®-1L76i | SLB®-IL82 | SP®-2331** |
|---------------------------------|------------|-------------|----------|-------|-------|------------|----------|------------------|--------------|-----------|-----------|-------------|------------|------------|-----------|------------|
| Non-Polar Primary (1°) Column   | •          | •           | •        | •     | •     | •          | •        |                  |              |           |           |             |            |            |           |            |
| Polar Secondary (2°) Column     |            |             |          |       |       |            |          | •                |              | •         | •         | •           |            | •          | •         |            |
| Polar Primary (1°) Column       |            |             |          |       |       |            |          | •                | •            | •         | •         | •           | •          | •          | •         | •          |
| Non-Polar Secondary (2°) Column | •          | •           |          |       |       | •          |          |                  |              |           |           |             |            |            |           |            |

<sup>\*\*</sup> Indicated columns are available in Intuvo format

### **Chiral Applications**

|   | Astec <sup>®</sup> CHIRALDEX <sup>®</sup> TA | Astec® CHIRALDEX® PN | Astec <sup>®</sup> CHIRALDEX <sup>®</sup> DP | Astec <sup>®</sup> CHIRALDEX <sup>®</sup> BP | Astec <sup>®</sup> CHIRALDEX <sup>®</sup> DM | Supelco® DEX 325 | Supelco® DEX 225 | Astec <sup>®</sup> CHIRALDEX <sup>®</sup> PM | Supelco® DEX 110 | Supelco® DEX 120** | Astec <sup>®</sup> CHIRALDEX <sup>®</sup> DA | Astec <sup>®</sup> CHIRALDEX <sup>®</sup> PH | a-Cyclodextrins | β-Cyclodextrins | γ-Cyclodextrins |
|---|--|----------------------|--|--|--|------------------|------------------|--|------------------|--------------------|--|--|-----------------|-----------------|-----------------|
| Oxygen containing analytes in the form of alcohols, ketones, acids, aldehydes, and lactones; halogenated compounds                                  | •  |                      |  |  |  |                  |                  |  |                  |                    |  |  |                 |                 |                 |
| Lactones and aromatic amines; epoxides; styrene oxide   |  | •                    |  |  |  |                  |                  |  |                  |                    |  |  |                 |                 |                 |
| Aliphatic and aromatic amines;<br>aliphatic and some aromatic esters;<br>polar racemates  |  |                      | •  |  |  |                  |                  |  |                  |                    |  |  |                 |                 |                 |
| Amino acids; amines; furans   |  |                      |  | •  |  |                  |                  |  |                  |                    |  |  |                 |                 |                 |
| Aliphatic, olefinic, and aromatic enantiomers   |  |                      |  |  | •  | •                | •                |  |                  |                    |  | •  |                 |                 |                 |
| Terpenes and tertiary amines  |  |                      |  |  |  |                  |                  | •  | •                | •                  |  |  |                 |                 |                 |
| Heterocyclic amines   |  |                      |  |  |  |                  |                  |  |                  |                    | •  |  |                 |                 |                 |
| Xylenes; menthols; cresols;<br>substituted phenols; substituted<br>benzenes; epoxide enantiomers  |  |                      |  |  |  |                  |                  |  |                  |                    |  |  | •               |                 |                 |
| Acids; alcohols; amines; diols;<br>esters; ethers; halohydrocarbons;<br>hydrocarbons; ketones; positional<br>isomers; silanes; terpenes; terpineols |  |                      |  |  |  |                  |                  |  |                  |                    |  |  |                 | •               |                 |
| a-BHC; carvone; carboxylic acids; methamphetamine   |  |                      |  |  |  |                  |                  |  |                  |                    |  |  |                 |                 | •               |

 $<sup>\</sup>begin{tabular}{ll} ** Indicated columns are available in Intuvo format \\ \end{tabular}$ 

### **General Purpose (non-MS) Applications**

|                           | Equity®-1** | SPB®-1** | Equity®-5** | SPB®-5** | SPB®-20** | Equity®-1701** | SPB®-35 | SPB®-50 | SPB®-225 | PAG | SUPELCOWAX® 10** | SP®-2330 | SP®-2380** | SP®-2340** |
|---------------------------|-------------|----------|-------------|----------|-----------|----------------|---------|---------|----------|-----|------------------|----------|------------|------------|
| Non-Polar Column          | •           | •        | •           | •        |           |                |         |         |          |     |                  |          |            |            |
| Intermediate Polar Column |             |          |             |          | •         | •              | •       | •       |          |     |                  |          |            |            |
| Polar Column              |             |          |             |          |           |                |         |         | •        | •   | •                |          |            |            |
| Highly Polar Column       |             |          |             |          |           |                |         |         |          |     |                  | •        | •          | •          |

<sup>\*\*</sup> Indicated columns are available in Intuvo format

## Reference substances for gas chromatography

Most of the high-purity products in our »reference substances for GC« range are completely synthetic in origin, which means they are largely free from isomers that are difficult to separate by GC. Their assay is generally greater than 90%, and is usually over 99.5 or 99.7%. Every pack includes a gas chromatogram under the appropriate test conditions. These reference substances can be used when identifying unknown compounds in a gas chromatogram, as standards in quantitative GC analysis, or in the characterization of GC column properties. Reference substances belonging to the hydrocarbon group are packed in pierceable ampoules; fatty acid methyl esters and other reference substances come in screw-capped glass vials.

### **Ordering information**

#### n-Hydrocarbons

| Product                 | Chain length | Assay [%]                    | Content /<br>Packaging | Cat. No. |
|-------------------------|--------------|------------------------------|------------------------|----------|
| Pentane                 | C5           | ≥99.5%                       | 10 mL, 50 mL           | 76870    |
|                         |              | ≥99.7%                       | 5 mL                   | 1.09719  |
| Hexane                  | C6           | ≥99.7%                       | 10 mL, 50 mL           | 52750    |
|                         |              | ≥99.7%                       | 5 mL                   | 1.09687  |
| Heptane                 | C7           | ≥99.8%                       | 5 mL, 50 mL            | 51730    |
|                         |              | ≥99.5%                       | 5 mL                   | 1.09686  |
| Octane                  | C8           | ≥99.7%                       | 5 mL, 50 mL            | 74820    |
|                         |              | ≥99.0%                       | 5 mL                   | 1.09716  |
| Nonane                  | C9           | ≥99.8%                       | 50 mL                  | 74250    |
| Decane                  | C10          | ≥99.8%                       | 5 mL, 25 mL            | 30540    |
|                         |              | ≥99.5%                       | 5 mL                   | 1.09603  |
| Undecane                | C11          | ≥99.8%                       | 5 mL, 25 mL            | 94000    |
|                         |              | ≥99.5%                       | 5 mL                   | 1.09794  |
| Dodecane                | C12          | ≥99.8%                       | 5 mL, 25 mL            | 44010    |
|                         |              | ≥98.5%                       | 5 mL                   | 1.09658  |
| Dodecane, TraceCERT®    |              | certified content on the CoA | 100 mg                 | 92064    |
| Tridecane               | C13          | ≥99.5%                       | 5 mL, 50 mL            | 91490    |
|                         |              | ≥99.3%                       | 5 mL                   | 1.09609  |
| Tetradecane             | C14          | ≥99.5%                       | 5 mL, 25 mL            | 87139    |
|                         |              | ≥99.0%                       | 5 mL                   | 1.09658  |
| Tetradecane, TraceCERT® |              | certified content on the CoA | 100 mg                 | 04003    |
| Pentadecane             | C15          | ≥99.8%                       | 5 mL                   | 1.09658  |
| Hexadecane              | C16          | ≥99.8%                       | 5 mL, 25 mL            | 52209    |
|                         |              | ≥99.5%                       | 5 mL                   | 1.09605  |
| Heptadecane             | C17          | ≥99.5%                       | 1 mL, 5 mL             | 51578    |
|                         |              | ≥99.3%                       | 5 mL                   | 1.09604  |
| Octadecane              | C18          | ≥98.5%                       | 1 g, 5 g               | 74691    |
| Nonadecane              | C19          | ≥99.5%                       | 1 g, 5 g               | 74158    |
| Eicosane                | C20          | ≥99.5%                       | 1 g, 5g                | 44818    |
| Heneicosane             | C21          | ≥99.5%                       | 1 g, 5 g               | 51523    |
| Docosane                | C22          | ≥99.5%                       | 1 g, 5 g               | 43942    |
| Tricosane               | C23          | ≥99.5%                       | 1 g, 5 g               | 91447    |
| Tetracosane             | C24          | ≥99.5%                       | 1 g, 5 g               | 87089    |
| Pentacosane             | C25          | ≥99.5%                       | 250 mg, 1 g            | 76493    |
| Hexacosane              | C26          | ≥98.0%                       | 250 mg, 1 g            | 52183    |
| Heptacosane             | C27          | ≥98.0%                       | 250 mg, 1 g            | 51559    |
| Octacosane              | C28          | ≥98.0%                       | 250 mg, 1 g            | 74684    |
| Nonacosane              | C29          | ≥98.0%                       | 250 mg, 1 g            | 74156    |
| Triacontane             | C30          | ≥98.0%                       | 250 mg, 1 g            | 90270    |
| Hentriacontane          | C31          | ≥98.0%                       | 250 mg, 1 g            | 51529    |

| Product           | Chain length | Assay [%] | Content /<br>Packaging | Cat. No. |
|-------------------|--------------|-----------|------------------------|----------|
| Dotriacontane     | C32          | ≥98.0%    | 100 mg                 | 44253    |
| Tritriacontane    | C33          | ≥98.0%    | 250 mg                 | 93435    |
| Tetratriacontane  | C34          | ≥98.0%    | 250 mg, 1g             | 88152    |
| Pentatriacontane  | C35          | ≥98.0%    | 250 mg                 | 76968    |
| Hexatriacontane   | C36          | ≥98.0%    | 1 g                    | 52919    |
| Heptatriacontane  | C37          | ≥97.0%    | 100 mg, 1 g            | 51848    |
| Tetracontane      | C40          | ≥98.0%    | 250 mg                 | 87086    |
| Dotetracosane     | C42          | ≥98.0%    | 100 mg                 | 44250    |
| Tetratetracontane | C44          | ≥98.0%    | 250 mg                 | 88144    |
| Hexatetracontane  | C46          | ≥98.0%    | 250 mg                 | 52913    |
| Octatetracosane   | C48          | ≥98.0%    | 100 mg                 | 74892    |
| Tetrapentacontane | C54          | ≥98.0%    | 500 mg                 | 87992    |

### **Aromatic, unsaturated & branched Hydrocarbons**

| Product                | Assay  | Packaging          | Cat. No. |
|------------------------|--------|--------------------|----------|
| Benzene                | ≥99.9% | 5 mL               | 1.09646  |
|                        | ≥99.9% | 5 mL, 50 mL        | 12540    |
| Toluene                | ≥99.7% | 5 mL               | 1.09768  |
|                        | ≥99.9% | 5 mL, 25 mL        | 89680    |
| o-Xylene               | ≥99.0% | 5 mL               | 1.09798  |
|                        | ≥99.0% | 5 mL, 50 mL        | 95660    |
| m-Xylene               | ≥99.3% | 5 mL               | 1.09797  |
|                        | ≥99.5% | 5 mL, 50 mL        | 95670    |
| p-Xylene               | ≥99.3% | 5 mL               | 1.09799  |
|                        | ≥99.5% | 5 mL, 50 mL        | 95680    |
| 2-Methylbutane         | ≥99.7% | 5 mL, 10 mL        | 59060    |
| 2-Methylpentane        | ≥99.5% | 5 mL, 50 mL        | 68310    |
| 2,2-Dimethylbutane     | ≥98.0% | 5 mL, 10 mL        | 39730    |
| 2,3-Dimethylbutane     | ≥99.5% | 5 mL, 10 mL        | 39760    |
| 3,4-Dimethylhexane     | ≥99.0% | 1 mL               | 40512    |
| 2,2,4-Trimethylpentane | ≥99.7% | 5 mL, 10 mL, 50 mL | 59030    |
| 1-Hexene               | ≥99.8% | 5 mL, 10 mL, 50 mL | 52930    |
| 1-Nonene               | ≥99.5% | 5 mL               | 74323    |
| 1-Decene               | ≥99.0% | 5 mL, 50 mL        | 30649    |
| 1-Dodecene             | ≥99.5% | 5 mL               | 44146    |
| 1-Tetradecene          | ≥99.8% | 5 mL, 25 mL        | 87187    |
| 1-Octadecene           | ≥99.5% | 1 mL, 5 mL         | 74738    |
| 4-Methyl-1-pentene     | ≥99.5% | 5 mL               | 68510    |
| 1,7-Hexadecadiene      | ≥90%   | 25 mg              | 52206    |

### **Hydrocarbon mixes and kits**

| Product Description   |   | Pck size   | Cat. No. |
|---|---|------------|----------|
| Hydrocarbon Test mix, TraceCERT®                              | C12, C13, C14, C15, C16, C17: varied concentrations in CHCl3  | 2 mL       | 48244    |
| n-Paraffin mix C5, C6, C7, C8                                 | varied concentrations, neat   | 500 mg     | 47100    |
| n-Paraffin Mix C7, C8, C9, C10                                | varied concentrations, neat   | 500 mg     | 47101    |
| n-Paraffin Mix C10, C12, C14, C16                             | varied concentrations, neat   | 500 mg     | 47102    |
| n-Paraffin Mix C18, C20, C22, C24, TraceCERT®                 | 2% (w/w) each component in octane   | 5 mL       | 47108    |
| Alkan standard solution C8-C20                                | ~40mg/L each, in hexane   | 1 mL, 5 mL | 04070    |
| Alkane standard solution C21-C40                              | 40mg/L each, in toluene   | 1 mL, 5 mL | 04071    |
| Alphatic Hydrocarbons Kit                                     | 34 individual ampules of neats plus 6 ampules of solutions  | 1 each     | 44575-U  |
| Florisil Applicability Test acc. to DIN EN ISO 9377-2/53, kit | Kit contains: Stearyl stearate Minal oil standard mixture type A and B (~10mg/mL in n-heptane) Extraction medium stock solution | 1 each     | 52462    |

| Product   | Description  | Pck size    | Cat. No. |
|---|--|-------------|----------|
| Extraction Medium Stock solution                            | Acc. to DIN EN ISO 9377-2                            | 100 mL      | 49574    |
| Mineral oil standard mixture Type A and B for EN ISO 9377-2 | 71 57 1  |             | 18602    |
| QC Standard solution for EN ISO 9377-2                      | Mineral oil type A and B ~5mg/L each in n-heptane    | 1 mL        | 51706    |
| Alkane standard mixture for performance tests of GC systems | C10 - C40 (all even), 50 mg/L each in n-heptane      | 2 mL        | 68281    |
| Hydrocarbon kit acc. to DIN EN 14039/ISO 16703              |  | 1 kit       | 56681    |
| Alkane Standard Mixture for performance tests of GC-systems | C10, C20 – C40 (all even), 50 mg/L each in n-heptane | 2 mL, 10 mL | 94234    |

### Fatty acid methyl esters (FAME)

| Designation                                | Empirical formula | Assay [%]                    | Content / Packaging | Cat. No. |
|--|-------------------|------------------------------|---------------------|----------|
| Methyl valerate                            | C5                | ≥99.8%                       | 1 mL, 5 mL          | 94560    |
| Methyl hexanoate                           | C6                | ≥99.8%                       | 1 mL, 5 mL          | 21599    |
| Methyl hexanoate, TraceCERT®               | C6                | certified content on the CoA | 50 mg               | 94776    |
| Methyl heptanoate                          | C7                | ≥99.8%                       | 1 mL, 5 mL          | 75218    |
| Methyl octanoate                           | C8                | ≥99.8%                       | 5 mL                | 21719    |
| Methyl octanoate, TraceCERT®               | C8                | certified content on the CoA | 50 mg               | 06934    |
| Methyl octanoate                           | C8                | ≥99.5%                       | 5 mL                | 1.09633  |
| Methyl nonanoate                           | C9                | ≥99.8%                       | 1 mL, 5 mL          | 76368    |
| Methyl decanoate                           | C10               | ≥99.5%                       | 5 mL                | 1.09637  |
| Methyl undecanoate                         | C11               | ≥99.0%                       | 1 mL, 5 mL          | 94118    |
| Methyl undecanoate, TraceCERT®             | C11               | certified content on the CoA | 100 mg              | 47147    |
| Methyl laurate                             | C12               | ≥99.5%                       | 5 mL                | 61689    |
| Methyl laurate                             | C12               | ≥99.0%                       | 5 mL                | 1.09693  |
| Methyl tridecanoate                        | C13               | ≥99.5%                       | 5 mL                | 91558    |
| Methyl myristate                           | C14               | ≥99.5%                       | 1 mL, 5 mL          | 70129    |
| Methyl myristate, TraceCERT®               | C14               | certified content on the CoA | 100 mg              | 55791    |
| Methyl myristate                           | C14               | ≥99.5%                       | 5 mL                | 1.09736  |
| Methyl pentadecanoate                      | C15               | ≥99.5%                       | 1 mL, 5 mL, 25 mL   | 76560    |
| Methyl palmitate                           | C16               | ≥99.0%                       | 1 g, 5 g            | 76159    |
| Methyl heptadecanoate                      | C17               | ≥99.0%                       | 1 g, 5 g            | 51633    |
| Methyl heptadecanoate, TraceCERT®          | C17               | certified content on the CoA | 100 mg              | 90606    |
| Methyl margarate                           | C17               | ≥99.0%                       | 5 mL                | 1.09754  |
| Methyl stearate                            | C18               | ≥99.5%                       | 1 g, 5 g            | 85769    |
| Methyl stearate, TraceCERT®                | C18               | certified content on the CoA | 100 mg              | 75533    |
| Methyl stearate                            | C18               | ≥99.0%                       | 5 mL                | 1.09602  |
| Methyl oleate (Methyl cis-9-octadecenoate) | C18-cis 9         | ≥99.0%                       | 1 mL, 5 mL          | 75160    |
| Methyl oleate (Methyl cis-9-octadecenoate) | C18-cis 9         | ≥96.0%                       | 5 mL                | 1.09743  |
| Methyl cis-11-octadecenoate                | C18-cis 11        | ≥99.0%                       | 100 mg              | 17264    |
| Methyl trans-vaccenate                     | C18-trans11       | ≥95.0%                       | 25 mg               | 07586    |
| Methyl nonadecanoate                       | C19               | ≥98.0%                       | 1 g, 5 g            | 74208    |
| Methyl arachidate                          | C20               | ≥99.0%                       | 1 g, 5 g            | 10941    |
| Methyl cis-11-eicosenoate                  | C20-cis 11        | ≥98.5%                       | 100 mg              | 17263    |
| Methyl cis,cis-11,14-eicosadienoate        | C20-cis,cis 11,14 | ≥98.0%                       | 100 mg              | 17272    |
| Methyl heneiscosanoate                     | C21               | ≥98.5%                       | 1 g                 | 51535    |
| Methyl behenate                            | C22               | ≥99.0%                       | 1 g, 5 g            | 11940    |
| Methyl tricosanoate                        | C23               | ≥99.0%                       | 250 mg, 1 g         | 91478    |
| Methyl tetracosanoate                      | C24               | ≥99.0%                       | 250 mg, 1 g         | 87115    |
| Methyl cis-15-tetracosenoate               | C24-cis 15        | ≥99.0%                       | 100 mg              | 17265    |
| Methyl pentacosanoate                      | C25               | ≥98.0%                       | 100 mg              | 76497    |
| Methyl hexacosanoate                       | C26               | ≥99.0%                       | 100 mg              | 52203    |
| Methyl octacosanoate                       | C28               | ≥98.0%                       | 250 mg              | 74701    |
|  |                   |                              |                     |          |

All *Trace*CERT® branded products are certified reference materials (CRMs), produced and certified in accordance with ISO/IEC 17025 and ISO 17034. The CRM is traceable to primary material from an NMI, e.g. NIST or MNIJ.


Certified content by quantitative NMR incl. uncertainty and expiry date are given on the certificate.

### Discover the GC workflow solutions

### PROFICIENCY TESTING

### SAMPLE PREPARATION AND COLLECTION





### Proficiency testing

- Supelco® Proficiency Testing solutions represents over 20 years of experience and expertise in providing PT studies worldwide.
- Our PT samples are manufactured to CRM grade. More than 20 000 PT samples a year are sent to over 2 500 labs worldwide.
- Our quality and services allow you to work smarter, enabling us to live in a safer and healthier world.



### Sample Preparation and Collection

- SPME Solid Phase Micro Extraction (Fibers, Holders, Accessories)
- SPE (Tubes, Manifolds, Accessories, QuEChERS, Bulk Adsorbents)
- Purge & Traps (Traps & Purge Vessels)
- Solvents (Suprasolv® FID/ECD, MS, HS)
- Milli-Q® water purification systems and VOC-Pak® Polisher
- Derivatisation Reagents

#### **Air Sampling**

- Thermal Desorption (Adsorbent Tubes & Accessories)
- Solvent Desorption (Adsorbent Tubes & Accessories)
- Whole Air (Sampling Bags & Bulbs, Sampler)





SPE Products

Solvents










SPME Products

Derivatisation Reagents

### **GC ANALYSIS**

### **DETECTION AND CALIBRATION**





### GC Analysis

- GC Columns for GC, GC/MS, Fast GC, GCxGC
- General Purpose Columns (e.g. SLB®-1ms, SLB®-5ms, SLB®-35ms, Equity®-1701, SUPELCOWAX®, Nukol™/FFAP)
- Special Application Columns (e.g. FAME, PAH, PCB, Dioxins, VOC, PLOT, Chiral)
- Ionic Liquid Columns (SLB®-IL i-Series, Watercol™)
- Packed Columns (SS & Glass)



**GC Columns** (Packed)



GC Columns (Capillary)

#### **Accessories**

- Liner & Septa
- Fittings, Ferrules & Column Connectors
- Flow Measurement
- · Vials & Syringes
- Gas Management & Purification
  - Gas Generators
  - Gas Purifiers
  - Plumbing (Valves, Fittings & Tubing)





**Gas Generators** 



Liner

Septa

### **Detection and Calibration**

- Certified Reference Materials (CRMs) & Reference Standards
- Neats & Solutions (for almost every application area)
- Matrix Standards
- Pharmacopoeia & Metrological Institute Standards



**Certified Reference Materials** & Reference Standards



| Notes |  |
|-------|--|
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |
|       |  |

| Notes |  |  |
|-------|--|--|
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |
|       |  |  |

### Supelco<sub>®</sub>

**Analytical Products** 

MilliporeSigma 400 Summit Drive Burlington, MA 01803

SigmaAldrich.com/GC

#### To place an order or receive technical assistance

Order/Customer Service: SigmaAldrich.com/order Technical Service: SigmaAldrich.com/techservice

Safety-related Information: SigmaAldrich.com/safetycenter

SigmaAldrich.com

