Millipore[®]

Preparation, Separation, Filtration & Testing Products

Now is the TIME

Shift from animal-based testing to *in vitro* testing.

EMDMillipore.com/PyroMAT

(11

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.



che Move

OMAT

Millipore

PyroMAT

KI


Now is the time! Meet regulation requirements with the PyroMAT[™] Monocyte Activation Test (MAT)

Used to detect a broad range of pyrogens in parenteral products such as pharmaceuticals, biopharmaceuticals or medical devices, the MAT gives an *in vitro* alternative to conventional animal testing in accordance with regulatory guidelines.

With **our new Mono-Mac-6 cell-line based MAT,** we provide a robust, reliable and standardized method for pyrogen detection.

PyroMAT[™]: the only cell-line based pyrogen test

Mimics the human immune reaction



The MAT kit is based on the Mono-Mac-6 (MM6) cell line and uses interleukin-6 (IL-6) as a read out.

Once these human monocytes come into contact with pyrogens from a contaminated sample, they produce cytokines such as interleukin-6 (IL-6), which can be detected with an enzyme-linked immunosorbent assay (ELISA).



With a Limit Of Detection (LOD) of 0.05 EU*/mL, PyroMAT[™] System is the only ready-to-use cell line based MAT kit on the market. A new robust and easy to perform pyrogen test.

*Endotoxin Units



Advantages of the MAT

Detection of a broad range of pyrogens: patient safety is ensured if the full range of pyrogens are tested. Like the Rabbit Pyrogen Test (RPT), MAT is effective for detection of both endotoxins and NEPs.

Extension of the range of products that can be tested: the most frequently applied methods, RPT and Bacterial Endotoxin Test (BET), are both limited in the product types they are able to test. The MAT offers more flexibility regarding its applications.

In vitro assay that mimics the human immune reaction: for a robust predictive model that reduces animal consumption.

Compliance with international regulations and guidelines: in line with ethical trends of industry and regulatory authorities to decrease the use of animal based testing.

Cell Line Benefits

Standardized reactivity and high sensitivity (LOD 0.05 EU/mL) Convenience of a ready-to-use cell line that prevents customers from laborious lab work and avoids the need for a cell culture laboratory.

Qualified cells: in addition to being cited in the international validation of MAT, Mono-Mac-6 cells are qualified for the detection of all types of pyrogens by testing the expression of all surface Toll-Like Receptors (TLRs).

PyroMAT™ System

All the advantages of the monocyte activation test combined with the benefits of using a cell line.

Ordering Information			
	Product Description	Store at	Catalog Number
PyroMAT™ kit	The kit contains reagents for the incubation with the cells and the interleukin IL-6 ELISA. Note: the PyroMAT™ cells are not included and should be ordered separately.	2-8 °C	Pyr0MATkit
PyroMAT™ cells	Contains the cryopreserved Mono-Mac-6 cells for one test (2 \times 1 mL)	-80 °C or lower	Pyr0MATcells
Reference standard endotoxin	International Reference Standard Endotoxin (RSE)—positive control of endotoxin, TLR-4 agonist	Lyophilized store at -20 °C, reconstituted at -40 °C	1.44161.0001
Non-Endotoxin Pyrogen control HKSA	Heat Killed Staphyloccocus Aureus (HKSA), a Gram-positive extra-cellular growing bacterium, TLR-2 agonist	Lyophilized store at 2-8 °C, reconstituted at -20 °C	MATHKSA
Non-Endotoxin Pyrogen control Flagellin	Flagellin, TLR-5 agonist	-20 °C	MATFLAGELLIN
Interleukin-6 control	For control reactions in the Interleukin IL-6 ELISA	2-8 °C	Pyr0MATIL6
Software protocols	To analyze results according to European Pharmacopoeia (Methods A, B and C). These analysis tools can be downloaded on EMDMillipore.com/PyroMAT		

All the lab material needed to perform the test is described in the user guide which is available on our website.



We provide additional services.

Learn more at EMDMillipore.com/ BioMonitoringServices

REQUEST INFORMATION

MilliporeSigma 400 Summit Drive Burlington, MA 01803



To Place an Order or Receive Technical Assistance

In the U.S. and Canada, call toll-free 1-800-521-8956 For other countries across Europe and the world, please visit: **EMDMillipore.com/offices** For Technical Service, please visit: **EMDMillipore.com/techservicee**

MilliporeSigma, the vibrant M, Millipore and PyroMAT are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources. © 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. PB1102ENUS 2017-06614 04/2018

VIII