Y-27632: Selective inhibitor of the Rho-associated protein kinase p160ROCK

The small GTPase Rho and its downstream effectors, the Rho-associated coiled-coil forming protein serine/threonine kinase (ROCK) family, have been implicated in various cellular functions including actin cytoskeleton organization, cell adhesion, cell motility, vascular and smooth muscle contraction and cytokinesis [1,2]. It has been suggested that Rho/ROCK inhibition may be a novel therapeutic target in the treatment of various cardiovascular diseases, such as hypertension and arteriosclerosis [2], as well as bronchial asthma [4], cancer [5] and Alzheimer's disease [6]. In recent animal experiments, Rho-kinase inhibitors were reported to suppress coronary artery spasm [3,7]. In addition, they have been shown to inhibit the development of coronary arteriosclerotic lesions and promote the regression of coronary vascular lesions in vivo [3,7]. Recent clinical studies have demonstrated the inhibitory effect of the Rho-kinase inhibitor HA-1077 (Prod. No. H-139) on coronary artery spasm in patients with vasospastic angina, and on exercise-induced myocardial ischemia in patients with stable effort angina [3,8].

Sigma-RBI is pleased to offer one of the best characterized ROCK inhibitors, Y-27632 (Prod. No. Y 0503), which selectively targets p160ROCK from the family of Rho-associated protein kinases [1]. Y-27632 inhibits human p160ROCK displaying a K_i value of 140 μM vs rat brain protein kinase C (26 mM) (PLC, Prod. No. P 7956), bovine cAMP-dependent protein kinase (25 mM) (Prod. No. C 8482), and chicken myosin light chain kinase (>250 mM) and exhibits better potency than the Rho kinase inhibitor HA-1077 and the PKA/PKC inhibitor H-7 (Prod. Nos. I 6891, I 7016) [1,9]. Y-27632 has also been shown to induce relaxation of rabbit aortic smooth muscle that has been contracted using phenylephrine (Prod. No. P 6126) displaying an IC₅₀ value of 700 μM [1]. In addition, it inhibits the contraction of pig coronary artery strips

and guinea-pig trachea induced by several agonists, including **histamine** (Prod. No. **H 7125**), **acetylcholine** (Prod. No. **A 5626**), **serotonin** (Prod. No. **H 9523**), **endothelin** (Prod. No. **E 7764**) and the thromboxane agonist **U-46619** (Prod. No. **D 8174**), exhibiting IC₅₀ values in the range of 0.3 - 1 mM [1]. In various hypertensive rat models, Y-27632 significantly dose-dependently decreased blood pressure (30 mg/kg p.o.) while exhibiting little or no effect on blood pressure in control rats [1].

These results suggest that Y-27632 should serve as a useful tool for investigating the Rho/ROCK pathway *in vivo* and for studying its role in normal physiology and disease.

Sold under license from Mitsubishi Pharma Corporation.

References

- 1. Uehata, M., et al., Nature, 389, 990-994 (1997)
- 2. Takahara, A., et al., Eur. J. Pharmacol., 460, 51-57 (2003).
- 3. Shimokawa, H.J., Cardiovasc. Pharmacol., 39, 319-327 (2002).
- 4. Nakahara, T., et al., Eur. J. Pharmacol., 389, 103-106 (2000).
- 5. Somlyo, A.V., et al., FASEB J., 17, 223-234 (2003).
- 6. Zhou, Y., et al., Science, 302, 1215-1217 (2003).
- 7. Kandabashi, T., et al., Arterioscler. Thromb. Vasc. Biol., 23, 2209-2214 (2003).
- 8. Masumoto, A., et al., Circulation, 105, 1545-1547 (2002).
- 9. Ishizaki, T., et al., Mol. Pharmacol., **57**, 976-983 (2000).

The Sigma-RBI eHandbook of Receptor Classification and Signal Transduction is now available online!

Receptor classification and signal transduction enzymes are just a click away...

The Updated and Expanded 4th Edition is now available online!

- Over 90 entries detailing the most current information on a range of receptors, ion channels, transporters and key cell signaling enzymes
- Each entry contains a peer-reviewed overview, 12 key references and a chart summarizing the classification of the subject receptors, ion channels, transporters or enzymes and detailing key ligands for their study

Visit the Sigma-RBI eHandbook of Receptor Classification and Signal Transduction online at sigma-aldrich.com/sigma-rbiehandbook

Also available in hard copy, 232 pages, spiral bound; to order your own copy (Prod. No. \$ 8188), call 800-325-3010 (US), fax 800-325-5052 (US) or contact your local Sigma-Aldrich office.



