

2007 Volume 2 Number 3



ENZYMES, KITS AND REAGENTS FOR ANALYSIS OF:

AGAROSE

ALGINIC ACID

CELLULOSE, LICHENEN AND GLUCANS

HEMICELLULOSE AND XYLAN

CHITIN AND CHITOSAN

**CHONDROITINS** 

DEXTRAN

HEPARANS

HYALURONIC ACID

INULIN

PEPTIDOGLYCAN

PECTIN

Cellulose, one of the most abundant biopolymers on earth, is a linear polymer of

 $\beta$ -(1-4)-D-glucopyranosyl units. Inter- and intra-chain hydrogen bonding is shown in red.

PULLULAN

STARCH AND GLYCOGEN

Complex Carbohydrate Analysis: Enzymes, Kits and Reagents



sigma-aldrich.com

# The Online Resource for Nutrition Research Products Only from Sigma-Aldrich





Designed to help you locate the chemicals and kits needed to support your work, the **Bioactive Nutrient Explorer** is an Internet-based tool created to aid medical researchers, pharmacologists, nutrition and animal scientists, and analytical chemists studying dietary plants and supplements.

The **Bioactive Nutrient Explorer** identifies the compounds found in a specific plant and arranges them by chemical family and class. You can also search for compounds having a similar chemical structure or for plants containing a specific compound.

The **Bioactive Nutrient Explorer** now includes a searchable database of plants listed by physiological activity in key areas of research, such as cancer, diabetes, metabolism and other disease or normal states. Plant Detail pages include common and Latin synonyms and display associated physiological activities, while Product Detail pages show the structure family and plants that contain the compound, along with comparative product information for easy selection. When you have found the product you need, a simple mouse click connects you to our easy online ordering system.

- New! Search for Plants Associated with Physiological Activity
- Locate Chemicals found in Specific Plants
- Identify Structurally Related Compounds



**Bioactive Nutrient Explorer** Helping Scientists Connect Bioactives to Botanicals sigma-aldrich.com/nutrition





2007 Volume 2 Number 3

### Table of Contents

### Complex Carbohydrate Analysis: Enzymes, Kits and Reagents

Agarase 2
Alginate Lyase 2
Cellulose, Lichenen and Glucan Degrading
Enzymes2-4
Cellulase
Driselase
β-Glucanase4
Laminarinase4
Lyticase4
Hemicellulose and Xylan Degrading
Enzymes5
Hemicellulase5
Xylanase5
Chitin and Chitosan Degrading Enzymes. 6-7
Chitinase6
Chitosanase7
Chondroitinases7-8
Chondroitinase ABC
Chondroitinase AC
Chondroitinase C
Dextranases
Heparinases
Heparinase I
Heparinase II 10
Heparinase III10
Hyaluronidases 10-11
Inulinase 11
Peptidoglycan Degrading Enzymes 12-13
Lysozyme 12
Mutanolysin13
Pectin Degrading Enzymes14-16
Pectinase15
Pectinesterase16
Pectolyase16
Pullulanases17
Starch and Glycogen Degrading Enzymes.18
α-Amylase19
β-Amylase
Amyloglucosidase 21
α-Glucosidase22
Kits for Carbohydrate Analysis23-24
Total Dietary Fiber23
Starch23
Glucose24
Fructose24
Complex Carbohydrates25-29

### Introduction

Complex carbohydrates compose the most abundant class of biopolymers on earth. Because of their structural and functional diversity, they have found applications in biomedical, nutritional, textile, cosmetic and countless other industries.

The complex carbohydrates of the extracellular matrices such as hyaluronic acid and chondroitin sulfate are finding utility in antinflammatory and cell proliferation applications. *In vivo*, hyaluronic acid forms a coating around chondrocytes in articular cartilage and together with the proteoglycan, aggrecan, is responsible for the uptake and retention of water. The two major glycan components of aggrecan are chondroitin sulfate and keratan sulfate. Hyaluronic acid may also interact with cell surface receptors, such as CD44, involved in lymphocyte activation. The degradation products of hyaluronic acid may also interact with Toll-like receptors in macrophages. Heparan sulfate is commonly found as a component of cell surface proteoglycans. It is also found in the extracellular matrix. Heparan sulfate appears to have a broad range of biological functions including regulation of thrombosis, growth factor signaling, cell proliferation, adhesion and mobility. Depending on its morphology, location and ligands, heparan sulfate may inhibit or promote metastasis. Heparan sulfate is known to bind several protein ligands. Most notably, its binding affinity with antithrombin has been extensively utilized in the form of the anticoagulant, heparin.

Dextrans also help to decrease vascular thrombosis. By binding to the endothelium, platelets and red blood cells, dextrans impart an electronegative environment in the blood vessel resulting in a reduction of red blood cell aggregation and platelet adhesion to the vascular endothelium. *In vivo*, dextran solutions have also been used for blood volume expansion. Conversely, chitosan has the ability to induce clot formation. It is used in wound healing, particularly as a coating for bandages. Chitin may also aid in wound healing by accelerating collagen production. Chitosan is also used to enhance plant growth and may help plants resist fungal infection.

Chitin and starch are used as binders in the paper, dye, textiles and adhesives industries. Chitin and chitosan are also used as filtration aids, particularly in the waste water treatment industry. Chitosan aids in particulate aggregation as well as removal of phosphorus, metals, and grease from waste water. Modified agarose, chitin, starch and dextrans have been manipulated to produce media with controlled pore size for chemical separations. Beaded forms of cross-linked agarose and dextrans are the components of size exclusion, ion exchange and affinity chromatography media. Agarose and soluble starch are commonly used as electrophoresis media.

In the food industry, starches, aglinates, agarose, chitins, chitosans and pectins are used as gelling, thickening and encapsulating agents. Pectins and inulins are common components of dietary fiber supplements and may help to increase nutrient uptake. Pullulan is a common component of edible films.

### The Enzyme Explorer

### Your Comprehensive Source for Products and Technical Resources for Glycobiology sigma-aldrich.com/enzymeexplorer

### Carbohydrate Analysis

- · Complex Carbohydrate and Polysaccharide Analysis
- Proteoglycan and Glycoprotein Analysis

### Carbohydrate Metabolism

- Carbohydrate Metabolite Library
- Enzymes Involved in Carbohydrate Metabolism
- Metabolic Pathway Charts and Animations


### Enzymatic-Based Kits for the Quantitation of Carbohydrates

- Total Dietary Fiber
- Starch
- Glucose
- Fructose
- Sucrose

For additional technical information including literature citations pertaining to the content in this publication, visit the Enzyme Explorer's "Enzymatic Carbohydrate Analysis Resource"

sigma-aldrich.com/enzymeexplorer The Enzyme Research Resource

For Hazard Information and other information please refer to the Sigma Biochemicals, Reagents and Kits for Life Science Research Catalog or sigma-aldrich.com



Agarose is the principal neutral gelling component of agar extracted from algae. Agarose is a complex range of polysaccharide chains composed of alternating  $\alpha$ -(1-3)-D-galactosyl- $\beta$ -(1-4)-anhydro-L-galactosyl units.

### Agarase

Agarase catalyzes the hydrolysis of 1,3- $\beta$ -D-galactosidic linkages in agarose, giving the tetramer as the predominant product.

### Agarase from Pseudomonas atlantica

Agarose 3-glycanohydrolase [37288-57-6] E.C. 3.2.1.81


### lyophilized powder, activity: 1,000-3,000 units/mg solid

Contains phosphate buffer salts. May contain bovine serum albumin to standardize protein content

One unit will produce 1.0 µg of reducing sugar (measured as p-galactose) from agar per min at pH 6.0 at 40 °C.

A6306-1KU	1,000 units
A6306-5KU	5,000 units

# **Alginate Lyase**



### Alginic Acid

Composed of blocks of  $\beta$ -D-(1-4) mannuronic acid homopolymeric regions (MMMM...),  $\alpha$ -L-(1-4)-guluronic acid (GGGG...) homopolymeric regions, and alternating coplolymer regions of  $\beta$ -D-(1-4) mannuronic acid -  $\alpha$ -L-(1-4)-guluronic acid (GMGMGM...). Bacterial alginic acid can be acetylated at the 2 or 3 positions on mannuronic acid.

Alginate lyase cleaves at the  $\beta$ -(1-4)-D-mannuronic bonds residues to yield oligosaccharides with 4-deoxy- $\alpha$ -L-erythro-hex-4-enopyranuronosyl groups at their non-reducing terminus.

### Alginate Lyase from Flavobacterium sp.

### E.C. 4.2.2.3 powder


### activity: ≥10,000 units/g solid

Add 0.15 mL of enzyme solution (1 un/mL) to 4.5 mL of 0.1% sodium alginate (pH 6.3). Incubate at 37 °C for 30 minutes. Terminate reaction by addition of 4.65 mL of 0.1 N NaOH.


One unit will produce an increase the A235 nm of 1.0 per minute per mL of sodium alginate solution at pH 6.3 at 37  $^\circ\mathrm{C}$ 

2-8°C

# Cellulose, Lichenan and Glucan Degrading Enzymes




Polymer of β-(1-4)-p-glucopyranosyl units




### Lichenen

Repeating linear polymer of two  $\beta\text{-}(1\text{-}4)\text{-}D\text{-}glucopyranosyl and one }\beta\text{-}(1\text{-}3)\text{-}D\text{-}glucopyranosyl unit.}$ 



### Cerial β-Glucan

Polymer of  $\beta$ -(1-4)-p-glucopyranosyl units occuring as predominantly as cellotriose and cellotetraose separated by single  $\beta$ -(1-3)-p-glucopyranosyl units. Cross-linking can occur within the consecutive cellotriose regions.



Yeast β-Glucan

Polymer of  $\beta\text{-}(1\text{-}3)\text{-}D\text{-}glucopyranosyl units with branching at }\beta\text{-}(1\text{-}6)\text{-}D\text{-}glucopyranosyl units.}$ 

2

sigma-aldrich.com/biofile

Cellulase catalyzes the endohydrolysis of 1,4- $\beta$ -D-glucosidic linkages in cellulose, lichenin and cereal  $\beta$ -D-glucans

### Cellulase from Aspergillus sp.

### Carezyme® 1000L

[9012-54-8] E.C. 3.2.1.4

### activity: $\geq$ 1000 U/g

Produced by submerged fermentation of a genetically modified *Aspergillus* microorganism

A product of Novozyme Corp.

2-8°C
-------

C2605-50ML	50 mL
C2605-250ML	250 mL

### Cellulase from Aspergillus niger

1,4-(1,3:1,4)-β-D-Glucan 4-glucano-hydrolase [9012-54-8] E.C. 3.2.1.4 EC No. 2327344

### powder, activity: ≥0.3 units/mg solid

One unit will liberate 1.0  $\mu mole$  of glucose from cellulose in 1 hr at pH 5.0 at 37  $^{\circ} C$  (2 hr incubation time).

### 2-8°C

C1184-5KU	5,000 units
C1184-25KU	25,000 units
C1184-100KU	100,000 units

### Cellulase from Trichoderma reesei ATCC 26921

1,4-(1,3:1,4)-β-D-Glucan 4-glucano-hydrolase [9012-54-8] E.C. 3.2.1.4 Celluclast<sup>®</sup> 1.5L

### ▶ aqueous solution, activity: ≥700 U/g

A product of Novozyme Corp.

2-8°C

C2730-50ML

### ▶ lyophilized powder, activity: ≥1 unit/mg solid

One unit will liberate 1.0  $\mu mole$  of glucose from cellulose in 1 hr at pH 5.0 at 37  $^{\circ} C$  (2 hr incubation time).

2-8°C

C8546-2.5KU	2,500 units
C8546-5KU	5,000 units
C8546-10KU	10,000 units

### Cellulase from Trichoderma viride

1,4-(1,3:1,4)-β-D-Glucan 4-glucano-hydrolase [9012-54-8] E.C. 3.2.1.4 EC No. 2327344

### plant cell culture tested, activity: 3-10 units/mg solid Composition: protein ~50% (biuret); contains lactose and glucose

One unit will liberate 1.0  $\mu mole$  of glucose from cellulose in 1 hr at pH 5.0 at 37  $^{\circ} C$  (2 hr incubation time).

### 2-8°C

C1794-5KU	5,000 units
C1794-10KU	10,000 units

# crude powder, activity: 3-10 units/mg solid Composition: protein ~50% (biuret)

One unit will liberate 1.0  $\mu mole$  of glucose from cellulose in 1 hr at pH 5.0 at 37  $^{\circ} C$  (2 hr incubation time).

### 2-8°C

C9422-5KU	5,000 units
C9422-10KU	10,000 units

### ▶ Onozuka RS, powder, activity: ≥5,000 units/g solid

### Manufactured by Yakult

One unit will liberate 1.0  $\mu mole$  of glucose from cellulose in 1 hr at pH 5.0 at 37  $^{\circ} C$  (2 hr incubation time).

2-8°C	
C0615-1G	1 g

### Driselase

### Driselase from Basidiomycetes sp.

[85186-71-6]

50 mL

### powder, Protein: ~15%

Crude powder containing laminarinase, xylanase and cellulase.

D9515-1G	1 g
D9515-5G	5 g
D9515-25G	25 g

Ů

IGMA

# Cellulose, Lichenan and Gluan Degradins Enzymes

Laminarinase catalyzes the endohydrolysis of 1,3- or 1,4-linkages in  $\beta$ -D-glucans when the glucose residue whose reducing group is involved in the linkage to be hydrolyzed is itself substituted at C-3.

### β-Glucanase from Aspergillus niger

[9074-98-0] E.C. 3.2.1.6 EC No. 2329802

### BioChemika, powder, dark-brown, activity: ~1 unit/mg

One Unit corresponds to the amount of enzyme which will release 1  $\mu$ moL of reducing sugar equivalents (expressed as glucose) per minute at pH 5.0 and 55 °C, using  $\beta$ -D-glucan (Catalog No. 49102) as substrate

Γ	2-	8°	С

49101-100MG	100 mg
49101-500MG	500 mg

### β-1,3-D-Glucanase from Helix pomatia

[9044-93-3] E.C. 3.2.1.39 EC No. 2329273

### BioChemika, powder, light beige, activity: 0.5-1.5 units/mg

One Unit corresponds to the amount of enzyme which liberates 1  $\mu mol$  glucose from laminarin (Catalog No. 61340) per minute at pH 5.0 and 37 °C.

Improved filterability of wines by enzymic decomposition of carbohydrate-containing colloids<sup>1</sup>; Induction of hydrolases as a defense reaction against pathogens, review<sup>2</sup>

Lit. cited: 1. Wucherpfennig, K., and Dietrich, H. , Weinwirtschaft 118, 598 (1982)

2. Boller, T. , UCLA Symp. Mol. Cell Biol., New Ser. 22, 247 (1985)

49103-10MG	10 mg
49103-50MG	50 mg

### Laminarinase from Penicillium sp.

endo-1,3(4)-β-glucanase; 1,3-(1,3:1,4)-β-D-Glucan 3(4)-glucanohydrolase [62213-14-3] E.C. 3.2.1.6

### lyophilized powder, activity: 5-10 units/mg protein

Lyophilized powder containing acetate buffer salts

composition

Protein ~70% (biuret)

One unit will liberate 1.0 mg of reducing sugar (measured as glucose) from laminarin per min at pH 5.0 at 37  $^\circ \rm C.$ 

Contains cellulase and  $\alpha$ -amylase

2-8°C

### L9259-25UN

### Laminarinase from Trichoderma sp.

endo-1,3(4)- $\beta$ -glucanase; 1,3-[1,3;1,4]- $\beta$ -D-Glucan 3(4)-glucanohydrolase [62213-14-3] E.C. 3.2.1.6

### powder, activity: 100-400 units/g solid

Contains chitinase activity

One unit will liberate 1.0 mg of reducing sugar (measured as glucose) from laminarin per min at pH 5.0 at 37  $^\circ\mathrm{C}.$ 

Contains cellulase and  $\alpha$ -amylase

### 2-8°C

Ð

m / b i o f i l

sigma-aldrich.co

L5272-5UN	5 units
L5272-25UN	25 units

### Lyticase

Lyticase hydrolyzes poly- $\beta$ -(1-3)-glucose such as yeast cell wall  $\beta$ -glucan.

### Lyticase from Arthrobacter luteus

### [37340-57-1]

Yeast cells are difficult to disrupt because the cell walls may form capsules or resistant spores. DNA can be extracted from yeast by using lysing enzymes such as lyticase, chitinase, zymolase, and gluculase to induce partial spheroplast formation: spheroplasts are subsequently lysed to release DNA. Lyticase is preferred to digest cell walls of yeast and generate spheroplasts from fungi for transformation. Reported to be useful for lysis of *Ashbya, Candida, Debaryomyces, Eremothecium, Endomyces, Hansenula, Hanseniaspora, Kloeckera, Kluyveromyces, Lipomyces, Saccharomycopsis, Saccharomycodes, and Schwanniomyces* species.

One unit will produce a  $\Delta A800$  of 0.001 per min at pH 7.5 at 25 °C, using a suspension of yeast as substrate in a 3 mL reaction mixture.

# $\blacktriangleright$ lyophilized powder, activity: ${\geq}2,000$ units/mg protein, Protein: ${\geq}20\%$

Partially purified, lyophilized powder containing potassium phosphate buffer salts and stabilizers

-20°C

L2524-10KU	10,000 units
L2524-25KU	25,000 units
L2524-50KU	50,000 units
L2524-200KU	200,000 units

### ▶ lyophilized powder, activity: ≥200 units/mg solid

-20°C

L4025-25KU	25,000 units
L4025-50KU	50,000 units
L4025-100KU	100,000 units
L4025-250KU	250,000 units
L4025-1MU	1,000,000 units

### ▶ partially purified powder, activity: ≥2,000 units/mg protein

Partially purified powder containing ammonium sulfate and stabilizer

### composition

25 units

Protein ~20% (biuret)

2-8°C	
L5263-25KU	25,000 units
L5263-50KU	50,000 units
L5263-200KU	200,000 units

Enzymes

### Lyticase from Oerskovia xanthineolytica

### [37340-57-1]

# recombinant, expressed in *Escherichia coli*, lyophilized powder

Purified recombinant  $\beta\text{-}(1,3)\text{-}glucanase preparation that is protease-free. Vial of <math display="inline">\geq\!500$  units.


One unit will produce a  $\Delta A800$  of 0.001 per min at pH 7.5 at 25 °C, using a suspension of yeast as substrate in a 3 mL reaction mixture.

An exceptionally stable enzyme preparation with very low levels of nucleic acid and nuclease contamination.

### -20°C

L4276-1VL	1 vial

# Hemicellulose and Xylan Degrading Enzymes



Polymer of  $\beta$ -(1-4)-D-xylopyranosyl units

Xylan hemicelluloses are a group of plant-derived heteropolysaccharides associated with cellulose and lignin. The most common hemicelluloses are: xylan, glucuronoxylan, arabinoxylan, glucomannan and xyloglucan. In angiosperms, the principal hemicellulose component, xylan, is a polymer of  $\beta$ -(1-4)- p-xylopyranose. In arabinoxylan, branching occurs at the C2 & C3 positions with  $\alpha$ -L-arabinofunaose. Glucuronoxylan, also found in angiosperms, has the xylan backbone with 4-0 methylglucuronic acid branching. In addition, arabinose branching as well as acetylation may be present. Gymnosperms contain glucomannans comprised primarily of p-mannosyl and p-glucosyl residues.

### Hemicellulase

### Hemicellulase from Aspergillus niger

[9025-56-3]

# powder, activity: 0.3-3.0 units/mg solid (using a $\beta$ -galactose dehydrogenase system and locust bean gum as substrate)

An undefined mixture of glycolytic enzymes usually containing xylanase, mananase and other activities.

Contains lactose as standardization of activity

One unit will produce a relative fluidity change of 1 per 5 minutes using locust bean gum as substrate at pH 4.5 at 40  $^\circ\mathrm{C}$ 

-20°C

H2125-150KU

150,000 units

### **Xylanase**

Xylanase catalyzes the endohydrolysis of  $\beta$ -(1-4)-D-xylosidic linkages in xylans yielding various  $\beta$ -(1-4)-D-xylooligosaccharides.

### Xylanase from Thermomyces lanuginosus

**Pentopan Mono BG**<sup>®</sup> [37278-89-0]

# powder, activity: ≥2500 units/g, recombinant, expressed in Aspergillus oryzae

Purified endo (1,4)-β-xylanase from *Thermomyces lanuginosus*. Produced by submerged fermentation of a genetically modified *Aspergillus oryzae* microorganism.

A product of Novozyme Corp.

2-8°C

X2753-10G	10 g
X2753-50G	50 g

### Xylanase from *Trichoderma viride*

endo-1,4- $\beta$ -Xylanase; 1,4- $\beta$ -D-Xylanxylanohydrolase [9025-57-4] E.C. 3.2.1.8 EC No. 2534397

### **lyophilized powder, activity: 100-300 units/mg protein** Contains sorbitol and sodium acetate buffer salts

composition

Protein ~50% (biuret)

### 2-8°C


X3876-250UN	250 units
X3876-1KU	1,000 units

Ú

IGMA

# **Chitin and Chitosan Degrading Enzymes**

### Chitinase



Polymer of  $\beta$ -(1-4)-N-Acetyl-D-glucosamine units

Chitinases have been detected in many organisms, including bacteria, fungi, plants, invertebrates and vertebrates. Chitinases are broadly classified as endo- and exochitinases. The endochitinase activity is defined as the random cleavage at internal points in the chitin chain. The exochitinase activity is defined as the progressive action starting at the non reducing end of chitin with the release of chitobiose or N-acetylglucosamine units. Chitobiosidase and N-acetyl-β-glucosaminidase are considered exochitinases. The combination of endo- and exochitinases results in a synergistic increase in the chitinolytic activity.

### Chitinase from Serratia marcescens

Chitodextrinase; Poly(1,4-β-[2-acetamido-2-deoxy-D-glucoside]) glycanohydrolase [9001-06-3] E.C. 3.2.1.14

lyophilized powder, activity: 400-1,200 units/g solid

Lyophilized powder containing phosphate buffer salts composition

Protein 20-40% (biuret)

One unit will liberate 1.0 mg of N-acetyl-D-glucosamine from chitin per hour at pH 6.0 at 25  $^\circ\rm C$  in a 2 hour assay.

### -20°C

C7809-1UN	1 unit
C7809-5UN	5 units
C7809-10UN	10 units

### Chitinase from Streptomyces griseus

Chitodextrinase; Poly(1,4-β-[2-acetamido-2-deoxy-D-glucoside]) glycanohydrolase

[9001-06-3] E.C. 3.2.1.14 EC No. 2325787

# lyophilized powder (Essentially salt free), activity: 200-800 units/g solid

One unit will liberate 1.0 mg of N-acetyl-D-glucosamine from chitin per hour at pH 6.0 at 25  $^\circ\rm C$  in a 2 hour assay.

### -20°C

C6137-5UN	5 units
C6137-25UN	25 units
C6137-50UN	50 units

### Chitinase from Trichoderma viride

N-acetyl- $\beta$ -glucosaminidase and chitodextrinase E.C. 3.2.1.14 and 3.2.1.52

### lyophilized powder, activity: ≥600 units/g solid

The chitinolytic enzymes from *T. viride* are a mixture of extracellular chitinolytic enzymes, which exhibit exo- and endochitinase activities including N-acetyl- $\beta$ -glucosaminidase and chitobiosidase.


One unit will liberate 1.0 mg of N-acetyl-D-glucosamine from chitin per hour at pH 6.0 at 25 °C in a 2 hour assay.

C8241-25UN

Enzymes

SIGMA

### Chitosanase



Polymer of  $\beta$ -(1-4)-D-glucosamine units

Chitosanase catalyzes the endohydrolysis of  $\beta$ -(1-4)-linkages between p-glucosamine (GlcN-GlcN) residues in chitosan. The enzyme from Streptomyces has been reported to also hydrolyze the GlcNAc-GlcN linkage in partially acetylated chitosan.

### Chitosanase from Streptomyces sp.

Chitosan N-acetylglucosaminohydrolase [51570-20-8] E.C. 3.2.1.132

# buffered aqueous glycerol solution, activity: $\geq$ 15 units/mg protein

Solution in 50% glycerol containing 100 mM sodium acetate, pH 5.0  $\,$ 

One unit will liberate 1.0  $\mu mole$  of reducing sugar (measured as p-glucosamine equivalents) from chitosan per minute at pH 5.5 at 37 °C.

### -20°C

10 units

### Chitosanase from Streptomyces griseus

Chitosan N-acetylglucosaminohydrolase [51570-20-8] E.C. 3.2.1.132

# lyophilized powder, activity: >50 units/mg protein (Bradford)


Lyophilized powder containing potassium phosphate buffer salts. Purified by chromatography

One unit will release 1  $\mu mole$  of glucosamine from chitosan per min at pH 5.0 at 37  $^{\circ} C$ 

chitinase	<1.0	unit/mg protein
-20°C		


C9830-10UN	10 units
------------	----------

### **Chondroitinases**




### Chondroitin Sulfate A

Alternating copoly (β-glucuronic acid-(1-3))-N-acetyl-β-galactosamine-4-sulfate



### Chondroitin Sulfate B

Alternating copoly ( $\beta$ -Iduronic acid-(1-3))-N-acetyl- $\beta$ -galactosamine-4-sulfate



Alternating copoly (β-glucuronic acid-(1-3))-N-acetyl-β-galactosamine-6-sulfate

sigma-aldrich.com/biofiles

SIGMA

# Chondroitinases

### **Chondroitinase ABC**

Chondroitinase ABC catalyzes the eliminative degradation of polysaccharides containing  $\beta$ -(1-4)-D-hexosaminyl and  $\beta$ -(1-3)-D-glucuronosyl or  $\alpha$ -(1-3)-L-iduronosyl linkages to disaccharides containing 4-deoxy- $\beta$ -D-gluc-4-enuronosyl groups. It acts on chondroitin 4-sulfate, chondroitin 6-sulfate, dermatan sulfate, and acts slowly on hyaluronate. Initial rates of degradation of chondroitin sulfate B, chondroitin, and hyaluronic acid were 40%, 20%, and 2%, respectively, that of chondroitin sulfate A and chondroitin sulfate C.

### Chondroitinase ABC from Proteus vulgaris

Chondroitin ABC Lyase [9024-13-9] E.C. 4.2.2.4 EC No. 2327779

### lyophilized powder, activity: 50-250 units/mg protein (usingchondroitin sulfate C as substrate)

Contains potassium phosphate buffer salts and stabilizer; BSA-free.

### composition, protein ~10% (Lowry) Affinity purified

Packages based on chondroitinase C

One unit will liberate 1.0 µmole of 2-acetamido-2-deoxy-3-O-( $\beta$ -D-gluc-4-ene-pyranosyluronic acid)-4-O-sulfo-D-galactose from chondroitin sulfate A or 1.0 µmole of 2-acetamido-2-deoxy-3-O-( $\beta$ -D-gluc-4-ene-pyranosyluronic acid)-6-O-sulfo-D-galactose from chondroitin sulfate C per min at pH 8.0 at 37 °C.

-20°C	
C3667-5UN	5 units
C3667-10UN	10 units

protease ..... essentially free

### **Chondroitinase AC**

Chondroitinase AC is an eliminase that degrades chondroitin sulfates A and C, but not chondroitin sulfate B. The enzyme cleaves, via an elimination mechanism, sulfated and non-sulfated polysaccharide chains containing  $\beta$ -(1-4) and  $\beta$ -(1-3) linkages between hexosamines and glucuronic acid residues. The reaction yields oligosaccharide products, mainly disaccharides, containing unsaturated uronic acids that can be detected by UV spectroscopy at 232 nm. The enzyme shows approximately equal activity with chondroitin sulfates A and C, while the activity observed with chondroitin sulfate B is approximately 7% of this value. This activity is most likely due to the presence of chondroitin sulfates A and C (10%) in the chondroitin sulfate B.

### Chondroitinase AC from Flavobacterium heparinum

Chondroitin AC lyase [9047-57-8] E.C. 4.2.2.5

### lyophilized powder, activity: 0.5-1.5 units/mg solid

(using chondroitin sulfate A as substrate, also cleaves chondroitin sulfate C)  $% \left( {{\mathcal{L}}_{{\rm{s}}}^{\rm{c}}} \right)$ 

Contains potassium phosphate buffer salts and BSA as stabilizer. composition, protein ~15% (Lowry)

One unit will cause a  $\Delta A_{_{232}}$  of 1.0 per minute due to the release of unsaturated disaccharide from chondroitin sulfate A at pH 7.3 at 37 °C. Reaction volume: 3.1 mL (light path 1 cm).

Glycosaminoglycan (GAG) degradation enzymes ..... may contain trace amount

5 units

### C2780-5UN

### Chondroitinase C

Chondroitinase C cleaves chondroitin sulfate C producing tetrasaccharide plus an unsaturated 6-sulfated disaccharide (delta Di-6S). It also cleaves hyaluronic acid producing unsaturated nonsulfated disaccharide ( $\Delta$  Di-OS). Chondroitin sulfate A is also degraded producing oligosaccharides and delta Di-6S, but not delta Di-4S. Chondroitinase C cleaves the GalNAc bond of the pentasaccharides or hexasaccharides derived from the linkage region of chondroitin sulfate chains and tolerates sulfation of the C-4 or C-6 of the GalNAc residue and C-6 of the Gal residues, as well as 2-O-phosphorylation of the Xyl residue. In contrast, it does not act on a GalNAc–GlcA linkage when attached to a 4-O-sulfate Gal residue.

### Chondroitinase C from Flavobacterium heparinum

Chondroitin C lyase [60184-91-0] E.C. 4.2.2.—


-20°C

### lyophilized powder, activity: ≥200 units/mg solid

One unit will form 0.1  $\mu mole$  of unsaturated uronic acid per hr at pH 8.0 at 25  $^\circ C$  using chondroitin sulfate C as substrate.

C0954-75UN	75 units

### Dextranases



Dextran is composed of approximately 95%  $\alpha$ -(1-6)-D-linkages. The remaining  $\alpha$ -(1-3) linkages account for the branching of dextran. Conflicting data on the branch lengths implies that the average branch length is less than three glucose units. However, other methods indicate branches of greater than 50 glucose units exist. Native dextran has been found to have a molecular weight (MW) in the range of 9 million to 500 million Da. Lower MW dextrans will exhibit slightly less branching and have a more narrow range of MW distribution. Dextrans with MW greater than 10,000 Da behave as if they are highly branched. As the MW increases, dextran molecules attain greater symmetry. Dextrans with MW of 2,000 to 10,000 Da exhibit the properties of an expandable coil. At MWs below 2,000 Da dextran is more rod-like.

Φ

biofil

ر ۲

с С

sigma-aldrich.

### Dextranase

Dextranase catalyzes the endohydrolysis of  $\alpha\text{-}(1\text{-}6)\text{-}D\text{-}glucosidic linkages in dextran.}$ 

### Dextranase from Chaetomium erraticum

1,6-α-D-Glucan 6-glucanohydrolase E.C. 3.2.1.11

A product of Novozymes Corp.

### solution

A fungal dextranase produced by submerged fermentation of *Chaetomium erraticum*.

Stable in the pH range of 3-7 and at temperatures up to approx. 70 °C. For most applications, the preferred conditions are pH 5-6 and a temperature of 50-60 °C.

2-8°C
-------

D0443-50ML	50 mL
D0443-250ML	250 mL

### Dextranase from Penicillium sp.

1,6-α-D-Glucan 6-glucanohydrolase

E.C. 3.2.1.11

One unit will liberate 1.0  $\mu mole$  of isomaltose (measured as maltose) per min at pH 6.0 at 37  $^\circ C$ , using dextran as substrate.

### Iyophilized powder, activity: 400-800 units/mg protein composition, protein 35% (Lowry)

2-8°C

D8144-500UN	500 units
D8144-1KU	1,000 units

### Iyophilized powder, activity: 100-250 units/mg protein

Partially purified, lyophilized powder

composition, protein ~25% (Lowry)

2-8°C

D4668-500UN	500 units
D4668-1KU	1,000 units

### Iyophilized powder, activity: 10-25 units/mg solid


Crude

2-8°C

D5884-5KU	5,000 units
D3004-3KU	5,000 units

### **Heparinases**

### Heparinase Specificities



Heparan and heparin glycosaminoglycans (GAGs) are complex heterogeneous mixtures of repeating disaccharide units consisting of a uronic acid (D-glucuronic or L-iduronic acid) and D-glucosamine or N-acetyl-D-glucosamine. Various degrees of sulfation occur (at O and/or N) on each monosacchiride unit, ranging from zero to tri-sulfation. In general, heparan is less sulfated than heparin.

Heparinase selectively cleaves sulfated glycans containing  $\alpha$ -(1-4)-glycosidic linkages between the glucosamine and uronic acid residues in the heparin polymer. The cleavage proceeds via an elimination reaction, resulting in the formation of oligosaccharides containing unsaturated uronic acid residues (double bond between C4 and C5). These cleavage products can be detected by UV spectroscopy (232 nm). The three forms of heparinase (I, II, and III) have varying substrate specificities.

### Heparinase I

Heparinase I cleaves heparin and heparan sulfate (relative activity about 3:1) at the linkages between hexosamines and O-sulfated iduronic acids, yielding mainly disaccharides. The enzyme also cleaves the antithrombin III binding pentasaccharide domain in the heparin molecule.

### Heparinase I from Flavobacterium heparinum

Heparinase; Heparin Iyase I [9025-39-2] E.C. 4.2.2.7

# Lyophilized powder stabilized with approx. 25% bovine serum albumin, activity: 200-600 units/mg solid

mol wt 42.8 kDa

One unit will form 0.1  $\mu mole$  of unsaturated uronic acid per hr at pH 7.5 at 25 °C. One International Unit (I.U.) is equivalent to approx. 600 Sigma units.

-20°C	

H2519-50UN	50 units
H2519-100UN	100 units
H2519-250UN	250 units

Ú

**GMA** 

### Heparinases

### Heparinase II

Heparinase II cleaves heparan sulfate, and to a lesser extent heparin (relative activity about 2:1), at the  $\alpha$ -(1-4) linkages between hexosamines and uronic acid residues (both glucuronic and iduronic), yielding mainly disaccharides.

### Heparinase II from Flavobacterium heparinum

Heparin lyase II

[149371-12-0]

Lyophilized powder stabilized with approx. 25% bovine serum albumin, lyophilized powder, activity: 100-300 units/mg solid mol wt 84.1 kDa

One unit will form 0.1  $\mu$ mole of unsaturated uronic acid per hr at pH 7.0 at 25 °C. One International Unit (I.U.) is equivalent to approx. 600 Sigma units.

### -20°C

H6512-10UN	10 units
H6512-25UN	25 units
H6512-100UN	100 units

### Heparinase III

Heparinase III cleaves at the  $\alpha$ -(1-4) linkages between hexosamine and glucuronic acid residues in heparan sulfate, yielding mainly disaccharides. The enzyme is not active towards heparin. Sulfation at the 6-position of glucosamine inhibits cleavage by heparinase III,

### Heparinase III from Flavobacterium heparinum

Heparin Lyase III; Heparitinase from *Flavobacterium heparinum*; Heparitinase I [37290-86-1] E.C. 4.2.2.8

# Lyophilized powder stabilized with approx. 25% bovine serum albumin, activity: 200-600 unit/mg solid


Heparin-degrading lyase that recognizes heparin sulfate proteoglycan as its primary substrate.

mol wt 70.8 kDa

One unit will form 0.1 µmole of unsaturated uronic acid per hr at pH 7.5 at 25 °C. One International Unit (I.U.) is equivalent to approx. 600 Sigma units.

H8891-5UN	5 units
H8891-10UN	10 units
H8891-50UN	50 units

### **Hyaluronidases**



### Hyaluronic Acid

Composed of alternating residues of  $\beta$ -D-(1-3) glucuronic acid and  $\beta$ -D-(1-4)-N-acetylglucosamine

The mammalian hyaluronidases (EC 3.2.1.35) cleave hyaluronic acid and similar glycosaminoglycans by hydrolysis. The enzyme from *Streptomyces* (EC 4.2.2.1) is a lyase that catalyzes cleavage by an elimination reaction yielding a 4-deoxy-4,5-unsaturated oligosaccharides. Its specificity towards chondroitins and other glycosaminoglycans is unclear.

### Mammalian Hyaluronidase

The mammalian glycolytic hyaluronidase (EC 3.2.1.35) catalyzes the random hydrolysis of the 1-4 bond between N-acetyl-D-glucosamine and D-glucuronic acid in hyaluronic acid. It also hydrolyzes  $\beta$ -(1-4)-D-glycosidic linkages between N-acetyl-galactosamine or N-acetylgalactosamine sulfate and glucuronic acid in chondroitin sulfates A and C, and dermatan.

### Hyaluronidase from bovine testes

Hyaluronate 4-glycanohydrolase; Hyaluronoglucosaminidase [37326-33-3] E.C. 3.2.1.35 EC No. 2534643

mol wt ~55 kDa (four subunits of 14 kDa each)

One unit is based on the change in absorbance at 600 nm (change in turbidity) of a USP reference standard hyaluronidase which is assayed concurrently with each lot.

### Type I-S, lyophilized powder, activity: 400-1000 units/mg solid

-20°C

H3506-100MG	100 mg
H3506-500MG	500 mg
H3506-1G	1 g
H3506-5G	5 g

Type IV-S, lyophilized powder (essentially salt-free), activity: 750-1500 units/mg solid

-20°C

H3884-50MG	50 mg
H3884-100MG	100 mg
H3884-500MG	500 mg
H3884-1G	1 g

Enzymes

SIGNA

### Type VIII, lyophilized powder, activity: ~300 units/mg

Prepared from sterile filtered solution of Type I-S.

-20°C

H3757-100MG 100 mg

### Type VI-S, lyophilized powder, activity: 3,000-15,000 units/mg solid

Chromatographically purified

-20°C

200	
H3631-3KU	3,000 units
H3631-15KU	15,000 units
H3631-30KU	30,000 units

### Hyaluronidase from sheep testes

Hyaluronate 4-glycanohydrolase: Hyaluronoglucosaminidase [37326-33-3] E.C. 3.2.1.35 EC No. 2534643

mol wt 55 kDa

One unit is based on the change in absorbance at 600 nm (change in turbidity) of a USP reference standard hyaluronidase which is assayed concurrently with each lot.

### ► Type V, lyophilized powder, activity: ≥1,500 units/mg solid 20°C

H6254-500MG	500 mg
H6254-1G	1 g

### ▶ Type II, lyophilized powder, activity: ≥300 units/mg

Lyophilized powder containing lactose

H2126-100MG	100 mg
H2126-500MG	500 mg
H2126-1G	1 g
H2126-5G	5 g

### ▶ Type III, lyophilized powder, activity: ≥500 units/mg

Lyophilized powder containing 20-50% lactose

-20°C

H2251-100MG	100 mg
H2251-500MG	500 mg
H2251-1G	1 g
H2251-5G	5 g

# Hyaluronate Lyase (*Streptomyces Hyaluronidase*)

### Hyaluronidase from Streptomyces hyalurolyticus

Hyaluronate Lyase from *Streptomyces hyalurolyticus* [9001-54-1] E.C. 4.2.2.1 EC No. 2326141

### lyophilized powder

Hyaluronate lyase cleaves hyaluronic acid at the  $\beta$ -D-GalNAc-(1-4)- $\beta$ -D-GlcA bond, yielding 3-(4-deoxy- $\beta$ -D-gluc-4-enuronosyl)-N-acetyl-D-glucosamine tetra- and hexasaccharides. Unlike other hyaluronidases, this enzyme is specific for hyaluronic acid and is inactive with chondroitin and chondroitin sulfate.<sup>1</sup>


Lit. cited: 1. Ohya, T., and Kaneko, Y., Biochim. Biophys. Acta 198, 607 (1970)

-20°C

H1136-1AMP

### Inulinase

### **Inulinase Specificity**



### Inulin

Inulins are fructan oligosaccharides composed  $\alpha$ -D-glucopyranosyl-[ $\beta$ -(2-1) D-fructofuranosyl-D-fructofuranosides. Inulins can generally contain 2 to 140 fructose units.

Inulinase	from	Aspergil	lus niger

I2017-250ML

Fructozyme L<sup>™</sup>; Inulase E.C. 3.2.1.7 Inulinase catalyzes endohydrolysis of β-(2-1)-D-fructosidic linkages in inulin. A product of Novozyme Corp. [2=8°C] [2017-50ML 50 mL

s i g


ma-aldrich.com/biofiles

SIGMA

1 amp

# Peptidoglycan Degrading Enzymes

### Lysozyme and Mutanolysin Specificity



(direct cross link in gram negative)

### Peptidoglycan

Polymer of  $\beta$ -(1-4)-N-Acetyl-D-glucosamine units. Alternating residues are modified to form N-acetylmuramic acid with the addition of lactate to form branching links to the tetrapeptide.

### Lysozyme

### Lysozyme from chicken egg white

Mucopeptide N-acetylmuramoylhydrolase; Muramidase [12650-88-3] E.C. 3.2.1.17

Lysozyme hydrolyzes  $\beta$ -(1-4) linkages between N-acetylmuramic acid and N-acetyl-D-glucosamine residues in peptidoglycan and between Nacetyl-D-glucosamine residues in chitodextrin. Grampositive cells are quite susceptible to this hydrolysis as their cell walls have a high proportion of peptidoglycan. Gram-negative bacteria are less susceptible due to the presence of an outer membrane and a lower proportion of peptidoglycan. However, these cells may be hydrolyzed in the presence of EDTA that chelates metal ions in the outer bacterial membrane.

The enzyme is active over a broad pH range (6.0 to 9.0). At pH 6.2, maximal activity is observed over a wider range of ionic strengths (0.02 to 0.100 M) than at pH 9.2 (0.01 to 0.06 M).

Used to prepare spheroplasts.

Single-chain mol wt 14.7 kDa

One unit will produce a  $\Delta A_{450}$  of 0.001 per min at pH 6.24 at 25 °C, using a suspension of *Micrococcus lysodeikticus* as substrate, in a 2.6 mL reaction mixture (1 cm light path).

### Iyophilized powder, Protein: 95% activity: 50,000 upits/mg r

### ~95%, activity: ~50,000 units/mg protein

Dialyzed and lyophilized, containing buffer salts as sodium acetate and sodium chloride

3× crystallized

-20°C	
L6876-1G	1 g
L6876-5G	5 g
L6876-10G	10 g
L6876-25G	25 g
L6876-100G	100 g

# Iyophilized powder, Protein: ~95%, activity: ~50,000 units/mg protein (E1%282)

### Features and Benefits

• Highly purified by repeated crystallization and dialysis

• Each lot is use-tested for isolation of plasmid DNA from *E. coli* essentially salt-free

3× crystallized

L7651-1G	1 g
L7651-5G	5 g
L7651-10G	10 g
L7651-25G	25 g
L7651-100G	100 g
L7651-25G	2!

### aseptically filled, Lyophilized powder

Prepared from L6876

-20°C

7773-50MG	50 mg

Enzymes

Φ

sigma-aldrich.com/biofil

### Lysozyme chloride form from chicken egg white

Mucopeptide N-acetylmuramoylhydrolase; Muramidase [9066-59-5] E.C. 3.2.1.17

### Grade VI, activity: ~60,000 units/mg protein

Enzyme which breaks down the cell walls of bacteria; used to prepare spheroplasts.

Lyophilized powder containing sodium chloride and sodium acetate

mol wt ~14.3 kDa

### Composition protein ~90%

### 3× Crystallized

One unit will produce a  $\Delta A_{450}$  of 0.001 per min at pH 6.24 at 25 °C, using a suspension of *Micrococcus lysodeikticus* as substrate, in a 2.6 mL reaction mixture (1 cm light path).

### -20°C

L2879-1G	1 g
L2879-5G	5 g
L2879-25G	25 g

### Lysozyme from human milk

Mucopeptide N-acetylmuramoylhydrolase; Muramidase [12671-19-1] E.C. 3.2.1.17

### lyophilized powder, activity: $\geq$ 100,000 units/mg protein

Lyophilized powder containing sodium phosphate and sodium chloride

### composition protein ~10%

One unit will produce a  $\Delta A_{450}$  of 0.001 per min at pH 6.24 at 25 °C, using a suspension of *Micrococcus lysodeikticus* as substrate, in a 2.6 mL reaction mixture (1 cm light path).

L6394-25KU	25,000 units
L6394-100KU	100,000 units

### Lysozyme from human neutrophils

Mucopeptide N-acetylmuramoylhydrolase; Muramidase [9001-63-2] E.C. 3.2.1.17

# ${\geq}95\%$ (SDS-PAGE), lyophilized powder, activity: ${\geq}100{,}000$ units/mg protein

Lyophilized from 50 mM sodium acetate, pH 6.0, with 100 mM  $\ensuremath{\mathsf{NaCl}}$ 

One unit will produce a  $\Delta A_{450}$  of 0.001 per min at pH 6.24 at 25 °C, using a suspension of *Micrococcus lysodeikticus* as substrate, in a 2.6 mL reaction mixture (1 cm light path).

-20°C

L8402-.1MG

0.1 mg

### **Mutanolysin**

# Mutanolysin from *Streptomyces globisporus* ATCC 21553

### [55466-22-3]

Mutanolysin is an N-acetylmuramidase. Like lysozyme, it is a muralytic enzyme that cleaves the  $\beta$ -N-acetylmuramyl-(1–4)-N-acetylglucosamine linkage of the bacterial cell wall polymer peptidoglycanpolysaccharide. Its carboxy terminal moieties are involved in the recognition and binding of unique cell wall polymers. Mutanolysin lyses *Listeria* and other Gram-positive bacteria such as *Lactobacillus* and *Lactococcus*.

Provides gentle cell lysis for the isolation of easily degradable biomolecules and RNA from bacteria. It has been used in the formation of spheroplasts for isolation of DNA.

### mol wt 23 kDa

One unit will produce a  $\Delta A_{600}$  of 0.01 per minute at pH 6.0 at 37 °C in a 1 mL volume using a suspension of *Streptococcus faecalis* cell wall as substrate.

### Iyophilized powder, activity:

### $\geq$ 4000 units/mg protein (biuret),

Chromatographically purified

Lyophilized powder containing  $\mathsf{Ficoll}^{\circledast}$  and sodium succinate buffer salts

-20°C

M9901-1KU	1,000 units
M9901-5KU	5,000 units
M9901-10KU	10,000 units
M9901-50KU	50,000 units

### ▶ aseptically filled, lyophilized powder, activity: ≥4000 units/mg protein (biuret)

Lyophilized powder containing  $\mathsf{Ficoll}^{\circledast}$  and sodium succinate buffer salts

Prepared from M9901

M4782-5KU

-20°C


5,000 units

SIGMA

### \_\_\_\_\_


### **Pectin Degrading Enzymes**

### **Pectinase and Pectinesterase Specificities**




Homogalacturonan


Poly- $\alpha$ -(1-4)-p-galacturonic acid backbone with random-partial methylation and acetylation



### Rhamnogalacturonan I

Alternating  $\alpha$ -(1-2)-L-rhamnosyl- $\alpha$ -(1-4)-D-galacturonosyl backbone with two types of branching composed of arabinofuranose or galactose oligomers





Pectins are complex branched heteropolysaccharides primarily containing an  $\alpha$ -(1-4) polygalacturonic acid backbone which can be randomly acetylated and methylated. Three different pectins have been isolated from plant cell walls. Homogalacturonans are composed of the simple  $\alpha$ -(1-4) polygalacturonic acid backbone. Substituted homogalacturonans are modifications of this backbone with  $\beta$ -D-xylose branching at C3, or apiofuranose substitutions in the backbone with  $\beta$ -D-Apiosyl-(1,3')- $\beta$ -D-Apiose branching. Rhamnogalacturonan I contains alternating  $\alpha$ -(1-4) galacturonosyl and  $\alpha$ -(1-2) rhamnosyl residues, with primarily oligo  $\alpha$ -(1-3) arabinose and oligo  $\beta$ -(1-4) galactose branching. Rhamnogalacturonan II is composed of the simple  $\alpha$ -(1-4) polygalacturonic acid backbone with complex branching composed of up to 11 different monosaccharide types.

14

SIGMA

A

### Pectinase

Pectinase catalyzes the random hydrolysis of  $\alpha$ -(1-4)-D-galactosiduronic linkages in pectin and other galacturonans.

### Pectinase from Aspergillus aculeatus

### Pectinex Ultra SPL®

### aqueous solution, activity: $\geq$ 26,000 units/mL

Highly active pectolytic enzyme preparation produced by a selected strain of *Aspergillus aculeatus* 

A product of Novozyme Corp.

2-8°C

P2611-50ML	50 mL
P2611-250ML	250 mL

### Pectinase from Aspergillus niger

### Pectinex® 3XL

### aqueous solution

Pectolytic enzyme preparation produced from a selected strain of *Aspergillus niger*: contains mainly pectintranseliminase, polygalacturonase, and pectinesterase and small amounts of hemicellulases and cellulases.

A product of Novozyme Corp.

2-8°C

P2736-50ML	50 mL
P2736-250ML	250 mL

### Pectinase solution from Aspergillus niger

Polygalacturonase solution from Aspergillus niger; Poly- $(1,4-\alpha-D-galacturonide)$  glycanohydrolase

[9032-75-1] E.C. 3.2.1.15

Used in plant protoplast preparation to digest cell wall prior to organelle isolation.

Solution in 40% glycerol

One unit will liberate 1.0  $\mu$ mole of galacturonic acid from polygalacturonic acid per min at pH 4.0 at 25 °C.

### aqueous glycerol solution, activity:

 $\geq$ 5 units/mg protein (Lowry)

2-8°0	

P4716-5KU	5,000 units	
P4716-10KU	10,000 units	
P4716-25KU	25,000 units	
P4716-100KU	100,000 units	

# $\blacktriangleright$ plant cell culture tested, aqueous glycerol solution, activity: ${\geq}5$ units/mg protein (Lowry)

### 2-8°C

P0690-10KU	10,000 units
P0690-25KU	25,000 units

### Pectinase from Rhizopus sp.

Macerozyme R-10; Polygalacturonase; Poly-(1,4- $\alpha$ -D-galacturonide) glycanohydrolase

[9032-75-1] E.C. 3.2.1.15 EC No. 2328856

Used in plant protoplast preparation to digest cell wall prior to organelle isolation.

One unit will liberate 1.0  $\mu$ mole of galacturonic acid from polygalacturonic acid per min at pH 4.0 at 25 °C.

### powder, activity: 400-800 units/g solid

Crude source of pectinase activity, also containing cellulase and hemicellulase activities.

-20°C

P2401-500UN	500 units
P2401-1KU	1,000 units
P2401-5KU	5,000 units

 plant cell culture tested, crude powder activity: 400-800 units/g solid

-20°C

P4300-1KU	1,000 units
P4300-5KU	5,000 units

### Visit the

# Enzyme Explorer Assay Library

Features over 600 detailed procedures for measuring enzyme activity and related metabolites. The Library is the result of over ten years of in-house process development by Sigma-Aldrich scientists.

sigma-aldrich.com/enzymeexplorer

# ma-aldrich.com/biofiles

SIGMA

s i g

### Pectin Degrading Enzymes

Pectinesterase catalyzes the hydrolysis of the methyl esters of pectin to yield pectate and methanol.

### Pectinesterase from orange peel

Pectin methylesterase; Pectin pectylhydrolase [9025-98-3] E.C. 3.1.1.11 EC No. 2328070

### lyophilized powder, activity: $\geq$ 150 units/mg protein Contains (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> and sodium chloride

Composition protein 20-50% (biuret)

One unit will release 1.0 microequivalent of acid from pectin per min at pH 7.5 at 30  $^\circ\text{C}.$ 

Protein determined by biuret.

### 2-8°C

P5400-1KU

### **Pectolyase**

Pectolyase catalyzes the eliminative cleavage of  $\alpha$ -(1- 4)-Dgalacturonan methyl ester to give oligosaccharides with 4-deoxy-6-O-methyl- $\alpha$ -D-galact-4-enuronosyl groups at their non-reducing ends.

### Pectolyase from Aspergillus japonicus

### E.C. 3.2.1.15

Reported to contain two types of pectinase, endopolygalacturonase (EC 3.2.1.15), endo-pectin lyase (EC 4.2.2.10) and a maceration stimulating factor.

Used in plant protoplast preparation to digest cell wall prior to organelle isolation.

Lyophilized powder containing lactose

One unit will liberate 1.0  $\mu mole$  of galacturonic acid from polygalacturonic acid per min at pH 5.5 at 25 °C.

### Iyophilized powder, activity: ≥0.3 units/mg solid

2-8°C

1,000 units

P3026-100MG	100 mg
P3026-250MG	250 mg
P3026-1G	1 g

▶ plant cell culture tested, lyophilized powder activity: ≥0.3 unit/mg solid

Composition protein ~60% (Lowry)

2-8°C	
P5936-100MG	100 mg
P5936-250MG	250 mg
P5936-1G	1 g

# The Enzyme Explorer

**The Protease Finder** 

Cleave proteins exactly where you want to with Sigma's new Protease Finder. The Protease Finder will identify the protease needed to cleave a specific peptide sequence at your desired location.

### Simple to use:


- Select either Endo- or Exoproteolytic cleavage.
- Enter your protein sequence into the positional boxes
   Submit the request to instantly receive the protease(s) capable of the cleavage

sigma-aldrich.com/proteasefinder



# sigma-aldrich.com/biofiles

### **Pullulanase Specificity**



### **Pullulan** Linear polymer of $\alpha$ -(1-6)-linked maltotriose units

Pullulanase catalyzes the hydrolysis of  $\alpha$ -(1-6)-D-glucosidic linkages in pullulan (a linear polymer of  $\alpha$ -(1-6)-linked maltotriose units), and, similar to isoamylase, in amylopectin and glycogen, and the  $\alpha$ - and  $\beta$ - limit dextrins of amylopectin and glycogen.

### Pullulanase from Bacillus acidopullulyticus

Pullulan 6-glucano-hydrolase [9075-68-7] E.C. 3.2.1.41 Promozyme<sup>®</sup> 400 L

### aqueous solution, $\geq$ 400 units/mL

Heat-stable debranching enzyme obtained from a selected strain of *Bacillus acidopullulyticus*, and belongs to the group of debranching enzymes known as pullulanases.

One unit is defined as the amount of enzyme which hydrolyzes pullulan, liberating reducing carbohydrate with a reducing power equivalent to 1.0  $\mu$ mole glucose per minute at pH 5.0 and 40 °C. density ...... 1.25 g/mL, 25 °C

A product of Novozymes Corp.

2 %

2-00	
P2986-50ML	50 mL
P2986-250ML	250 mL

### Pullulanase from Klebsiella pneumoniae

Amylopectin 6-gluconohydrolase; Limit dextrinase [9075-68-7] E.C. 3.2.1.41

One unit will liberate 1.0  $\mu mole$  of maltotriose (measured as glucose) from pullulan per min at pH 5.0 at 25 °C.

### Iyophilized powder, activity: 10-30 units/mg protein

Lyophilized powder containing potassium phosphate buffer salts and stabilizer

### Composition protein ~10% (Lowry)

-20°C	
P1067-100UN	
P1067-250UN	

ammonium sulfate suspension, activity:
 ≥5 units/mg protein (biuret)

Suspension in 3.2 M (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> solution, pH 6.2

Highly purified by a modification of ion exchange chromatography.

P5420-100UN	100 units
P5420-250UN	250 units

s i g

ma-aldrich.com/biofiles


Ú

IGMA

100 units

250 units

# Starch and Glycogen Degrading Enzymes



 Amylopectin and Glycogen

 Polymers of α-(1-4)-D-glucopyranosyl units. Amylopectin has approximately 4% α-(1-6)-D-glucopyranosyl branching.

 Glycogen tends to have a higher frequency of branches of shorter length than amylopectin.

### α-Amylase

 $\alpha$ -Amylase catalyzes the endohydrolysis of  $\alpha$ -(1-4)-D-glucosidic linkages in polysaccharides containing three or more  $\alpha$ -(1-4)-linked D-glucose units.

### α-Amylase from Bacillus licheniformis

1,4-α-D-Glucan-glucanohydrolase [9000-85-5] E.C. 3.2.1.1

# lyophilized powder, activity: 500-1,500 units/mg protein, ~95% (SDS-PAGE)

Lyophilized powder containing potassium

### Composition protein ~70% (BCA)

One unit will liberate 1.0 mg of maltose from starch in 3 min at pH 6.9 at 20 °C.

2-8°C	(etable)
A4551-100MG	100 mg
A4551-1G	1 g

### α-Amylase from Aspergillus oryzae

1,4-α-d-Glucan-glucanohydrolase

[9001-19-8] E.C. 3.2.1.1 EC No. 2325881

# Iyophilized powder, activity: 150-250 units/mg protein (biuret) Crude

One unit will liberate 1.0 mg of maltose from starch in 3 min at pH 6.9 at 20  $^\circ\text{C}.$ 

### -20°C

A6211-250KU	250,000 units
A6211-1MU	1,000,000 units
A6211-5MU	5,000,000 units

### ▶ Fungamyl<sup>®</sup> 800L aqueous solution, activity: ≥0.8 units/g

A product of Novozyme Corp.

Γ	2-8°C	1

A8220-50ML	50 mL
A8220-250ML	250 mL

### Taka-Diastase from Aspergillus oryzae

 $\alpha\text{-Amylase; }1,4\text{-}\alpha\text{-}D\text{-}Glucan-glucanohydrolase; Taka-Amylase A [9001-19-8] E.C. 3.2.1.1 EC No. 2325881$ 

### BioChemika, powder, activity: ~100 units/mg

One unit corresponds to the amount of enzyme which liberates 1  $\mu mol$  maltose per minute at pH 6.0 and 25 °C (starch according to Zulkowsky, Catalog No. 85642, as substrate).

2-8°	С

86247-25G	25 g
86247-100G	100 g

### α-Amylase from Bacillus sp.

1,4- $\alpha$ -D-Glucan-glucanohydrolase from *Bacillus* sp. [9000-90-2] E.C. 3.2.1.1

### ▶ powder, activity: ≥400 units/mg protein (Lowry) Contains starch as an extender.

One unit will liberate 1.0 mg of maltose from starch in 3 min at pH 6.9 at 20  $^\circ\text{C}.$ 

### -20°C

A6814-1MU	1,000,000 units
A6814-5MU	5,000,000 units
A6814-25MU	25,000,000 units

### ▶ Duramyl<sup>®</sup>, liquid, activity: ≥300 units/g

Protein-engineered  $\alpha$ -amylase produced by submerged fermentation of a genetically modified species of *Bacillus*. A product of Novozyme Corp.

2-6 0	
A7720-50ML	50 mL
A7720-250ML	250 mL

### Type II-A, lyophilized powder,

### activity: 1,500-3,000 units/mg protein (biuret)

mol wt 50-55 kDa by SDS-PAGE

4× crystallized

One unit will liberate 1.0 mg of maltose from starch in 3 min at pH 6.9 at 20  $^\circ\text{C}.$ 

-20°C

2 % C

A6380-25MG	25 mg
A6380-100MG	100 mg
A6380-250MG	250 mg
A6380-500MG	500 mg
A6380-1G	1 g

### α-Amylase from Bacillus amyloliquefaciens

BAN<sup>™</sup> 240L; 1,4-α-D-Glucan glucanohydrolase [9000-85-5] E.C. 3.2.1.1. EC No. 2325656

### liquid, activity: ≥250 units/g

This enzyme is active at high temperatures (70-90 °C).

mol wt 55 kDa

One unit is the amount of enzyme which dextrinizes 5.26 g dry starch per hour under standard conditions. A product of Novozyme Corp.

A7595-50ML	50 mL
A7595-250ML	250 mL

Ú

**G**MA

# Starch and Glycogen Degrading Enzymes

1,4-α-D-Glucan-glucanohydrolase [9000-85-5] E.C. 3.2.1.1 Termamyl<sup>®</sup> 120

### Type XII-A, saline solution, activity: 500-1,000 units/mg protein (biuret)

Aqueous solution containing approx. 15% sodium chloride and 25% sucrose

One unit will liberate 1.0 mg of maltose from starch in 3 min at pH 6.9 at 20 °C.

Reported to be heat stable at temperatures as high as ~90 °C. A product of Novozyme Corp.

### 2-8°C

Enzymes

A3403-500KU	500,000 units
A3403-1MU	1,000,000 units
A3403-5MU	5,000,000 units

### suitable for determination of starch (Kit STA-20)

### 2-8°C

A4582-.5ML

### α-Amylase, heat-stable

α-Amylase; 1,4-αD-Glucan-glucanohydrolase [9000-85-5] E.C. 3.2.1.1

### solution, For use in Total Dietary Fiber Assay, TDF-100A 2-8°C

A330	6-10M	L			10 mL
		~	 		

### α-Amylase from barley malt

 $1,4-\alpha$ -D-Glucan-glucanohydrolase [9000-90-2] E.C. 3.2.1.1

### Type VIII-A, powder

 $\alpha$ -amylase activity:  $\geq 1$  unit/mg solid  $\beta$ -amylase activity:  $\geq 1$  unit/mg solid contains lactose as standardization of activity Package size based on α-amylase activity One unit will liberate 1.0 mg of maltose from starch in 3 min at pH 6.9 at 20 °C.

### 2-8°C

A2771-10KU	10 KU
A2771-50KU	50 KU

### α-Amylase from human pancreas

[9000-90-2] E.C. 3.2.1.1 EC No. 2325656

### lyophilized powder, activity: ≥100 units/mg protein

Lyophilized from Tris buffer containing NaCl and CaCl<sub>2</sub>. purified by 3× crystallization

Prepared by modified method of Levitzki et al.

One unit will liberate 1.0 mg of maltose from starch in 3 min at pH 6.9 at 20 °C.

### -20°C

A9972-100UG

### α-Amylase from human saliva

1,4-α-D-Glucan-glucanohydrolase

### [9000-90-2] E.C. 3.2.1.1 EC No. 2325656

composition Protein ~10% (biuret)

One unit will liberate 1.0 mg of maltose from starch in 3 min at pH 6.9 at 20 °C.

### Type XIII-A, lyophilized powder, activity:

300-1,000 units/mg protein Lyophilized powder containing  $(NH_4)_2SO_4$  and sodium citrate. -20°C

A1031-1KU	1,000 units
A1031-5KU	5,000 units

### > Type IX-A, lyophilized powder, activity: 1,000-3,000 units/mg protein

Lyophilized powder containing (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> and sodium citrate Chromatographically purified

-20°C

0.5 mL

A0521-100UN	100 units
A0521-500UN	500 units
A0521-2.5KU	2,500 units

### $\alpha$ -Amylase from porcine pancreas

E.C. 3.2.1.1

Molecular Weight: 51-54 kDa.

 $\alpha$ -Amylase isolated from porcine pancreas is a glycoprotein. It is a single polypeptide chain of approximately 475 residues containing two SH groups and four disulfide bridges and a tightly bound Ca<sup>2+</sup> necessary for stability. Chloride ions are necessary for activity and stability. The pH range for activity is 5.5 to 8.0, with the pH optimum at 7.

One unit will liberate 1.0 mg of maltose from starch in 3 min at pH 6.9 at 20 °C.

### > Type VI-B, activity: 10-30 units/mg solid

Contains lactose

Package size based on  $\alpha$ -amylase activity

2-8°C

A3176-500KU	500,000 units
A3176-1MU	1,000,000 units
A3176-2.5MU	2,500,000 units
A3176-5MU	5,000,000 units
A3176-10MU	10,000,000 units

### Type I-A, DFP Treated, saline suspension, activity: 700-1400 units/mg protein

Suspension in 2.9 M NaCl solution containing 3 mM CaCl<sub>2</sub> DFP treated. 2× crystallized

2-8°C

100 µg

A6255-10MG	10 mg
A6255-25MG	25 mg
A6255-100MG	100 mg

Φ

### Type VII-A, DFP treated, ammonium sulfate suspension,

### activity: ≥500 units/mg protein

Suspension in 3.2 M  $(NH_4)_2SO_4$ , pH 6.1 2-8°C

A2643-10MG	10 mg
A2643-50MG	50 mg

# Type I-A, PMSF treated, saline suspension, activity: 700-1400 units/mg protein

Suspension in 2.9 M NaCl solution containing 3 mM CaCl\_2. 2× crystallized

2-8°C

A4268-25MG	25 mg
A4268-100MG	100 mg

### β-Amylase

 $\beta$ -Amylase catalyzes the exo-hydrolysis of  $\alpha$ -(1-4)-D-glucosidic linkages in polysaccharides resulting in the successive liberation of maltose units from the non-reducing ends of the chains.

### β-Amylase from barley

1,4-α-D-Glucan maltohydrolase [9000-91-3] E.C. 3.2.1.2 EC No. 2325661

### Type II-B, activity: 20-80 units/mg protein (biuret)

Crude

One unit will liberate 1.0 mg of maltose from starch in 3 min at pH 4.8 at 20  $^\circ\text{C}.$ 

### 2-8°C

A7130-10KU	10,000 units
A7130-50KU	50,000 units
A7130-250KU	250,000 units

### β-Amylase from sweet potato

1,4-α-D-Glucan maltohydrolase [9000-91-3] E.C. 3.2.1.2 EC No. 2325661

# Type I-B, ammonium sulfate suspension, activity: 750-1,000 units/mg protein

Crystalline suspension in 2.3 M (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> One unit will liberate 1.0 mg of maltose from starch in 3 min at pH 4.8 at 20 °C.

### 2-8°C

A7005-10KU	10,000 units
A7005-25KU	25,000 units
A7005-50KU	50,000 units
A7005-100KU	100,000 units

### Amyloglucosidase

Amyloglucosidase catalyzes the hydrolysis of terminal  $\alpha$ -(1-4)linked D-glucose residues successively from the non-reducing ends of maltooligo-and polysaccharides with release of  $\beta$ -D-glucose. Most forms of the enzyme can rapidly hydrolyze  $\alpha$ -(1-6)-Dglucosidic bonds when the next bond in the sequence is 1,4- and some preparations of this enzyme hydrolyze 1,6- and  $\alpha$ -(1-3)-Dglucosidic bonds in other polysaccharides.

### Amyloglucosidase from Aspergillus niger

Exo-1,4- $\alpha$ -glucosidase; 1,4- $\alpha$ -D-Glucan glucohydrolase; Glucoamylase [9032-08-0] E.C. 3.2.1.3 EC No. 2328772

### ▶ lyophilized powder, activity: ≥80 units/mg protein (biuret)

Lyophilized powder containing less than 0.02% glucose One unit will liberate 1.0 mg of glucose from starch in 3 min at pH 4.5 at 55  $^\circ C.$ 

### -20°C

A7420-5MG	5 mg
A7420-25MG	25 mg
A7420-100MG	100 mg

### ▶ AMG<sup>™</sup> 300L, ≥300 units/mL

A7095-50ML 50 ml	L
2-8°C	_
A product of Novozymes Corp.	
density ~1.2 g/mL, 25 °	С
aqueous solution	
Stabilized with glucose	
Stabilized with alucose	

### ▶ aqueous glucose solution activity: ≥5000 units/mL

Solution in 1 M glucose containing 0.5% sodium benzoate as preservative

One unit will liberate 1.0 mg of glucose from starch in 3 min at pH 4.5 at 55  $^\circ\text{C}.$ 

2-8°C	
A3042-50ML	50 mL

### > ammonium sulfate suspension, activity:

 $\geq$ 40 units/mg protein Suspension in 3.2 M (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> solution, pH approx. 6.0

One unit will liberate 1.0 mg of glucose from starch in 3 min at pH 4.5 at 55  $^\circ\text{C}.$ 

# 2-&\*C A1602-25MG 25 mg A1602-100MG 100 mg

Ú

**G**MA

g

# Starch and Glycogen Degrading Enzymes

### Amyloglucosidase from Candida tsukubaensis

Exo-1,4- $\alpha$ -glucosidase; 1,4- $\alpha$ -p-Glucan glucohydrolase; Glucoamylase [9032-08-0] E.C. 3.2.1.3

# ammonium sulfate suspension, activity: 50-150 units/mg protein

An acid-stable amyloglucosidase, maintaining high activity at pH values down to 2.5.

Suspension in 3.2 M (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>, pH approx. 5.5.

One unit will liberate 1.0 mg of glucose from starch in 3 min at pH 4.5 at 55  $^\circ\mathrm{C}.$ 

### 2-8°C

A2330-25UN	25 units
A2330-100UN	100 units

### Amyloglucosidase from *Rhizopus* sp.

Exo-1,4- $\alpha$ -glucosidase; 1,4- $\alpha$ -D-Glucan glucohydrolase; Glucoamylase [9032-08-0] E.C. 3.2.1.3 EC No. 2328772

### activity: ≥40,000 units/g solid Lyophilized salt-free powder

### Composition protein ~35% (biuret)

One unit will liberate 1.0 mg of glucose from starch in 3 min at pH 4.5 at 55  $^\circ\mathrm{C}.$ 

# A9228-1G 1g

### α-Glucosidase

 $\alpha$ -Glucosidase catalyzes the hydrolysis of terminal  $\alpha$ -(1-4)-linked D-glucose residues successively from the non-reducing ends of maltooligo-and to a lesser extent polysaccharides with release of  $\beta$ -D-glucose. Most forms of the enzyme can slowly hydrolyze  $\alpha$ -(1-6)-D-glucosidic bonds.

### α-Glucosidase from Bacillus stearothermophilus

 $\alpha\text{-D-Glucoside glucohydrolase; Maltase}$  [9001-42-7] E.C. 3.2.1.20

### lyophilized powder, activity: ≥50 units/mg protein

Lyophilized powder containing potassium phosphate buffer salt

One unit will liberate 1.0 µmole of D-glucose from p-nitrophenyl  $\alpha\text{-}D\text{-}glucoside$  per min at pH 6.8 at 37 °C.

Protein determined by biuret.

 $\beta\text{-Glucosidase},\,\alpha\text{-galactosidase}$  and  $\beta\text{-galactosidase}$  ......  ${\leq}0.1\%$  [2-8°C]

G3651-250UN

250 units

### $\alpha$ -Glucosidase from rice

 $\alpha\text{-}\text{D-}Glucoside glucohydrolase; Maltase [9001-42-7] E.C. 3.2.1.20$ 

# Type V, ammonium sulfate suspension, activity: 40-80 units/mg protein

Suspension in 2.8 M  $(NH_4)_2SO_4$  solution One unit will convert 1.0 µmole of maltose to 2.0 µmoles of p-glucose per min at pH 4.0 at 37 °C. Protein determined by biuret. [2-8°C]

100 units

G9259-100UN

### α-Glucosidase from Saccharomyces cerevisiae

 $\alpha\text{-}\text{D-}\text{Glucosidase};$   $\alpha\text{-}\text{D-}\text{Glucoside}$  glucohydrolase; Maltase from yeast [9001-42-7] E.C. 3.2.1.20 EC No. 2326047

For the determination of  $\alpha\text{-}amylase$  and the synthesis of various 1'-osucrose and 1-O-fructose esters

Protein determined by biuret.

# ▶ recombinant, expressed in unspecified host, lyophilized powder, activity: ≥125 units/mg protein

Lyophilized powder containing potassium phosphate buffer salt pH 7.15 and approx. 70% lactose

One unit will liberate 1.0  $\mu mole$  of D-glucose from p-nitrophenyl  $\alpha\text{-}Dglucoside$ 

per min at pH 6.8 at 37 °C.

	2-8°C
--	-------

G0660-750UN	750 units

# ▶ Type I, lyophilized powder, activity: $\geq$ 10 units/mg protein (using p-nitrophenyl $\alpha$ -D-glucoside as substrate.)

contains phosphate buffer salts and EDTA as balance composition

Protein ~50%

Sold on basis of p-nitrophenyl  $\alpha$ -p-glucoside units.

One unit will liberate 1.0 µmole of D-glucose from p-nitrophenyl  $\alpha$ -D-glucoside per min at pH 6.8 at 37 °C.

G5003-100UN	100 units
G5003-1KU	1,000 units

SIGMA's sigma-aldrich.com/biofile

# Kits for Carbohydrate Analysis

Sigma manufactures several unique enzymatic-based kits for the quantitation of important carbohydrates. These kits utilize spectrophotometric, and gravimetric detection making them easyto-use, yielding high sensitivity, and consistent results.

### Total Dietary Fiber Assay Kit, Cat. No. TDF100A-1KT

### sufficient for ~100 assays

For the determination of total dietary fiber. Uses a combination of enzymatic and gravimetric methods to analyze samples of dried, defatted foods to determine soluble fiber, protein, and ash content. This procedure is based on the method published by AOAC.<sup>1</sup>

### Reference:

<sup>1</sup>Official Methods of Analysis, 16th ed., AOAC, Arlington, VA, Vol. II, Sec. 45.4.07, Method 985.29, 1105 (1997).

### **Total Dietary Fiber Assay Procedure**

Heat stable α-Amylase, incubation at pH 6.0, 15 min., 95 °C

### $\downarrow$

Protease incubation at pH 7.5, 30 min., 60 °C

 $\downarrow$ 

Amyloglucosidase incubation at pH 4.5, 30 min., 60 °C





### Dietary Fiber, Total, Assay Control Kit, Cat. No. TDFC10-1KT


sufficient for ~10 assays


Set of 6 standards for use as internal controls in conjunction with the Total Dietary Fiber Assay Kit (TDF100A)

### Starch (GO/P) Assay Kit, Cat. No. STA20-1KT

### sufficient for 20 assays

For the quantitative, enzymatic determination of starch in food and other materials. The hydrolysis of starch to glucose is catalyzed by  $\alpha$ -amylase and amyloglucosidase. Glucose is oxidized to gluconic acid and hydrogen peroxide by glucose oxidase. Hydrogen peroxide reacts with *o*-dianisidine in the presence of peroxidase to form a colored product. Oxidized *o*-dianisidine reacts with sulfuric acid to form a more stable colored product. The intensity of the pink color measured at 540 nm is proportional to the original glucose concentration.





Visit the Enzyme Explorer Protease Inhibitor Guide

The on-line guide to broad spectrum and high specificity inhibition strategies, reagents and cocktails

sigma-aldrich.com/enzymeexplorer


Ú

IGMA

### Starch (HK) Assay Kit, Cat. No. SA20-1KT

### sufficient for 20 assays


For the quantitative, enzymatic determination of native starch in food and other materials. The hydrolysis of starch to glucose is catalyzed by amyloglucosidase. Glucose is phosphorylated by hexokinase. Glucose-6-phosphate is then oxidized to 6-phosphogluconate in the presence of NAD in a reaction catalyzed by glucose 6-phosphate dehydrogenase. The increase in absorbance at 340 nm is directly proportional to the glucose concentration.



### Glucose (GO) Assay Kit, Cat. No. GAGO20-1KT

### sufficient for 20 assays


For the quantitative, enzymatic determination of glucose in food and other materials. Glucose is oxidized to gluconic acid and hydrogen peroxide by glucose oxidase. Hydrogen peroxide reacts with *o*-dianisidine in the presence of peroxidase to form a colored product. Oxidized *o*-dianisidine reacts with sulfuric acid to form a more stable colored product. The intensity of the pink color measured at 540 nm is proportional to the original glucose concentration.



### Glucose (HK) Assay Kit, Cat. No. GAHK20-1KT

### sufficient for 20 assays


For the quantitative, enzymatic determination of glucose in food and other materials. Glucose is phosphorylated by hexokinase to form glucose 6-phosphate. Glucose 6-phosphate is then oxidized to 6-phospho gluconate in the presence of NAD by glucose 6-phosphate dehydrogenase. During this oxidation, an equimolar amount of NAD is reduced to NADH. The consequent increase in absorbance at 340 nm is directly proportional to glucose concentration.



### Fructose Assay Kit, Cat. No. FA20-1KT

### sufficient for 20 assays

For the quantitative, enzymatic determination of fructose in food and other materials. Fructose is phosphorylated by ATP using hexokinase. Fructose 6-phosphate is then converted to glucose 6-phosphate by phosphoglucose isomerase. Glucose 6-phosphate is then oxidized to 6-phosphogluconate in the presence of NAD by glucose 6-phosphate dehydrogenase. During this oxidation, an equimolar amount of NAD is reduced to NADH. The consequent increase in absorbance at 340 nm is directly proportional to fructose concentration.



Kits

SIGMA

The following table contains selected polysaccharides, glycosaminoglycans, and related reagents for use as standards, glycan detection and measurement of enzymatic activity.

Product Name	Description	Cat. No.
Alginic acid ammonium calcium salt Macrocystis pyrifera (kelp)	-	A7253-100G
Alginic acid from brown algae		A7003-100G
		A7003-250G
		A7003-500G
		A7003-1KG
Alginic acid sodium from salt brown algae	for immobilization of microorganisms	71238-50G
5		71238-250G
		71238-1KG
Alginic acid sodium from salt brown algae	Low viscosity	A2158-100G
		A2158-250G
		A2158-500G
		A2158-1KG
Alginic acid sodium from salt brown algae	Medium viscosity	A2033-100G
ignio dola sociati non salt storm digao	includin historicy	A2033-250G
		A2033-2500G
		A2033-300G
Nalais said sadium from salt brown alass	Lourvisessity elect call autors tested, pourder	
Alginic acid sodium from salt brown algae	Low viscosity, plant cell culture tested, powder	A0682-100G
Amylopectin potato starch	•	A8515-25G
Annula a shin. Anuna	Average Calestone	A8515-100G
Amylopectin Azure	Amylase Substrate	A4640-1G
		A4640-5G
		A4640-25G
		A4640-50G
Amylose from potato	Essentially Free of Amylopectin	A0512-250MG
		A0512-1G
		A0512-5G
		A0512-25G
Amylose–Remazol Brilliant Blue R	Amylase Substrates	A3508-1G
		A3508-5G
Chitin from crab shells	suitable for analysis of chitinase, purified powder	C9752-250MG
		C9752-1G
		C9752-5G
Chitin from crab shells	practical grade, powder	C7170-100G
		C7170-1KG
Chitin from crab shells	practical grade, coarse flakes	C9213-500G
		C9213-1KG
Chitin azure	Chitinase Substrate	C3020-100MG
		C3020-1G
Chitosan from crab shells	≥75% deacetylated	C3646-10G
	-	C3646-25G
		C3646-100G
		C3646-500G
Chondroitin disaccharide ∆di-OS sodium salt		C3920-5MG
		C3920-10MG
Chondroitin disaccharide ∆di-4S sodium salt		C4045-5MG
Shenarettin diadeendride zur Folgodiditt Bait		C4045-10MG
Chondroitin disaccharide ∆di-6S sodium salt		C4045-10MG C4170-5MG
enonaronan disacchande Zur-os soulditt salt	-	
		C4170-25MG
Chondroitin disaccharide ∆di-UA-2S sodium salt	-	C5820-1MG
Chondroitin 6-sulfate sodium salt from shark cartilage	~90%, balance is chondroitin sulfate A	C4384-250MG
		C4384-1G
		C4384-5G
		C4384-25G
Chondroitin sulfate B sodium salt	from porcine intestinal mucosa, $\geq$ 90%, lyophilized powder	C3788-25MG
		C3788-100MG
Nonoclonal Anti-Chondroitin Sulfate antibody produced in mouse	clone CS-56, ascites fluid	C80352ML
		C80355ML
Chondroitin sulfate A sodium salt from bovine trachea	cell culture tested	C9819-5G
		C9819-25G

sigma-aldrich.com/biofiles

**SIGMA** 

Φ
—
4
0
۵
~
Ε
0
U
•
ے
U
5
σ
-
a
a
Ε
σ
S

Ś

Product Name	Description	Cat. No.
Curdlan from Alcaligenes faecalis		C7821-5G
Dammar Resin	-	30424-250G
		30424-1KG
Dextran from Leuconostoc mesenteroides	average mol wt 9,000-11,000	D9260-10G
		D9260-50G
		D9260-100G
		D9260-500G
Dextran from Leuconostoc mesenteroides	average mol wt 64,000-76,000	D4751-10G
		D4751-50G
		D4751-100G
		D4751-500G
		D4751-1KG
Dextran from Leuconostoc mesenteroides	average mol wt 35,000-45,000	D1662-10G
		D1662-50G
		D1662-100G
		D1662-500G
Dextran from Leuconostoc mesenteroides	average mol wt 425,000-575,000	D1037-50G
		D1037-100G
		D1037-500G
Dextran Leuconostoc mesenteroides	average mol wt 100,000-200,000	D1037-500G
Sextrait Leuconosioc mesenteroides	average mor wit 100,000-200,000	D4876-30G
		D4876-500G
		D4876-1KG
Dextran from Leuconostoc mesenteroides	industrial grade, average mol wt 5,000,000-40,000,000	D5501-100G
		D5501-500G
		D5501-1KG
Dextran from Leuconostoc mesenteroides	average mol wt ~2,000,000	D5376-100G
		D5376-500G
Dextran solution from Leuconostoc mesenteroides	20 % (w/w) (Autoclaved)	D8802-25ML
		D8802-50ML
3-D-Glucan from barley	powder	G6513-50M0
-		G6513-100M
		G6513-500N
		G6513-1G
		G6513-5G
Glucan from baker's yeast (S. cerevisiae)		G5011-25M0
Sidean non bakers yeast (s. eerewside)		G5011-100N
Glycogen from bovine liver		G0885-1G
Siycogen noni bovine liver	-	G0885-1G G0885-5G
		G0885-10G
Norman Organista la Compacta (alla de la del		G0885-25G
Slycogen Crepidula fornicata (slipper limpet)	•	G1633-5G
Glycogen from Mytilus edulis (Blue mussel)		G1508-5G
		G1508-25G
Slycogen from oyster		G1765-5MG
		G1765-10M0
		G1765-25M0
Glycogen from oyster		G8751-5G
		G8751-25G
		G8751-100G
Glycogen from rabbit liver		G8876-500N
		G8876-1G
		G8876-5G
		G8876-10G
Glycogen azure	from rabbit liver suitable for substrate for a amulase	G5510-1G
Siycol chitosan	from rabbit liver, suitable for substrate for $\alpha$ -amylase	
	$\geq$ 60% (colloidal titration), crystalline	G7753-500M
Sybor of itestant		G7753-1G

Product Name	Description	Cat. No.
Guar		G4129-250G
		G4129-500G
		G4129-1KG
Gum accroides		G9627-500G
Gum arabic from acacia tree	reagent grade	G9752-500G
		G9752-1KG
Heparan sulfate sodium salt from bovine kidney		H7640-1MG
		H7640-5MG
		H7640-10MG
Heparan sulfate fast-moving fraction sodium salt from porcine		H9902-1MG
intestinal mucosa	$\geq$ 90% (electrophoresis)	H9902-5MG
Heparan sulfate proteoglycan	≥400 µg/mL glycosaminoglycan, sterile-filtered	H47771MG
Heparin sodium salt from porcine intestinal mucosa	mol wt 4,000-6,000 Da	H8537-50MG
		H8537-100MG
		H8537-250MG
		H8537-1G
Heparin sodium salt from porcine intestinal mucosa	average mol wt ~3,000	H3400-50MG
		H3400-100MG
		H3400-250MG
		H3400-250MG
Heparin sodium salt from porcine intestinal mucosa	Grade I-A, activity: ~170 USP units/mg	H3393-10KU
hepann souidh san nom porche intestinar nucosa	Glade PA, activity. ~170 03F drifts/ring	
		H3393-25KU
		H3393-50KU H3393-100KU
		H3393-250KU
		H3393-500KU
		H3393-1MU
Heparin-benzalkonium	activity: ~60 USP units/mg	H7280-1G
Heparin, deaminated sodium salt from porcine intestinal mucosa	Low molecular weight mono-aldehyde, heparin activity: >75 USP units/mg	H7405-250MG
		H7405-1G
Heparin disaccharide I-A sodium salt	α-δUA-2S-[1→4]-GlcNAc-6S	H9517-1MG
Heparin disaccharide I-H sodium salt	α-δUA-2S-[1→4]-Glc-6S	H8892-1MG
Heparin disaccharide I-S sodium salt	$(\alpha - \Delta UA - 2S - [1 \rightarrow 4] - GlcNS - 6S)$	H9267-1MG
Heparin disaccharide II-H sodium salt	$(\alpha \cdot \delta \cup A \cdot [1 \rightarrow 4] \cdot GlcN \cdot 6S)$	H9017-1MG
Heparin disaccharide III-H sodium salt	$(\alpha \cdot \delta \cup A - S2 - [1 \rightarrow 4] - GlcN)$	H9142-1MG
Heparin disaccharide III-S sodium salt	$(\alpha - \Delta UA - 2S - [1 \rightarrow 4] - GlcNS)$	H9392-1MG
Heparin disaccharide IV-A sodium salt	$(\alpha - \Delta UA - [1 \rightarrow 4] - GlcNAc)$	H08955MG
Heparin disaccharide IV-H ≥95%	$\alpha$ - $\Delta$ UA-[1 $\rightarrow$ 4]-GlcN	H9276-1MG
Hyaluronan biotin sodium salt	>97%, soluble powder	B1557-5MG
Hyaluronic acid potassium salt from human umbilical cord	suitable as substrate for hyaluronidase	H1504-50MG
		H1504-100MG
		H1504-500MG
		H1504-1G
Hyaluronic acid potassium salt from human umbilical cord	Highly polymerized	H1751-500MG
Hyaluronic acid sodium salt from bovine vitreous humor		H7630-10MG
		H7630-50MG
Hyaluronic acid sodium salt from rooster comb		H5388-100MG
		H5388-250MG
		H5388-1G
Hyaluronic acid sodium salt from Streptococcus equi		53747-1G
. ,		53747-10G
Hyaluronic acid sodium salt from Streptococcus zooepidemicus		H9390-1G
Hyaluronic acid disaccharide $\delta$ DiHA sodium salt	≥95%	H9649-1MG
Inulin from chicory		I2255-10G
nain non chicory		12255-10G
		I2255-100G
		I2255-1KG

sigma-aldrich.com/biofiles

SIGMA.

Product Name	Description	Cat. No.
nulin from dahlia tubers	Mr <sup>2</sup> 5000	13754-25G
		I3754-100G
		13754-1666
nulin_FITC	from dahlia tuber	F3272-1G
Lichenan from Cetraria islandica	practical grade, powder	L6133-250MG
		L6133-1G
Pectin from apple	meets USP testing specifications	P8471-100G
		P8471-500G
Pectin from citrus peel	Galacturonic acid: ≥74.0%	P9135-100G
		P9135-500G
		P9135-1KG
Pectin, esterified from citrus fruit	extent of labeling: ~90% esterified	P9561-5G
		P9561-25G
Pectin, esterified potassium salt from citrus fruit	extent of labeling: ~60% esterified	P9436-5G
		P9436-25G
		P9436-50G
Pectin, esterified potassium salt from citrus fruit	extent of labeling: ~30% esterified	P9311-5G
ectin, esterned potassium sait nom citrus nuit	extent of labeling. ~3078 estenned	
		P9311-25G
Peptidoglycan from Bacillus subtilis	•	69554-10MG
Peptidoglycan from Micrococcus luteus	•	53243-10MG
Peptidoglycan from Saccharomyces cerevisiae	• ·	72789-10MG
Peptidoglycan from Staphylococcus aureus		77140-10MG
		77140-25MG
Peptidoglycan from Streptomyces sp.		79682-10MG
Pullulan from Aureobasidium pullulans	suitable for substrate for pullulanase	P4516-1G
		P4516-5G
		P4516-25G
Stachvosa hydrate from Stachve tuharifara	>02%	
Stachyose hydrate from Stachys tuberifera	≥98%	S4001-10MG
Stachyose hydrate from Stachys tuberifera	≥98%	S4001-10MG S4001-100M0
Stachyose hydrate from Stachys tuberifera	≥98%	S4001-10MG S4001-100M0 S4001-500M0
Stachyose hydrate from <i>Stachys tuberifera</i>	≥98%	S4001-10MG S4001-100M0 S4001-500M0 S4001-1G
		S4001-10MG S4001-100M0 S4001-500M0 S4001-1G S4001-5G
Stachyose hydrate from <i>Stachys tuberifera</i>	≥98%	S4001-10MG S4001-100M0 S4001-500M0 S4001-1G
		S4001-10MG S4001-100M0 S4001-500M0 S4001-1G S4001-5G
Starch		S4001-10MG S4001-100M0 S4001-500M0 S4001-1G S4001-5G 85645-100G
Starch	from potatoes, for electrophoresis	S4001-10MG S4001-100M0 S4001-500M0 S4001-1G S4001-5G 85645-100G 85645-1KG
Starch Starch Starch from corn	from potatoes, for electrophoresis from potatoes	S4001-10MG S4001-100M0 S4001-500M0 S4001-1G S4001-5G 85645-100G 85645-1KG 85650-1KG
Starch Starch Starch from corn	from potatoes, for electrophoresis from potatoes Standard for Starch Assay Kits SA-20 and STA-20	\$4001-10MG \$4001-100M0 \$4001-500M0 \$4001-5G \$5645-100G 85645-1KG 85650-1KG \$5296-5G
Starch Starch Starch from corn	from potatoes, for electrophoresis from potatoes Standard for Starch Assay Kits SA-20 and STA-20	\$4001-10MG \$4001-100Md \$4001-500Md \$4001-5G \$4001-5G \$5645-100G \$5645-1KG \$5650-1KG \$5296-5G \$4180-100G
Starch Starch Starch from corn Starch from corn	from potatoes, for electrophoresis from potatoes Standard for Starch Assay Kits SA-20 and STA-20 practical grade	\$4001-10MG \$4001-100MG \$4001-500MG \$4001-5G 85645-100G 85645-1KG 85650-1KG \$5296-5G \$4180-100G \$4180-500G \$4180-1KG
Starch Starch Starch from corn	from potatoes, for electrophoresis from potatoes Standard for Starch Assay Kits SA-20 and STA-20	\$4001-10MG \$4001-500M0 \$4001-500M0 \$4001-5G \$5645-100G \$5645-1KG \$55296-5G \$4180-100G \$4180-500G \$4180-500G \$4180-1KG \$4126-2KG
Starch Starch from corn Starch from corn	from potatoes, for electrophoresis from potatoes Standard for Starch Assay Kits SA-20 and STA-20 practical grade Unmodified regular corn starch containing approx. 73% amylopectin and 27% amylose	\$4001-10MG \$4001-500M0 \$4001-500M0 \$4001-5G 85645-100G 85645-10G 85645-1KG \$5296-5G \$4180-100G \$4180-500G \$4180-500G \$4126-2KG \$4126-5KG
Starch Starch from corn Starch from corn	from potatoes, for electrophoresis from potatoes Standard for Starch Assay Kits SA-20 and STA-20 practical grade	\$4001-10MG \$4001-500M0 \$4001-500M0 \$4001-5G 85645-100G 85645-10G 85645-1KG \$5296-5G \$4180-100G \$4180-500G \$4180-1KG \$4126-2KG \$4126-5KG \$5651-100G
Starch Starch from corn Starch from corn	from potatoes, for electrophoresis from potatoes Standard for Starch Assay Kits SA-20 and STA-20 practical grade Unmodified regular corn starch containing approx. 73% amylopectin and 27% amylose	\$4001-10MG \$4001-500M0 \$4001-500M0 \$4001-5G 85645-100G 85645-10G 85645-1KG \$5296-5G \$4180-100G \$4180-500G \$4180-500G \$4126-2KG \$4126-5KG \$5651-100G
itarch itarch itarch from corn itarch from corn	from potatoes, for electrophoresis from potatoes Standard for Starch Assay Kits SA-20 and STA-20 practical grade Unmodified regular corn starch containing approx. 73% amylopectin and 27% amylose	\$4001-10MG \$4001-500M0 \$4001-500M0 \$4001-5G 85645-100G 85645-1KG 85650-1KG \$4180-500G \$4180-500G \$4180-500G \$4126-2KG \$4126-5KG \$5651-100G \$5651-500G
Starch Starch from corn Starch from corn	from potatoes, for electrophoresis from potatoes Standard for Starch Assay Kits SA-20 and STA-20 practical grade Unmodified regular corn starch containing approx. 73% amylopectin and 27% amylose	\$4001-10MG \$4001-100M \$4001-500M \$4001-5G \$5645-100G \$5645-10G \$5650-1KG \$5296-5G \$4180-100G \$4180-500G \$4180-500G \$4126-2KG \$4126-5KG \$5651-100G
Starch Starch from corn Starch from corn	from potatoes, for electrophoresis from potatoes Standard for Starch Assay Kits SA-20 and STA-20 practical grade Unmodified regular corn starch containing approx. 73% amylopectin and 27% amylose	\$4001-10MG \$4001-500M0 \$4001-500M0 \$4001-5G 85645-100G 85645-1KG 85645-1KG \$5296-5G \$4180-100G \$4180-500G \$4180-1KG \$4126-2KG \$4126-5KG \$5651-100G \$5651-500G
Starch Starch from corn Starch from corn Starch from corn Starch from potato	from potatoes, for electrophoresis from potatoes Standard for Starch Assay Kits SA-20 and STA-20 practical grade Unmodified regular corn starch containing approx. 73% amylopectin and 27% amylose	\$4001-10MG \$4001-500M0 \$4001-500M0 \$4001-5G 85645-100G 85645-1KG 85650-1KG \$4180-500G \$4180-500G \$4180-500G \$4126-2KG \$4126-5KG \$5651-100G \$5651-500G \$5651-2KG
Starch Starch from corn Starch from corn Starch from corn Starch from potato	from potatoes, for electrophoresis from potatoes Standard for Starch Assay Kits SA-20 and STA-20 practical grade Unmodified regular corn starch containing approx. 73% amylopectin and 27% amylose for electrophoresis	\$4001-10MG \$4001-100MG \$4001-50 85645-100G 85645-1KG 85650-1KG \$5296-5G \$4180-100G \$4180-500G \$4180-1KG \$4126-2KG \$4126-2KG \$5651-100G \$5651-500G \$5651-2KG
itarch itarch itarch from corn itarch from corn itarch from corn	from potatoes, for electrophoresis from potatoes Standard for Starch Assay Kits SA-20 and STA-20 practical grade Unmodified regular corn starch containing approx. 73% amylopectin and 27% amylose for electrophoresis	\$4001-10MG \$4001-100M \$4001-50 85645-100G 85645-1KG 85645-1KG \$5296-5G \$4180-100G \$4180-100G \$4180-500G \$4180-1KG \$4126-2KG \$4126-5KG \$5651-100G \$5651-1KG \$5651-2KG \$5651-5KG \$2630-100G
itarch itarch from corn itarch from corn itarch from corn itarch from potato itarch from potato	from potatoes, for electrophoresis from potatoes Standard for Starch Assay Kits SA-20 and STA-20 practical grade Unmodified regular corn starch containing approx. 73% amylopectin and 27% amylose for electrophoresis Soluble	\$4001-10MG \$4001-100M \$4001-50 85645-100G 85645-1KG 85645-1KG \$5296-5G \$4180-100G \$4180-100G \$4180-1KG \$4126-2KG \$4126-2KG \$5651-100G \$5651-2KG \$5651-2KG \$5651-5KG \$2630-100G \$2630-500G \$2630-1KG
itarch itarch from corn itarch from corn itarch from corn itarch from potato itarch from potato	from potatoes, for electrophoresis from potatoes Standard for Starch Assay Kits SA-20 and STA-20 practical grade Unmodified regular corn starch containing approx. 73% amylopectin and 27% amylose for electrophoresis	\$4001-10MG \$4001-100M \$4001-50 85645-100G 85645-1KG 85645-1KG \$5296-5G \$4180-100G \$4180-100G \$4180-1KG \$4126-2KG \$4126-2KG \$5651-100G \$5651-1KG \$5651-2KG \$5651-2KG \$2630-100G \$2630-100G \$2630-1KG \$2004-500G
itarch itarch from corn itarch from corn itarch from corn itarch from potato itarch from potato itarch from potato	from potatoes, for electrophoresis         from potatoes         Standard for Starch Assay Kits SA-20 and STA-20         practical grade         Unmodified regular corn starch containing approx. 73% amylopectin and 27% amylose         for electrophoresis         Soluble         Soluble	\$4001-10MG \$4001-100M \$4001-50 85645-100G 85645-1KG 85645-1KG \$5296-5G \$4180-100G \$4180-100G \$4180-1KG \$4126-2KG \$4126-2KG \$5651-100G \$5651-2KG \$5651-2KG \$5651-5KG \$2630-100G \$2630-1KG \$2004-500G \$2004-1KG
itarch itarch from corn itarch from corn itarch from corn itarch from potato itarch from potato itarch from potato	from potatoes, for electrophoresis from potatoes Standard for Starch Assay Kits SA-20 and STA-20 practical grade Unmodified regular corn starch containing approx. 73% amylopectin and 27% amylose for electrophoresis Soluble	\$4001-10MG \$4001-100M \$4001-50 85645-100G 85645-1KG 85645-1KG 85650-1KG \$5296-5G \$4180-100G \$4180-100G \$4180-1KG \$4126-2KG \$4126-2KG \$5651-100G \$5651-2KG \$5651-2KG \$5651-5KG \$2630-100G \$2630-1KG \$2004-500G \$2004-1KG \$4251-2KG
Starch Starch from corn Starch from corn Starch from corn Starch from potato Starch from potato Starch from potato Starch from potato Starch from potato	from potatoes, for electrophoresis         from potatoes         Standard for Starch Assay Kits SA-20 and STA-20         practical grade         Unmodified regular corn starch containing approx. 73% amylopectin and 27% amylose         for electrophoresis         Soluble         Soluble         Powder	S4001-10MG S4001-100M S4001-50 B5645-100G B5645-1KG B5645-1KG S5296-5G S4180-100G S4180-100G S4180-1KG S4126-2KG S4126-2KG S5651-100G S5651-2KG S5651-2KG S5651-5KG S2630-100G S2630-1KG S2630-1KG S2004-500G S2004-1KG S4251-2KG
Starch Starch from corn Starch from corn Starch from corn Starch from potato Starch from potato Starch from potato Starch from potato Starch from potato	from potatoes, for electrophoresis         from potatoes         Standard for Starch Assay Kits SA-20 and STA-20         practical grade         Unmodified regular corn starch containing approx. 73% amylopectin and 27% amylose         for electrophoresis         Soluble         Soluble	\$4001-10MG \$4001-100MG \$4001-50 85645-100G 85645-1KG 85645-1KG 85650-1KG \$5296-5G \$4180-100G \$4180-100G \$4180-100G \$4180-1KG \$4126-2KG \$4126-2KG \$5651-10KG \$5651-2KG \$5651-2KG \$2630-100G \$2630-100G \$2630-1KG \$2004-500G \$2004-1KG \$4251-2KG \$4251-2KG \$4251-2KG
Starch Starch from corn Starch from corn	from potatoes, for electrophoresis         from potatoes         Standard for Starch Assay Kits SA-20 and STA-20         practical grade         Unmodified regular corn starch containing approx. 73% amylopectin and 27% amylose         for electrophoresis         Soluble         Soluble         Powder	\$4001-10MG \$4001-100MG \$4001-50 85645-100G 85645-1KG 85645-1KG 85650-1KG \$5296-5G \$4180-100G \$4180-100G \$4180-100G \$4180-1KG \$4126-2KG \$4126-5KG \$5651-10KG \$5651-2KG \$5651-2KG \$2630-100G \$2630-100G \$2630-1KG \$2004-500G \$2004-1KG \$4251-2KG
Starch Starch Starch from corn Starch from corn Starch from corn Starch from potato	from potatoes, for electrophoresis         from potatoes         Standard for Starch Assay Kits SA-20 and STA-20         practical grade         Unmodified regular corn starch containing approx. 73% amylopectin and 27% amylose         for electrophoresis         Soluble         Soluble         Powder	\$4001-10MG           \$4001-100MG           \$4001-500MG           \$4001-5G           85645-100G           85645-1KG           85650-1KG           \$5296-5G           \$4180-100G           \$4180-100G           \$4180-500G           \$4180-100G           \$4126-2KG           \$5651-100G           \$5651-2KG           \$2630-100G           \$2630-100G           \$2630-100G           \$2630-10KG           \$2004-500G           \$2004-500G           \$2004-1KG           \$4251-2KG           \$4251-2KG           \$7260-100G
Starch Starch from corn Starch from corn Starch from corn Starch from potato Starch from potato Starch from potato Starch from potato Starch from potato	from potatoes, for electrophoresis         from potatoes         Standard for Starch Assay Kits SA-20 and STA-20         practical grade         Unmodified regular corn starch containing approx. 73% amylopectin and 27% amylose         for electrophoresis         Soluble         Soluble         Powder	\$4001-10MG \$4001-100MG \$4001-500MG \$4001-5G 85645-100G 85645-1KG 85650-1KG \$5296-5G \$4180-100G \$4180-100G \$4180-500G \$4180-1KG \$4126-2KG \$4126-2KG \$5651-10KG \$5651-2KG \$5651-2KG \$2630-100G \$2630-1KG \$2004-500G \$2004-1KG \$4251-2KG \$4251-2KG \$4251-2KG \$4251-2KG \$7260-100G \$7260-500G
tarch tarch tarch tarch from corn tarch from corn tarch from corn tarch from potato tarch from rice	from potatoes, for electrophoresis         from potatoes         Standard for Starch Assay Kits SA-20 and STA-20         practical grade         Unmodified regular corn starch containing approx. 73% amylopectin and 27% amylose         for electrophoresis         Soluble         Soluble         -	\$4001-10MG           \$4001-10MG           \$4001-16           \$4001-50           85645-100G           85645-1KG           85645-1KG           85650-1KG           \$5296-5G           \$4180-100G           \$4180-100G           \$4180-100G           \$4126-2KG           \$5651-100G           \$5651-100G           \$5651-100G           \$5651-100G           \$5651-100G           \$5651-2KG           \$2630-100G           \$2630-100G           \$2630-100G           \$2630-100G           \$2630-100G           \$2630-100G           \$2630-100G           \$2630-100G           \$2630-100G           \$2640-100G           \$2760-100G           \$7260-100G           \$7260-100G           \$7260-100G           \$7260-100G           \$7260-100G

Product Name	Description	Cat. No.
Starch from wheat	Purified	\$2760-500G
Starch from Azure	Potato starch covalently linked with Remazol Brilliant Blue	\$7629-1G
		\$7629-5G
		\$7629-25G
Starch from Azure	Insoluble corn starch covalently linked with Remazol Brilliant Blue	\$7776-1G
		\$7776-5G
Starch, from soluble	ACS reagent	\$9765-100G
		\$9765-250G
		\$9765-500G
		S9765-1KG
Xylan from beechwood		X4252-10G
		X4252-25G
		X4252-100G
Xylan from birch wood	Xylose residues: ≥90%	X0502-10G
		X0502-25G
		X0502-100G
Xylan from birch wood		95588-10G
		95588-25G
		95588-100G
Xylan from oat spelts	Xylose: ≥70%	X0627-10G
		X0627-25G
		X0627-100G
Xylan from oat spelts		95590-10G
		95590-50G





Your gateway to products, services and more for Life Science Research.

Read current and previous issues and register to receive future BioFiles issues.

Visit us at sigma-aldrich.com/biofiles

sigma-aldrich.com/biofiles

**SIGMA** 

### Argentina

SIGMA-ALDRICH DE ARGENTINA S.A. Free Tel: 0810 888 7446 Tel: (+54) 11 4556 1472 Fax: (+54) 11 4552 1698

### Australia

SIGMA-ALDRICH PTY LTD. Free Tel: 1800 800 097 Free Fax: 1800 800 096 Tel: (+61) 2 9841 0555 Fax: (+61) 2 9841 0500

### Austria

SIGMA-ALDRICH HANDELS GmbH Tel: (+43) 1 605 81 10 Fax: (+43) 1 605 81 20

### Belaium

SIGMA-ALDRICH NV/SA. Free Tel: 0800 14747 Free Fax: 0800 14745 Tel: (+32) 3 899 13 01 Fax: (+32) 3 899 13 11

### Brazil

SIGMA-ALDRICH BRASIL LTDA. Free Tel: 0800 701 7425 Tel: (+55) 11 3732 3100 Fax: (+55) 11 5522 9895

### Canada

SIGMA-ALDRICH CANADA LTD. Free Tel: 1800 565 1400 Free Fax: 1800 265 3858 Tel: (+1) 905 829 9500 Fax: (+1) 905 829 9292

### China

SIGMA-ALDRICH (SHANGHAI) TRADING CO. LTD. Free Tel: 800 819 3336 Tel: (+86) 21 6141 5566 Fax: (+86) 21 6141 5567

### **Czech Republic**

SIGMA-ALDRICH S.R.O. Tel: (+420) 246 003 200 Fax: (+420) 246 003 291 Denmark SIGMA-ALDRICH DENMARK A/S Tel: (+45) 43 56 59 10 Fax: (+45) 43 56 59 05

SIGMA-ALDRICH FINLAND OY

SIGMA-ALDRICH CHIMIE S.à.r.I.

SIGMA-ALDRICH CHEMIE GmbH

Tel: (+358) 9 350 9250

Fax: (+358) 9 350 92555

Free Tel: 0800 211 408

Free Fax: 0800 031 052

Tel: (+33) 474 82 28 00

Fax: (+33) 474 95 68 08

Free Tel: 0800 51 55 000

Free Fax: 0800 64 90 000

Fax: (+49) 89 6513 1160

Tel: (+30) 210 994 8010

Fax: (+30) 210 994 3831

SIGMA-ALDRICH Kft

Tel: (+36) 1 235 9055

Fax: (+36) 1 235 9050

PRIVATE LIMITED

SIGMA-ALDRICH (O.M.) LTD.

Ingyenes zöld telefon: 06 80 355 355

Ingyenes zöld fax: 06 80 344 344

SIGMA-ALDRICH CHEMICALS

Bangalore: (+91) 80 6621 9600

New Delhi: (+91) 11 4165 4255

Hyderabad: (+91) 40 6684 5488

Bangalore: (+91) 80 6621 9650

New Delhi: (+91) 11 4165 4266

Hyderabad: (+91) 40 6684 5466

Mumbai: (+91) 22 2579 7589

Mumbai: (+91) 22 2570 2364

Tel: (+49) 89 6513 0

Finland

France

Germany

Greece

Hungary

India

Telephone

Fax

SIGMA-ALDRICH IRELAND LTD. Free Tel: 1800 200 888 Free Fax: 1800 600 222 Tel: (+353) 1 404 1900 Fax: (+353) 1 404 1910

Ireland

### Israel

SIGMA-ALDRICH ISRAEL LTD. Free Tel: 1 800 70 2222 Tel: (+972) 8 948 4100 Fax: (+972) 8 948 4200

### Italy

SIGMA-ALDRICH S.r.I. Numero Verde: 800 827018 Tel: (+39) 02 3341 7310 Fax: (+39) 02 3801 0737

### Japan SIGMA-ALDRICH JAPAN K.K. Tokyo Tel: (+81) 3 5796 7300

Tokyo Fax: (+81) 3 5796 7315 Korea SIGMA-ALDRICH KOREA Free Tel: (+82) 80 023 7111 Free Fax: (+82) 80 023 8111

### Tel: (+82) 31 329 9000 Fax: (+82) 31 329 9090 Malaysia SIGMA-ALDRICH (M) SDN. BHD

Tel: (+60) 3 5635 3321 Fax: (+60) 3 5635 4116

### Mexico

Order/Customer Service (800) 325-3010 • Fax (800) 325-5052

Technical Service (800) 325-5832 • sigma-aldrich.com/techservice

Development/Bulk Manufacturing Inquiries SAFC<sup>®</sup> (800) 244-1173

SAFC

SUPELCO<sup>®</sup>

SIGMA-ALDRICH QUÍMICA, S.A. de C.V. Free Tel: 01 800 007 5300 Free Fax: 01 800 712 9920 Tel: 52 722 276 1600 Fax: 52 722 276 1601

### The Netherlands

SIGMA-ALDRICH CHEMIE BV Free Tel: 0800 022 9088 Free Fax: 0800 022 9089 Tel: (+31) 78 620 5411 Fax: (+31) 78 620 5421

### New Zealand

SIGMA-ALDRICH NEW ZEALAND LTD. Free Tel: 0800 936 666 Free Fax: 0800 937 777 Tel: (+61) 2 9841 0555 Fax: (+61) 2 9841 0500

### Norway

SIGMA-ALDRICH NORWAY AS Tel: (+47) 23 17 60 60 Fax: (+47) 23 17 60 50

### Poland

SIGMA-ALDRICH Sp. z o.o. Tel: (+48) 61 829 01 00 Fax: (+48) 61 829 01 20

### Portugal

SIGMA-ALDRICH QUÍMICA, S.A. Free Tel: 800 202 180 Free Fax: 800 202 178 Tel: (+351) 21 924 2555 Fax: (+351) 21 924 2610

### Russia

SIGMA-ALDRICH RUS, LLC Tel: +7 (495) 621 6037 Fax: +7 (495) 621 5923

### Singapore SIGMA-ALDRICH PTE. LTD. Tel: (+65) 6779 1200 Fax: (+65) 6779 1822

### South Africa

SIGMA-ALDRICH SOUTH AFRICA (PTY) LTD. Free Tel: 0800 1100 75 Free Fax: 0800 1100 79 Tel: (+27) 11 979 1188 Fax: (+27) 11 979 1119

ISOTEC"

### Free Fax: 0800 378 785 Tel: (+44) 1747 833 000

Spain

Sweden

SIGMA-ALDRICH QUÍMICA, S.A.

Free Tel: 900 101 376

Free Fax: 900 102 028

Tel: (+34) 91 661 99 77

Fax: (+34) 91 661 96 42

Tel: (+46) 8 742 4200

Fax: (+46) 8 742 4243

Free Tel: 0800 80 00 80

Free Fax: 0800 80 00 81

Tel: (+41) 81 755 2828

Fax: (+41) 81 755 2815

United Kingdom

Free Tel: 0800 717 181

Switzerland

SIGMA-ALDRICH SWEDEN AB

SIGMA-ALDRICH CHEMIE GmbH

Fax: (+44) 1747 833 313 SAFC (UK) Free Tel: 0800 71 71 17 **United States** SIGMA-ALDRICH

SIGMA-ALDRICH COMPANY LTD.

### P.O. Box 14508 St. Louis, Missouri 63178 Toll-Free: 800 325 3010 Toll-Free Fax: 800 325 5052 Call Collect: (+1) 314 771 5750 Tel: (+1) 314 771 5765 Fax: (+1) 314 771 5757

Internet sigma-aldrich.com



### Accelerating Customers' Success through Leadership in Life Science, High Technology and Service

01948-504141

### World Headquarters 3050 Spruce St., St. Louis, MO 63103

(314) 771-5765 sigma-aldrich.com

The SIGMA-ALDRICH Group



©2007 Sigma-Aldrich Co. All rights reserved

©2007 Sigma-Aldrich Co. All rights reserved. SIGMA  $\frac{2}{8}$ , SAFC, SAFC, SIGMA-ALDRICH,  $\frac{2}{8}$  ISOTEC, ALDRICH,  $\frac{2}{9}$ , FLUKA,  $\frac{1}{9}$ , and SUPELCO are trademarks belonging to Sigma-Aldrich Co. and its affiliate Sigma-Aldrich Biotechnology, L.P. Riedel-de Haen<sup>®</sup> trademark under license from Riedel-de Haen GmbH. AMG 300L, BAN, Fructozyme L, are trademarks; Carezyme 1000L, Celluctast 1.5L, Pentopan Mono BG, Termamyl, Fungamyl, Duramyl, Promozyme, Pectinex Ultra SPL, are registered trademarks of Novozymes Corp. Ficoll is a registered trade-mark of Amersham Biosciences. Ficoll is a registered trademark of GE Healthcare. Sigma Fand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.

Riedel-de Haën\*



3050 Spruce Street • St. Louis, MO 63103 USA



///Fluka