

# THE REVOLUTION REVOLUTION NEEDEN ADDIA

The **ReadyStream® system** enter the new era of media preparation

> MilliporeSigma is the U.S. and Canada Life Science business of Merck KGaA, Darmstadt, Germany.

# **Millipore**®

Preparation, Separation, Filtration & Monitoring Products

# Media preparation DAY TOMORROW TODAY

### **Autoclave/Media Preparator**

### Classical culture media production from dehydrated culture media using autoclave or Media preparator





Time consumina



Need to clean lab equipment







Organizational effort to pre-warm media

High water and energy consumption

Significant investment & intensive

maintenance





High plastic waste generation

Non-sustainable transportation



XGDD

**X12D** bags 



VA

M



Require

substantial

storage space

### **ReadyStream® system**



Save space in your storage area and lab. Rethink how you prepare & dispense media with the revolutionary ReadyStream<sup>®</sup> system

# 100 L of preheated, ISO compliant Culture Media, prepared in only 20 minutes

# HOW is that possible?

Installation of the **ReadyStream**<sup>®</sup> bag and filter as well as instrument calibration in about 5 minutes



The **ReadyStream**<sup>®</sup> bags are irradiated, come pre-filled with our superior **GranuCult**<sup>®</sup> prime dehydrated culture media.

The disposable **ReadyStream**<sup>®</sup> filter set contains a 0.22 µm sterile filter allowing purification of 600 L ISO 11133 compliant water.

# 15 minutes for rehydrating the media bag



Water delivered through the **ReadyStream**<sup>®</sup> system rehydrates the media, generating a 10x concentrated media syrup.

By scanning the data matrix code on the consumables, the system tracks the media type, consumable's capacity, and shelf life.\*

### 100 L media ready to dispense

Thanks to the 10x media concentrate this bag delivers **100 L media**.\*\*

At the push of a button, the media is dispensed at the **required volume**, **temperature**, **concentration** and **flow rate** – all adjustable in presets.

\* Shelf life of **ReadyStream® TSB** is 10 days and for **ReadyStream® BPW** 5 days after first rehydration. \*\* **For ReadyStream® BPW**, smaller bags delivering 30 L of media are available too.

step

step

READY

# choose your future lab setup FROM 3 OPTIONS

# Standard ReadyStream® Setup - High Volume Users

- Simple setup, easy to install
- Can be used for any sample size but in particular suited for large sample sizes and dispensing volumes e.g., 3375 mL



# In-Line Connection -Connect it to any dilutor!

- Convenient connection to 1 or 2 dilutors at the same time
- Suited for small sample sizes, e.g., 25 g
- No change of connections required
- Addition of supplements possible



# Flex Connect -For separated media kitchen and sample prep locations

- Enhanced flexibility thanks to quick connection/disconnection
- Prepare media and dispense it independently from instrument location
  - Pre-filling of several bottles and temporary storage or delivery at different dilutor positions possible
    - In particular suited for small sample sizes
    - Addition of supplements possible

# Enhancing sustainability and operational Efficiency





# Energy & Water:

No need for resource intensive autoclaving and lab equipment cleaning, enabling significant savings in energy and water consumption compared to traditional media culture media production from dehydrated culture media.



# Time:

Efficiently prepare large quantities of preheated media at the push of a button, with the smart and innovative ReadyStream<sup>®</sup> process, saving significant time and effort compared to conventional culture media solutions.



# Waste:

Reduce waste compared to readyto-use media: 1 ReadyStream<sup>®</sup> bag delivers same media volume like 10 classic ready-to-use media bags.

The ReadyStream<sup>®</sup> platform is awarded with the green icon for greener alternative products.



More information on this program can be found here Data derived from a study investigating food testing processes from media preparation to primary dilution, executed by an independent external food testing (BAV INSTITUT GmbH, A Tentamus Company).

To access full study report scan or click


on the QR code.




1500 -62% -67% -35% -85% 1000 Traditional - Media Kitchen and Autoclave Ready-to-Use (RTU) Media Bags Media Preparator Instrument ReadyStream® System with Diluto ReadyStream® System Direct Dispensing Hands-on time (sec) per 375 g sample -24% 500 400 300 200 Media Preparator Inst Ready-to-Use (RTU) Media Bao ReadyStream<sup>®</sup> System with Dilutor Traditional - Media Kitchen and Autoclay


Total processing time (sec) per 375 g sample

ReadyStream<sup>®</sup> System Direct Dispensing



ReadyStream<sup>®</sup> System Direct Dispensing





# With the ReadyStream® system, you can

- in about 20 minutes prepare up to 100 L of media
- dispense pre-heated media of the volume, temperature and concentration you need
- easily store and carry the bags containing dehydrated media (less than 3 kg for the biggest bag)
- save energy, water and waste to meet your sustainability goals

And there's no need for autoclaving, bottles and washing, saving you time that surely can be well deployed for higher value tasks.

Don't wait to try our ReadyStream® system,

# **Ordering information**

Product	Packaging	Ordering No.
30 L ReadyStream <sup>®</sup> Media Bag GranuCult <sup>®</sup> Buffered Peptone Water	10x 30 L Media bags with dehydrated media*	5.74846.0030
100 L ReadyStream <sup>®</sup> Media Bag GranuCult <sup>®</sup> Buffered Peptone Water	3x 100L Media bags with dehydrated media*	5.74846.0100
100 L ReadyStream <sup>®</sup> Media Bag GranuCult <sup>®</sup> Tryptic Soy Broth	3x 100L Media bags with dehydrated media*	5.77000.0100
ReadyStream <sup>®</sup> Filter Set	5x Filter sets	5.74826.0001
ReadyStream <sup>®</sup> System	1x	RDYSTRM01
Female Connector for ReadyStream <sup>®</sup> Bag (for Dilutor Connection)	1x 5x	RDYFCON01 RDYFCON05

\*Bags are shipped with dry media which will be reconstituted with the ReadyStream<sup>®</sup> system. The 10x concentrated media allows dispensing of up to 30 L and 100 L, respectively, of 1x concentrated media with the different bag sizes.



To place an order or receive technical assistance in Europe, please call Customer Service:France: 0825 045 645Spain: 901 516 645 Option 1Germany: 069 86798021Switzerland: 0848 645 645Italy: 848 845 645United Kingdom: 0870 900 4645

For other countries across Europe, please call: +44 (0) 115 943 0840 Or visit: **SigmaAldrich.com/offices** For Technical Service visit: **SigmaAldrich.com/techservice** 

Millipore. Sigma-Aldrich. Supelco. Milli-Q. SAFC. BioReliance. © 2024 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, the vibrant M, BioReliance, Millipore, Milli-Q. SAFC, Sigma-Aldrich, Supelco, GranuCult, DiluCult and ReadyStream are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

MK\_BR8340EN Ver. 5.0 54070 04/2024