MUELLER-HINTON agar CLSI

For the determination of antibiotic susceptibility including sulfonamides by the disc-agar diffusion method.

IVD

In Vitro Diagnostic Medical Device -

For professional use only



Version 17-10-2008 Merck KGaA, 64271 Darmstadt

Droparat

Principle Microbiological

Microbiological method.

Mode of Action

The medium complies with the requirement of the Clinical and Laboratory Standards Institute (CLSI) and is manufactured to contain low concentrations of thymine and thymidine as well as appropriate levels of calcium and magnesium ions.

Thymine and thymidine concentrations are determined by the disc diffusion procedure with trimethoprim and sulfamethoxazole and Enterococcus faecalis ATCC 29212.

Calcium and/or magnesium concentrations are controlled by obtaining the correct zone diameters with aminoglycoside antibiotics and Pseudomonas aeruginosa ATCC 27853.

Typical Composition (g/litre)

Meat infusion 2.0; casein hydrolysate 17.5; starch 1.5; agar-agar 17.0..

See also General Instruction for Use "How to use Dehydrated Culture Media"

For MSDS, warnings and precautions see our website: www.merck-chemicals.com

Preparation

Suspend 38 g in 1 litre of demin. water and autoclave (15 min at 121 °C). If required, cool to 45-50 °C and add 5-10 % defibrinated blood, pour plates.

pH: 7.3 ± 0.2 at $25 \, ^{\circ}$ C.

Without blood, the plates are clear to opalescent and brownishyellow.

Storage

Usable up to the expiry date when stored $\,$ dry and tightly closed at +15 to +25 °C. Protect from light.

After first opening of the bottle the content can be used up to the expiry date when stored dry and tightly closed at +15 to +25 °C.

Specimen

e.g. Isolated bacteria from urine.

Clinical specimen collection, handling and processing, see general instructions of use.

Experimental Procedure

Carry out the sensitivity test acc. to CLSI.

Incubation for 24 h at 35 °C aerobically.

Literature

National Committee for Clinical Laboratory Standards. Approved Standard. Performance standards for antimicrobial disc susceptibility tests, 5th ed. Natrional Committee for Clinical Laboratory Standards, Villanova, Pa. (1993).

Ordering Information

Product	Ordering No.	Pack size
MUELLER-HINTON agar acc. to CLSI	1.05435.0500	500 g

Quality control

Test strains	Growth	Zone diameters within specifications
Escherichia coli ATCC 25922	good / very good	+
Staphylococcus aureus ATCC 25923	good / very good	+
Pseudomonas aeruginosa ATCC 27853	good / very good	+
Enterococcus faecalis ATCC 29212	good / very good	+