

media selection + Filter sizing = Ready to Run

We do the testing to identify the best filter for your custom media formulation.



Custom Media Screening

imMEDIAte ADVANTAGE™ (IA) for Cell Culture Media



Media Selection

Highest performing media is chosen



- Filter selection to meet your need
- Filter sizing

Titter recommendation example		
Batch Volume (L)	Mycoplasma and Bacterial Retention	Virus, Mycoplasma and Bacterial Retention
	Millipore Express® SHR-P (0.5/0.1 µm) filters	Viresolve® Barrier filters
< 100	XL 3 KHVEG03TT3	0.15 m ² VBKG015TC1
100 - 1000	XL 10 KHVEG10TT1	1.0 m ² VBKG100TC1
1000 - 2000	XL 10 KHVEG10TT1	2 x 1.0 m ² VBKG100TC1

Filter recommendation example



Performance delivered: custom media with the optimal filter



Ready for GMP Production

For more information, please contact the Tech Service group.

The life science business of Merck KGaA, Darmstadt, Germany sperates as MilliporeSigma in the J.S. and Canada. Millipore

Preparation, Separation, Filtration & Testing Products